

Technical Specification Sheet

Product Name: **ELASTOMERIC LIGATURE TIES** – are a high-performance elastic ligature with superior elastic retention.

Available in 49 colors.

Distributor: Orthodontal USA
1500 Weston Rd Suite 204, Weston, FL 33326, EE.UU

Phone: (+1) 9546016832
Email: sales@orthodontalusa.com

Product Name:	ELASTOMERIC LIGATURE TIES
Indication for use:	Ligature Tie
Compositions:	Polymers (Latex-Free)
Colors:	49 colors
Inner Diameter	0.060" ± 0.005"
Outer Diameter	0.125" ± 0.01"
Cross Section Thickness	0.060" ± 0.005"
Initial Extension Force	90 ± 10 gm
24 h Residual Force	50 ± 10%
Composition	Polyurethane Polyether Elastomer - 96-98 % FDA approved pigments 2-4%

Good Manufacturing Practices and is ISO 13485:2016 certified.