

Technical Specification Sheet

Product Name: **ELASTOMERIC POWER CHAIN** are a high performance elastic ligature with superior elastic retention capable of effecting superior anatomical tooth movement.

It has excellent tear strength minimizing chain breakage during placement and treatment.

Available in 46 colors.

Distributor: Orthodontal USA
1500 Weston Rd Suite 204, Weston, FL 33326, EE.UU

Phone: (+1) 9546016832
Email: sales@orthodontalusa.com

Product Name:	ELASTOMERIC POWER CHAIN
Indication for use:	Ligature Tie
Compositions:	Polymers (Latex-Free)
Colors:	49 colors
Inner Diameter	0.050" ± 0.005"
Link Lengths	Continuous 0.110", Short 0.132" and Long 0.165"
Cross Section Thickness	0.060" ± 0.005"
Initial Extension Force	250-300 g
Thickness	.027" ± 0.001"
Composition	Polyurethane Polyether Elastomer - 96-98 % FDA approved pigments 2-4%

Good Manufacturing Practices and is ISO 13485:2016 certified.