



MATHIEU NEEDLE HOLDER NARROW FINE TC - GOLD HANDLE

SKU:
N21150465F T/C

A-1 TECHNICAL OVERVIEW



FINE PRECISION JAWS

Slim, tapered jaws provide precise access and controlled handling during delicate orthodontic and dental procedures.



SPRING-ACTION HANDLE

Flexible spring-style handle provides smooth opening and closing action for comfortable, controlled instrument operation.

** All measurement are approximate*

www.orthodentalusa.com



MATHIEU NEEDLE HOLDER NARROW FINE TC - GOLD HANDLE

SKU:
N21150465F T/C

A-2 KEY FEATURES



PRECISION ENGINEERED

Designed for accurate and reliable measurements.



REUSABLE & DURABLE

Built to last with repeated sterilization and daily use.



EASY TO CLEAN & STERILIZE

Designed for easy cleaning and fully autoclavable, ensuring the highest standards of hygiene in your practice.



HIGH QUALITY STAINLESS STEEL

Made from premium stainless steel for superior strength and corrosion resistance.



AUTOCLAVABLE

Fully autoclavable and corrosion resistant for reliable sterilization.



** All measurement are approximate*



MATHIEU NEEDLE HOLDER NARROW FINE TC - GOLD HANDLE

SKU:
N21150465F T/C

A-3 TECHNICAL SPECIFICATIONS

SKU	N21150465F T/C
Product Name	Mathieu Needle Holder arrow Fine TC
Instrument Type	Surgical Needle Holder
Function	Holding and Controlling Surgical Suture Needles
Material	Premium Surgical Stainless Steel
Jaw Design	Narrow Fine Precision Jaws
Finish	Satin Finish
Insert Material	Tungsten Carbide (TC)
Handle Design	Mathieu Spring Handle with Locking Mechanism
Sterilization Temp	Up to 135°C (275°F)
Corrosion Resistant	Yes
Sterilization	Fully Autoclavable



** All measurement are approximate*