

 **MORELLI**

**the best
version of the
*Smile***



**Product
Catalog**



**QUALITY
GOVERNANCE AND
CONTINUITY
CORPORATE**

OUR HISTORY

Morelli Orthodontia was born in 1980 from an innovative idea, in the city of Sorocaba, in São Paulo. With more than 40 years of experience, Morelli has established itself as one of the largest manufacturers of orthodontic products in the world. Today, occupying an area of 20 thousand m², where monthly More than 20 million pieces are produced from a portfolio of more than 2,000 items.

Based on entrepreneurship, Morelli is a pioneer in Brazilian industry of orthodontic products. The company emerged in a conflictive scenario marked by political and economic instability. In the middle of this panorama, the innovative and entrepreneurial vision of our founders, allowed Morelli to meet the needs of the national orthodontic market. Today, with a network of more than 200 distributors in Brazil and presence in more than 30 countries, Morelli strengthens its presence in the market every day. The constant technological evolution of our products is one of the pillars of our philosophy, guaranteeing a high degree of quality and excellence.

A SOLID AND RELIABLE FUTURE

Morelli has demonstrated strength and confidence throughout more than four decades of growth, focusing on the long-term vision.

In 2020, it began a business succession process under the visionary leadership of its founder, Mr. Oraci Morelli. This process included in-depth studies on the challenges and opportunities of succession, with consulting specialized in succession management. **Corporate governance and business restructuring processes were implemented**, including the creation of an experienced Advisory Board to guide decisions.

At Morelli, excellence and continuity are non-negotiable pillars, and the commitment to clients and partners is very clear: the future is prepared with care and dedication, and the company is prepared to face it and shape it in favor of dentistry.

This is Morelli:

Continuous dedication, to always be Morelli.

Quick guide to Morelli Kits


Prescription	Model	Canines Premolars		Upper Canine Angle	1 case (20 pieces)		10 case (200 pieces)		100 case (2000 pieces)		Page
					22	18	22	18	22	18	
ROTH	Autoligado SLI	✓	✓	9°	10.14.900						11
	Autoligado SLP	✓	✓	9°	10.13.900						11
	Autoligado SLH	✓	✓	9°	10.23.900						11
	Ceramic SLI	✓	✓	9°	10.12.900						09
	Ceramic SLP	✓	✓	9°	10.16.900						09
	Ceramic SLH	✓	✓	9°	10.26.912						09
	Polysafira	✓	✓	9°	10.24.900						15
	Ceramic	✓	✓	9°	10.11.900						15
	Max	✓	✓	9°	10.15.903	10.15.902	10.15.933	10.15.932	10.15.923	10.15.922	21
	Max	✓	✓	13°	10.15.905	10.15.904	10.15.935	10.15.934	10.15.925	10.15.924	21
	MaX (Bidimensional)	✓	✓	9°	10.15.930						21
	Nano	✓	✓	9°	10.15.906						24
	Monobloc	✓	✓	13°	10.15.901						24
	Light			13°	10.10.907	10.10.906					25
	Light			9°	10.10.915						25
	Light	✓		13°	10.10.961	10.10.960					26
	Light	✓		9°	10.10.963						25
	Light	✓	✓	13°	10.10.979	10.10.978	10.10.996				26
	Light	✓	✓	9°	10.10.981		10.10.997				26
	Light (Biofuncional)	✓	✓	13°	10.10.920						27
	Standard			13°	10.10.901	10.10.900					28
	Standard			9°	10.10.911						28
	Standard	✓		13°	10.10.951	10.10.950					29
	Standard	✓		9°	10.10.953						29
	Standard	✓	✓	13°	10.10.969						30
	Standard	✓	✓	9°	10.10.971						29
	Composite	✓	✓	13°	10.17.901		10.17.932				15
EDGEWISE	Ceramic	✓	✓		10.18.901						16
	Slim Max				10.65.905	10.65.904	10.65.935	10.65.934			23
	Slim Max	✓			10.65.907	10.65.906	10.65.937	10.15.936			24
	Slim	✓			10.65.901	10.65.900					31
	Standard				10.30.901	10.30.900					31
M.B.T	Autoligado SLI	✓	✓		10.27.900						12
	Autoligado SLP	✓	✓		10.28.900						12
	Autoligado SLH	✓	✓		10.23.935						12
	Ceramic	✓	✓		10.32.900						16
	Polysafira	✓	✓		10.33.900						16
	Max		✓		10.35.906		10.35.936				23
	Max	✓	✓		10.35.905		10.35.935				23
	Light				10.35.904						27
	Light	✓	✓		10.35.903						28
	Light (Bidimensional)	✓			10.35.924						27
	Standard				10.35.901						30
	Standard	✓			10.35.902						30
CAPELOZZA	Max Padrão I	✓	✓		10.19.904		10.19.934				22
	Max Padrão II	✓	✓		10.19.905		10.19.935				22
	Max Padrão III	✓	✓		10.19.906		10.19.936				22
RICKETTS	Actual	✓	✓			10.25.901					32
	Actual Especial	✓	✓			10.25.900					32
	Slim	✓	✓			10.25.905					32

TABLE OF CONTENTS

Brackets	07
Aesthetic Brackets	13
Metal Brackets	17
Tubes	35
Accessories	43
Placas de Termoformado	53
Thermoformed Plates	55
Bands	59
Archs	65
Wires	73
Elastic	77
Functional Orthopedics	85
Headcaps and Tractioners	91
Instruments	95
Oral Hygiene	105
Compartments	109
Adhesivos	113


ORTHODONTIC SPECIFICATIONS

GUIDE




TORQUE

CERVICAL OU GENGIVAL




CERVICAL OU GENGIVAL

SLOT




ÂNGULO

CERVICAL OU GENGIVAL



CERVICAL OU GENGIVAL

IN / OUT (IN-SET / OFF-SET)




COMPRIMENTO




DISTAL OFFSET



CURVATURA DA BASE



ARCADA DENTÁRIA



Self-Ligating Brackets





CLIP EM NITIOLO

La aleación de níquel titanio tiene memoria elástica, pudiendo resistir deformaciones mecánicas durante todo el tratamiento. Debido a la memoria geométrica del material, el Clip siempre tenderá a volver a su posición original, asegurando una mejor lectura de torque.

ACABADO MATE

Para complementar la estética de la pieza, el Clip tiene un acabado mate, y al no ser un revestimiento, no se despega y permanece inalterable hasta el final del tratamiento.

MECANISMO SIMPLES Y PRÁCTICO

Sin necesidad de herramientas especiales o alicates.

SISTEMA ANTIBLOQUEO

Especialmente diseñado para evitar que el clip se bloquee debido a la acumulación de placa y sarro.



Video Ceramic SLI

MÁXIMA COMODIDAD

Por sus reducidas dimensiones y contornos redondeados.



ESTÉTICA SUPERIOR

Injectados en cerámica alúmina pura, los brackets están fabricados con la más moderna Tecnología CIM (Ceramic Injection Molding). Las propiedades translúcidas de la pieza garantizan una estética superior y una mejor adaptación a los colores dentales.

EXCELENTE ADAPTACIÓN Y RETENCIÓN

Con geometría optimizada, la base tiene curvaturas mesio-distal y cérvico-oclusal, que aseguran una mejor adaptación al diente. La base está bigranallada y silanizada, lo que garantiza una excelente retención mecánica y química.

DESLICE OPTIMO

Ranura con alivio de curvatura y extremos redondeados.

Morelli Ceramic SLI Prescription Roth with hooks in the canines (9°) and premolars



3.3mm	3.3mm	3.3mm	2.9mm	3.6mm	3.6mm	2.9mm	3.3mm	3.3mm	3.3mm	◀ Length
-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°	◀ Torque
0°	0°	+9°	+9°	+5°	+5°	+9°	+9°	0°	0°	◀ Angle
10.12.007	10.12.007	10.12.005	10.12.003	10.12.001	10.12.002	10.12.004	10.12.006	10.12.008	10.12.008	◀ Replacement .022" (2 UN)
15	14	13	12	11	21	22	23	24	25	
45	44	43	42	41	31	32	33	34	35	



1 CASE - 20 PIECE KIT
SLOT22 REF 10.12.900

3.0mm	3.0mm	3.3mm	2.7mm	2.7mm	2.7mm	2.7mm	3.3mm	3.0mm	3.0mm	◀ Length
-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°	◀ Torque
0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°	◀ Angle
10.14.014	10.14.012	10.12.011	10.12.009	10.12.009	10.12.009	10.12.009	10.12.010	10.14.013	10.14.015	◀ Replacement .022" (2 UN)

Morelli Ceramic SLP Prescription Roth with hook on canines (9th) and premolars



3.3mm	3.3mm	3.3mm	2.9mm	3.6mm	3.6mm	2.9mm	3.3mm	3.3mm	3.3mm	◀ Length
-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°	◀ Torque
0°	0°	+9°	+9°	+5°	+5°	+9°	+9°	0°	0°	◀ Angle
10.16.007	10.16.007	10.16.005	10.16.003	10.16.001	10.16.002	10.16.004	10.16.006	10.16.008	10.16.008	◀ Replacement .022" (2 UN)
15	14	13	12	11	21	22	23	24	25	
45	44	43	42	41	31	32	33	34	35	



1 CASE - 20 PIECE KIT
SLOT22 REF 10.16.900

3.0mm	3.0mm	3.3mm	2.7mm	2.7mm	2.7mm	2.7mm	3.3mm	3.0mm	3.0mm	◀ Length
-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°	◀ Torque
0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°	◀ Angle
10.13.014	10.13.012	10.16.011	10.16.009	10.16.009	10.16.009	10.16.009	10.16.010	10.13.013	10.13.015	◀ Replacement .022" (2 UN)

Morelli Ceramic SLH Prescription Roth with hook on canines (9th) and premolars



3.3mm	3.3mm	3.3mm	2.9mm	3.6mm	3.6mm	2.9mm	3.3mm	3.3mm	3.3mm	◀ Length
-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°	◀ Torque
0°	0°	+9°	+9°	+5°	+5°	+9°	+9°	0°	0°	◀ Angle
10.16.007	10.16.007	10.16.005	10.12.003	10.12.001	10.12.002	10.12.004	10.16.006	10.16.008	10.16.008	◀ Replacement .022" (2 UN)
15	14	13	12	11	21	22	23	24	25	
45	44	43	42	41	31	32	33	34	35	



1 CASE - 20 PIECE KIT
SLOT22 REF 10.26.912

3.0mm	3.0mm	3.3mm	2.7mm	2.7mm	2.7mm	2.7mm	3.3mm	3.0mm	3.0mm	◀ Length
-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°	◀ Torque
0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°	◀ Angle
10.13.014	10.13.012	10.16.011	10.12.009	10.12.009	10.12.009	10.12.009	10.16.010	10.13.013	10.13.015	◀ Replacement .022" (2 UN)

*Peças 4 e 5 inferiores metálicas para melhor resistência aos esforços de mastigação.
Reposição com 2 unidades são para as peças cerâmicas e 5 unidades para as peças metálicas.

BRACKET AUTOLIGABLE MORELLI

Nº1 LÍDER EN LA CATEGORÍA

La línea Morelli Autoligable fue especialmente desarrollada para promover una excelente mecánica de deslizamiento, con comodidad al paciente y sentido práctico al ortodoncista.

CLIP EM NITINOL

La aleación de níquel titanio tiene memoria elástica, pudiendo resistir deformaciones mecánicas durante todo el tratamiento. Debido a la memoria geométrica del material, el Clip siempre tenderá a volver a su posición original, asegurando una mejor lectura de torque

MECANISMO SIMPLES Y PRÁCTICO

Un sistema práctico que no requiere el uso de herramientas especiales o alicates. El explorador de tragamonedas hace que abrir el Clip sea una tarea simple y conveniente. Y para cerrar, basta con empujar el Clip

SISTEMA ANTIBLOQUEIO

El diseño del Clip fue especialmente diseñado para evitar que la acumulación de placa bacteriana provoque que el Clip se pegue. Facilitando el intercambio de arcos durante el tratamiento

DISEÑO REDONDEADO

Cantos redondeados y perfectamente pulidos, que garantizan comodidad al paciente.

SLOT CON ALIVIO DE CURVATURA Y EXTREMIDADES REDONDEADAS

Reducción de los efectos de unión y muescas, la fricción causada por la flexión del alambre se suaviza con los extremos redondeados y los relieves internos. Estas características garantizan una excelente mecánica de



RESISTENCIA SUPERIOR

Injectado por la más moderna tecnología MIM (Metal Injection Molding), su cuerpo único tiene alta precisión y resistencia.

EXCELENTE ADAPTACIÓN Y RETENCIÓN MECÁNICA

Con geometría optimizada y tratamiento de superficie, es la base más eficaz del mercado de brackets metálicos, pues cuenta con:

- Curvaturas méso-distal y cérvico-oclusal, garantizando una mejor adaptación
- Micro pines que maximizan el área de adhesión y el anclaje del adhesivo
- Bordes de contención que facilitan la remoción del excedente marginal, especialmente importante en los brackets autoligados, para evitar que el exceso de adhesivo fluya hacia arriba de la base y obstruya el Clip
- Chorreado para aumentar la rugosidad

Teste do Clip



Morelli SLI Prescription Roth (Self-Ligating Interactive) with hooks in the canines (9°) and premolars

3.0mm	3.0mm	3.0mm	2.9mm	2.9mm	2.9mm	2.9mm	3.0mm	3.0mm	3.0mm	◀ Length
-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°	◀ Torque
0°	0°	+9°	+9°	+5°	+5°	+9°	+9°	0°	0°	◀ Angle
10.14.007	10.14.007	10.14.005	10.14.003	10.14.001	10.14.002	10.14.004	10.14.006	10.14.008	10.14.008	◀ Replacement .022" (5 UN)
15	14	13	12	11	21	22	23	24	25	
45	44	43	42	41	31	32	33	34	35	
3.0mm	3.0mm	3.0mm	2.7mm	2.7mm	2.7mm	2.7mm	3.0mm	3.0mm	3.0mm	◀ Length
0°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°	◀ Torque
-22°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°	◀ Angle
10.14.014	10.14.012	10.14.010	10.14.009	10.14.009	10.14.009	10.14.009	10.14.011	10.14.013	10.14.015	◀ Replacement .022" (5 UN)

SLI
Self Ligating Interactive



1 CASO • KIT 20 PIEZA
SLOT22 REF 10.14.900

Morelli SLP Prescription Roth (Self-ligating Passive) with hooks in the canines (9°) and premolars

3.0mm	3.0mm	3.0mm	2.9mm	2.9mm	2.9mm	2.9mm	3.0mm	3.0mm	3.0mm	◀ Length
-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°	◀ Torque
0°	0°	+9°	+9°	+5°	+5°	+9°	+9°	0°	0°	◀ Angle
10.13.007	10.13.007	10.13.005	10.13.003	10.13.001	10.13.002	10.13.004	10.13.006	10.13.008	10.13.008	◀ Replacement .022" (5 UN)
15	14	13	12	11	21	22	23	24	25	
45	44	43	42	41	31	32	33	34	35	
3.0mm	3.0mm	3.0mm	2.7mm	2.7mm	2.7mm	2.7mm	3.0mm	3.0mm	3.0mm	◀ Length
-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°	◀ Torque
0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°	◀ Angle
10.13.014	10.13.012	10.13.010	10.13.009	10.13.009	10.13.009	10.13.009	10.13.011	10.13.013	10.13.015	◀ Replacement .022" (5 UN)

SLP
Self Ligating Passive



1 CASO • KIT 20 PIEZA
SLOT22 REF 10.13.900

Morelli SLH Prescription Roth (Self-ligating Hybrid) with hooks in the canines (9°) and premolars

SLP			SLI				SLP			◀ Length
3.0mm	3.0mm	3.0mm	2.9mm	2.9mm	2.9mm	2.9mm	3.0mm	3.0mm	3.0mm	◀ Torque
-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°	◀ Angle
0°	0°	+9°	+9°	+5°	+5°	+9°	+9°	0°	0°	◀ Replacement .022" (5 UN)
10.13.007	10.13.007	10.13.005	10.14.003	10.14.001	10.14.002	10.14.004	10.13.006	10.13.008	10.13.008	
15	14	13	12	11	21	22	23	24	25	
45	44	43	42	41	31	32	33	34	35	
SLP			SLI				SLP			◀ Length
3.0mm	3.0mm	3.0mm	2.7mm	2.7mm	2.7mm	2.7mm	3.0mm	3.0mm	3.0mm	◀ Torque
-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°	◀ Angle
0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°	◀ Replacement .022" (5 UN)
10.13.014	10.13.012	10.13.010	10.14.009	10.14.009	10.14.009	10.14.009	10.13.011	10.13.013	10.13.015	

SLH
Self Ligating Hybrid



1 CASO • KIT 20 PIEZA
SLOT22 REF 10.23.900

#TIP

Discover the sequence of arches for self-ligating brackets on the page 70.

SLI
Self Ligating Interactive




1 CASE - 20 PIECE KIT

SLOT22 REF 10.27.900

Morelli SLI Prescription M.B.T. (Interactive Self-ligating) with canine and premolar hooks

Length		3.0mm	3.0mm	3.0mm	2.9mm	2.9mm	2.9mm	2.9mm	3.0mm	3.0mm	3.0mm
Torque		-7°	-7°	0°	+10°	+17°	+17°	+10°	0°	-7°	-7°
Angle		0°	0°	+8°	+8°	+4°	+4°	+8°	+8°	0°	0°
Replacement .022" (5 UN)		10.27.007	10.27.007	10.27.005	10.27.003	10.27.001	10.27.002	10.27.004	10.27.006	10.27.008	10.27.008



Length		3.0mm	3.0mm	3.0mm	2.7mm	2.7mm	2.7mm	2.7mm	3.0mm	3.0mm	3.0mm
Torque		-17°	-12°	0°	-6°	-6°	-6°	-6°	0°	-12°	-17°
Angle		+2°	+2°	+3°	0°	0°	0°	0°	+3°	+2°	+2°
Replacement .022" (5 UN)		10.27.014	10.27.012	10.27.010	10.27.009	10.27.009	10.27.009	10.27.009	10.27.011	10.27.013	10.27.015

SLP
Self Ligating Passive




1 CASE - 20 PIECE KIT

SLOT22 REF 10.28.900

Morelli SLP Prescription M.B.T. (Passive Self-ligating) with canine and premolar hooks

Length		3.0mm	3.0mm	3.0mm	2.9mm	2.9mm	2.9mm	2.9mm	3.0mm	3.0mm	3.0mm
Torque		-7°	-7°	0°	+10°	+17°	+17°	+10°	0°	-7°	-7°
Angle		0°	0°	+8°	+8°	+4°	+4°	+8°	+8°	0°	0°
Replacement .022" (5 UN)		10.28.007	10.28.007	10.28.005	10.28.003	10.28.001	10.28.002	10.28.004	10.28.006	10.28.008	10.28.008



Length		3.0mm	3.0mm	3.0mm	2.7mm	2.7mm	2.7mm	2.7mm	3.0mm	3.0mm	3.0mm
Torque		-17°	-12°	0°	-6°	-6°	-6°	-6°	0°	-12°	-17°
Angle		+2°	+2°	+3°	0°	0°	0°	0°	+3°	+2°	+2°
Replacement .022" (5 UN)		10.28.014	10.28.012	10.28.010	10.28.009	10.28.009	10.28.009	10.28.009	10.28.011	10.28.013	10.28.015

SLH
Self Ligating Hybrid




1 CASE - 20 PIECE KIT

SLOT22 REF 10.23.935

Morelli SLH Prescription M.B.T. (Interactive Self-ligating) with canine and premolar hooks

		SLP			SLI				SLP		
Length		3.0mm	3.0mm	3.0mm	2.9mm	2.9mm	2.9mm	2.9mm	3.0mm	3.0mm	3.0mm
Torque		-7°	-7°	0°	+10°	+17°	+17°	+10°	0°	-7°	-7°
Angle		0°	0°	+8°	+8°	+4°	+4°	+8°	+8°	0°	0°
Replacement .022" (5 UN)		10.28.007	10.28.007	10.28.005	10.27.003	10.27.001	10.27.002	10.27.004	10.28.006	10.28.008	10.28.008



Length		3.0mm	3.0mm	3.0mm	2.7mm	2.7mm	2.7mm	2.7mm	3.0mm	3.0mm	3.0mm
Torque		-17°	-12°	0°	-6°	-6°	-6°	-6°	0°	-12°	-17°
Angle		+2°	+2°	+3°	0°	0°	0°	0°	+3°	+2°	+2°
Replacement .022" (5 UN)		10.28.014	10.28.012	10.28.010	10.27.009	10.27.009	10.27.009	10.27.009	10.28.011	10.28.013	10.28.015

SLP

SLI

SLP

Clip Test



Aesthetic Brackets



poly *Safira*



MORE THAN PURE AESTHETICS

- Round Center Hook
- Maximum Polishing
- Excellent Mechanical Retention
- Polycrystalline Sapphire
- Rounded Contours
- Pure Aesthetics
- Bilateral Positions

EVERYTHING YOU EXPECT FROM AN AESTHETIC BRACKET!

- Polycrystalline, translucent and highly aesthetic
- Comfort of the rounded and polished profile
- High glide at the ends of the slot
- Centralized sphere-shaped hooks
- Greater comfort for the patient and less replacement code
- Application of the same bracket in bilateral positions


Ceramic



COMFORT THAN PURE AESTHETICS

Morelli PolySafira Prescription Roth with hooks in the canines (9°) and premolars

3.2mm	3.2mm	3.2mm	2.9mm	3.6mm	3.6mm	2.9mm	3.2mm	3.2mm	3.2mm	◀ Length
-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°	◀ Torque
0°	0°	+9°	+9°	+5°	+5°	+9°	+9°	0°	0°	◀ Angle
10.24.007	10.24.007	10.24.005	10.24.003	10.24.001	10.24.002	10.24.004	10.24.006	10.24.007	10.24.007	◀ Replacement .022" (2 UN)



3.2mm	3.2mm	3.2mm	2.9mm	2.9mm	2.9mm	2.9mm	3.2mm	3.2mm	3.2mm	◀ Length
-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°	◀ Torque
0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°	◀ Angle
10.24.012	10.24.011	10.24.010	10.24.008	10.24.008	10.24.008	10.24.008	10.24.009	10.24.011	10.24.012	◀ Replacement .022" (2 UN)

Safira



1 CASE - 20 PIECE KIT
SLOT22 REF 10.24.900

Removal tip for ceramic brackets.



Ceramic




1 CASE - 20 PIECE KIT
SLOT22 REF 10.11.900

Check how the strength test is performed on Morelli ceramic pieces.



Morelli Ceramic Prescription Roth with hooks in the canines (9°) and premolars


3.2mm	3.2mm	3.2mm	2.9mm	3.6mm	3.6mm	2.9mm	3.2mm	3.2mm	3.2mm	◀ Length
-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°	◀ Torque
0°	0°	+9°	+9°	+5°	+5°	+9°	+9°	0°	0°	◀ Angle
10.11.007	10.11.007	10.11.005	10.11.003	10.11.001	10.11.002	10.11.004	10.11.006	10.11.007	10.11.007	◀ Replacement .022" (2 UN)



3.2mm	3.2mm	3.2mm	2.9mm	2.9mm	2.9mm	2.9mm	3.2mm	3.2mm	3.2mm	◀ Length
-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°	◀ Torque
0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°	◀ Angle
10.11.012	10.11.011	10.11.009	10.11.008	10.11.008	10.11.008	10.11.008	10.11.010	10.11.011	10.11.012	◀ Replacement .022" (2 UN)

Morelli Composite Prescription Roth with hooks in the canines (13°) and premolars

3.2mm	3.2mm	3.2mm	2.9mm	3.9mm	3.9mm	2.9mm	3.2mm	3.2mm	3.2mm	◀ Length
-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°	◀ Torque
0°	0°	+13°	+9°	+5°	+5°	+9°	+13°	0°	0°	◀ Angle
10.17.007	10.17.007	10.17.005	10.17.003	10.17.001	10.17.002	10.17.004	10.17.006	10.17.008	10.17.008	◀ Replacement .022" (10 UN)



3.2mm	3.2mm	3.2mm	2.9mm	2.9mm	2.9mm	2.9mm	3.2mm	3.2mm	3.2mm	◀ Length
-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°	◀ Torque
0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°	◀ Angle
10.17.014	10.17.012	10.17.010	10.17.009	10.17.009	10.17.009	10.17.009	10.17.011	10.17.013	10.17.015	◀ Replacement .022" (10 UN)

Composite



1 CASE - 20 PIECE KIT
SLOT22 REF 10.17.901

10 CASES - 200 PIECE KIT
SLOT22 REF 10.17.932

PRESCRIPTION M.B.T

poly *Safira*




1 CASE - 20 PIECE KIT

SLOT22 REF 10.33.900

Morelli PolySafira Prescription M.B.T. with hooks in the canines (9°) and premolars

Length	3.2mm	3.2mm	3.2mm	2.9mm	3.6mm	3.6mm	2.9mm	3.2mm	3.2mm	3.2mm
Torque	-7°	-7°	0°	+10°	+17°	+17°	+10°	0°	-7°	-7°
Angle	0°	0°	+8°	+8°	+4°	+4°	+8°	+8°	0°	0°
Replacement .022" (2 UN)	10.33.007	10.33.007	10.33.005	10.33.003	10.33.001	10.33.002	10.33.004	10.33.006	10.33.007	10.33.007



Length	3.2mm	3.2mm	3.2mm	2.9mm	2.9mm	2.9mm	2.9mm	3.2mm	3.2mm	3.2mm
Torque	-17°	-12°	0°	-6°	-6°	-6°	-6°	0°	-12°	-17°
Angle	+2°	+2°	+3°	0°	0°	0°	0°	+3°	+2°	+2°
Replacement .022" (2 UN)	10.33.013	10.33.011	10.33.009	10.33.008	10.33.008	10.33.008	10.33.008	10.33.010	10.33.012	10.33.014

Ceramic




1 CASE - 20 PIECE KIT

SLOT22 REF 10.32.900

Morelli Ceramic Prescription M.B.T. with hooks in the canines (9°) and premolars

Length	3.2mm	3.2mm	3.2mm	2.9mm	3.6mm	3.6mm	2.9mm	3.2mm	3.2mm	3.2mm
Torque	-7°	-7°	0°	+10°	+17°	+17°	+10°	0°	-7°	-7°
Angle	0°	0°	+8°	+8°	+4°	+4°	+8°	+8°	0°	0°
Replacement .022" (5 UN)	10.32.007	10.32.007	10.32.005	10.32.003	10.32.001	10.32.002	10.32.004	10.32.006	10.32.007	10.32.007



Length	3.2mm	3.2mm	3.2mm	2.9mm	2.9mm	2.9mm	2.9mm	3.2mm	3.2mm	3.2mm
Torque	-17°	-12°	0°	-6°	-6°	-6°	-6°	0°	-12°	-17°
Angle	+2°	+2°	+3°	0°	0°	0°	0°	+3°	+2°	+2°
Replacement .022" (5 UN)	10.32.013	10.32.011	10.32.009	10.32.008	10.32.008	10.32.008	10.32.008	10.32.010	10.32.012	10.32.014

PRESCRIPTION EDGEWISE

Ceramic




1 CASO - KIT 20 PIEZAS

REF 10.18.901

#TIP

Morelli Ceramic Prescription Edgewise with hooks in the canines and premolars

Length	3.2mm	3.2mm	3.2mm	2.9mm	3.6mm	3.6mm	2.9mm	3.2mm	3.2mm	3.2mm
Torque	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
Angle	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
Replacement .022" (10 UN)	10.18.013	10.18.013	10.18.013	10.18.012	10.18.011	10.18.011	10.18.012	10.18.013	10.18.013	10.18.013



Length	3.2mm	3.2mm	3.2mm	2.9mm	2.9mm	2.9mm	2.9mm	3.2mm	3.2mm	3.2mm
Torque	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
Angle	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
Replacement .022" (10 UN)	10.18.013	10.18.013	10.18.013	10.18.012	10.18.012	10.18.012	10.18.012	10.18.013	10.18.013	10.18.013

Metal Brackets





super kit
BRACKETS + TUBES

2000
parts
+
400
parts

BRACKETS *Max*
Prescription Roth
Slot 22

TUBES FOR GLUING
All Quadrants Roth
Prescription - Simple

- Composed of Max Line Brackets and Bonding Tubes, all produced using modern MIM (Metal Injection Molding) technology, which results in single-body parts.
- Rounded finish that provides more comfort to the patient;
- Highly adherent base that prevents detachment.
- Indicated for all types of dental malocclusions, especially the most severe ones;
- The best cost-benefit among Morelli bracket options;
- Allows the patient to choose different bandage colors for each maintenance!

10.15.413 Bráquete Roth Max .022"  14	10.15.409 Bráquete Roth Max .022"  13	10.15.402 Bráquete Roth Max .022"  12	10.15.400 Bráquete Roth Max .022"  11	10.15.401 Bráquete Roth Max .022"  21	10.15.403 Bráquete Roth Max .022"  22	10.15.411 Bráquete Roth Max .022"  23	10.15.415 Bráquete Roth Max .022"  24
10.15.413 Bráquete Roth Max .022"  14	20.11.211 Tubo Roth Simples .022"  16 17	SUPER KIT brackets + tubes Part code reference				20.11.212 Tubo Roth Simples .022"  27 26	10.15.415 Bráquete Roth Max .022"  24
10.15.426 Bráquete Roth Max .022"  46	20.11.221 Tubo Roth Simples .022"  47					20.11.222 Tubo Roth Simples .022"  37 36	10.15.428 Bráquete Roth Max .022"  34
10.15.422 Bráquete Roth Max .022"  44	10.15.418 Bráquete Roth Max .022"  43	10.15.416 Bráquete Roth Max .022"  42	10.15.416 Bráquete Roth Max .022"  41	10.15.416 Bráquete Roth Max .022"  31	10.15.416 Bráquete Roth Max .022"  32	10.15.420 Bráquete Roth Max .022"  33	10.15.424 Bráquete Roth Max .022"  34

Quantity of pieces per cavity in this package

Kit de Brackets 100 casos + 400 Tubos Prescripción Roth - Max .022" - Can. 9° Ang. - Gancho Can./Prés

2.9mm	2.9mm	2.9mm	2.9mm	2.9mm	2.5mm	3.2mm	3.2mm	2.5mm	2.9mm	2.9mm	2.9mm	2.9mm	2.9mm	◀ Length
-14°	-14°	-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°	-14°	-14°	◀ Torque
0°	0°	0°	0°	+9°	+9°	+5°	+5°	+9°	+9°	0°	0°	0°	0°	◀ Angle
14°	14°											14°	14°	◀ Distal
20.11.211	20.11.211	10.15.413	10.15.413	10.15.409	10.15.402	10.15.400	10.15.401	10.15.403	10.15.411	10.15.415	10.15.415	20.11.212	20.11.212	◀ Replacement .022" (10 UN)

17	16	15	14	13	12	11	21	22	23	24	25	26	27	
														100 CASES - 2400 PIECES KIT
														SLOT22 REF 10.15.723
47	46	45	44	43	42	41	31	32	33	34	35	36	37	

2.9mm	2.9mm	2.9mm	2.9mm	2.9mm	2.3mm	2.3mm	2.3mm	2.3mm	2.9mm	2.9mm	2.9mm	2.9mm	2.9mm	◀ Length
-30°	-30°	-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°	-30°	-30°	◀ Torque
0°	0°	0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°	0°	0°	◀ Angle
4°	4°											4°	4°	◀ Distal
20.11.221	20.11.221	10.15.426	10.15.422	10.15.418	10.15.416	10.15.416	10.15.416	10.15.416	10.15.420	10.15.424	10.15.428	20.11.222	20.11.222	◀ Replacement .022" (10 UN)

Kit de Brackets 100 casos + 400 Tubos Prescripción Roth - Max .022" - Can. 13° Ang. - Gancho Can./Prés

2.9mm	2.9mm	2.9mm	2.9mm	2.9mm	2.5mm	3.2mm	3.2mm	2.5mm	2.9mm	2.9mm	2.9mm	2.9mm	2.9mm	◀ Length
-14°	-14°	-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°	-14°	-14°	◀ Torque
0°	0°	0°	0°	+13°	+9°	+5°	+5°	+9°	+13°	0°	0°	0°	0°	◀ Angle
14°	14°											14°	14°	◀ Distal
20.11.211	20.11.211	10.15.413	10.15.413	10.15.405	10.115.402	10.15.400	10.15.401	10.15.403	10.15.407	10.15.415	10.15.415	20.11.212	20.11.212	◀ Replacement .022" (10 UN)

17	16	15	14	13	12	11	21	22	23	24	25	26	27	
														100 CASES - 2400 PIECES KIT
														SLOT22 REF 10.15.725
47	46	45	44	43	42	41	31	32	33	34	35	36	37	

2.9mm	2.9mm	2.9mm	2.9mm	2.9mm	2.3mm	2.3mm	2.3mm	2.3mm	2.9mm	2.9mm	2.9mm	2.9mm	2.9mm	◀ Length
-30°	-30°	-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°	-30°	-30°	◀ Torque
0°	0°	0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°	0°	0°	◀ Angle
4°	4°											4°	4°	◀ Distal
20.11.221	20.11.221	10.15.426	10.15.422	10.15.418	10.15.416	10.15.416	10.15.416	10.15.416	10.15.420	10.15.424	10.15.428	20.11.222	20.11.222	◀ Replacement .022" (10 UN)

LINE BRACKETS

Max

The Max Line brackets are constructed in a single body using the MIM (Metal Injection Molding) system, without welding, and have the strength of stainless steel.

- Patient comfort
- High performance in bonding



OPENING

Allows excess resin to escape and facilitates photopolymerization.

SANDBLA

Increased retention by 35%.

CERVICO-OCCLUSIONAL AND MESIODISTAL CURVATURE

The bases have been improved to obtain maximum adaptation to the enamel surface and, consequently, better bonding

MICRO RETENTIVE PINS

16% more bonding efficiency.



THE MAX LINE IS THE MOST COMPLETE IN PRESCRIPTIONS

The Max line offers the main prescriptions, and all are produced using the MIM (Metal Injection Molding) process.

When it comes to brackets, regardless of the production process, Morelli offers full availability and the largest diversity of prescriptions in Latin America.



ROTH

M.B.T.

CAPELOZZA

EDGEWISE

Morelli Max Prescription Roth with hooks in the canines (13°) and premolars

2.9mm	2.9mm	2.7mm	2.5mm	3.2mm	3.2mm	2.5mm	2.7mm	2.9mm	2.9mm	◀ Length
-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°	◀ Torque
0°	0°	+13°	+9°	+5°	+5°	+9°	+13°	0°	0°	◀ Angle
10.15.413	10.15.413	10.15.405	10.15.402	10.15.400	10.15.401	10.15.403	10.15.407	10.15.415	10.15.415	◀ Replacement .022" (10 UN)
10.15.313	10.15.313	10.15.305	10.15.302	10.15.300	10.15.301	10.15.303	10.15.307	10.15.315	10.15.315	◀ Replacement .018" (10 UN)
15	14	13	12	11	21	22	23	24	25	
45	44	43	42	41	31	32	33	34	35	
2.9mm	2.9mm	2.7mm	2.2mm	2.2mm	2.2mm	2.2mm	2.7mm	2.9mm	2.9mm	◀ Length
-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°	◀ Torque
0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°	◀ Angle
10.15.426	10.15.422	10.15.418	10.15.416	10.15.416	10.15.416	10.15.416	10.15.420	10.15.424	10.15.428	◀ Replacement .022" (10 UN)
10.15.326	10.15.322	10.15.318	10.15.316	10.15.316	10.15.316	10.15.316	10.15.320	10.15.324	10.15.328	◀ Replacement .018" (10 UN)

Morelli Max Prescription Roth with hooks in the canines (9°) and premolars

2.9mm	2.9mm	2.7mm	2.5mm	3.2mm	3.2mm	2.5mm	2.7mm	2.9mm	2.9mm	◀ Length
-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°	◀ Torque
0°	0°	+9°	+9°	+5°	+5°	+9°	+9°	0°	0°	◀ Angle
10.15.413	10.15.413	10.15.409	10.15.402	10.15.400	10.15.401	10.15.403	10.15.411	10.15.415	10.15.415	◀ Replacement .022" (10 UN)
10.15.313	10.15.313	10.15.309	10.15.302	10.15.300	10.15.301	10.15.303	10.15.311	10.15.315	10.15.315	◀ Replacement .018" (10 UN)
15	14	13	12	11	21	22	23	24	25	
45	44	43	42	41	31	32	33	34	35	
2.9mm	2.9mm	2.7mm	2.2mm	2.2mm	2.2mm	2.2mm	2.7mm	2.9mm	2.9mm	◀ Length
-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°	◀ Torque
0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°	◀ Angle
10.15.426	10.15.422	10.15.418	10.15.416	10.15.416	10.15.416	10.15.416	10.15.420	10.15.424	10.15.428	◀ Replacement .022" (10 UN)
10.15.326	10.15.322	10.15.318	10.15.316	10.15.316	10.15.316	10.15.316	10.15.320	10.15.324	10.15.328	◀ Replacement .018" (10 UN)

Morelli Max Two-dimensional Prescription Roth with hooks in the canines (9°) and premolars

2.9mm	2.9mm	2.9mm	2.5mm	3.2mm	3.2mm	2.5mm	2.9mm	2.9mm	2.9mm	◀ Length
-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°	◀ Torque
0°	0°	+9°	+9°	+5°	+5°	+9°	+9°	0°	0°	◀ Angle
10.15.413	10.15.413	10.15.409	10.15.302	10.15.300	10.15.301	10.15.303	10.15.411	10.15.415	10.15.415	◀ Replacement .022" (10 UN)
15	14	13	12	11	21	22	23	24	25	
45	44	43	42	41	31	32	33	34	35	
2.9mm	2.9mm	2.9mm	2.3mm	2.3mm	2.3mm	2.3mm	2.9mm	2.9mm	2.9mm	◀ Length
-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°	◀ Torque
0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°	◀ Angle
10.15.426	10.15.422	10.15.418	10.15.316	10.15.316	10.15.316	10.15.316	10.15.420	10.15.424	10.15.428	◀ Replacement .022" (10 UN)

Max



1 CASE - 20 PIECE KIT

SLOT22 REF 10.15.905

SLOT18 REF 10.15.904

10 CASES - 200 PIECE KIT

SLOT22 REF 10.15.935

SLOT18 REF 10.15.934

100 CASES - 2000 PIECE KIT

SLOT22 REF 10.15.925

SLOT18 REF 10.15.924

Max



1 CASE - 20 PIECE KIT

SLOT22 REF 10.15.903

SLOT18 REF 10.15.902

10 CASES - 200 PIECE KIT

SLOT22 REF 10.15.933

SLOT18 REF 10.15.932

100 CASES - 2000 PIECE KIT

SLOT22 REF 10.15.923

SLOT18 REF 10.15.922

Max
Bidimensional



1 CASE - 20 PIECE KIT

SLOT18 REF 10.15.930

SLOT22

Max



1 CASE - 20 PIECE KIT

SLOT22 REF 10.19.904

10 CASES - 200 PIECE KIT

SLOT22 REF 10.19.934

Morelli Max Prescription Capellozza with hooks in the canines and premolars

Length		2.9mm	2.9mm	2.9mm	2.5mm	3.2mm	3.2mm	2.5mm	2.9mm	2.9mm	2.9mm
Torque		-7°	-7°	-5°	+3°	+7°	+7°	+3°	-5°	-7°	-7°
Angle		0°	0°	+8°	+9°	+5°	+5°	+9°	+8°	0°	0°
Replacement .022" (10 UN)		10.19.107	10.19.107	10.19.105	10.19.103	10.19.101	10.19.102	10.19.104	10.19.106	10.19.108	10.19.108
		15	14	13	12	11	21	22	23	24	25
		45	44	43	42	41	31	32	33	34	35

Length		2.9mm	2.9mm	2.9mm	2.3mm	2.3mm	2.3mm	2.3mm	2.9mm	2.9mm	2.9mm
Torque		-22°	-17°	-11°	-1°	-1°	-1°	-1°	-11°	-17°	-22°
Angle		+2°	+2°	+5°	+2°	+2°	+2°	+2°	+5°	+2°	+2°
Replacement .022" (10 UN)		10.19.115	10.19.113	10.19.111	10.19.109	10.19.109	10.19.110	10.19.110	10.19.112	10.19.114	10.19.116

Max



1 CASE - 20 PIECE KIT

SLOT22 REF 10.19.905

10 CASES - 200 PIECE KIT

SLOT22 REF 10.19.935

Morelli Max Prescription Capellozza Standard II with hooks in the canines and premolars

Length		2.9mm	2.9mm	2.9mm	2.5mm	3.2mm	3.2mm	2.5mm	2.9mm	2.9mm	2.9mm
Torque		-7°	-7°	-5°	+3°	+7°	+7°	+3°	-5°	-7°	-7°
Angle		0°	0°	+5°	+9°	+5°	+5°	+9°	+5°	0°	0°
Replacement .022" (10 UN)		10.19.107	10.19.107	10.19.117	10.19.103	10.19.101	10.19.102	10.19.104	10.19.118	10.19.108	10.19.108
		15	14	13	12	11	21	22	23	24	25
		45	44	43	42	41	31	32	33	34	35

Length		2.9mm	2.9mm	2.9mm	2.3mm	2.3mm	2.3mm	2.3mm	2.9mm	2.9mm	2.9mm
Torque		-22°	-17°	-11°	+4°	+4°	+4°	+4°	-11°	-17°	-22°
Angle		+2°	+2°	+5°	+2°	+2°	+2°	+2°	+5°	+2°	+2°
Replacement .022" (10 UN)		10.19.115	10.19.113	10.19.111	10.19.119	10.19.119	10.19.120	10.19.120	10.19.112	10.19.114	10.19.116

Max



1 CASE - 20 PIECE KIT

SLOT22 REF 10.19.906

10 CASES - 200 PIECE KIT

SLOT22 REF 10.19.936


Morelli Max Prescription Capellozza Standard III with hooks in the canines and premolars

Length		2.9mm	2.9mm	2.9mm	2.5mm	3.2mm	3.2mm	2.5mm	2.9mm	2.9mm	2.9mm
Torque		-7°	-7°	-2°	+10°	+14°	+14°	+10°	-2°	-7°	-7°
Angle		0°	0°	+11°	+9°	+5°	+5°	+9°	+11°	0°	0°
Replacement .022" (10 UN)		10.19.107	10.19.107	10.19.125	10.19.123	10.19.121	10.19.122	10.19.124	10.19.126	10.19.108	10.19.108
		15	14	13	12	11	21	22	23	24	25
		45	44	43	42	41	31	32	33	34	35

Length		2.9mm	2.9mm	2.9mm	2.3mm	2.3mm	2.3mm	2.3mm	2.7mm	2.9mm	2.9mm
Torque		-22°	-17°	-11°	-6°	-6°	-6°	-6°	-11°	-17°	-22°
Angle		+2°	+2°	0°	0°	0°	0°	0°	0°	+2°	+2°
Replacement .022" (10 UN)		10.19.115	10.19.113	10.19.128	10.19.127	10.19.127	10.19.127	10.19.127	10.19.129	10.19.114	10.19.116

Morelli Max Prescripción MBT con ganchos en premolares

2.9mm	2.9mm	2.9mm	2.5mm	3.2mm	3.2mm	2.5mm	2.9mm	2.9mm	2.9mm	◀ Longitud
-7°	-7°	-7°	+10°	+17°	+17°	+10°	-7°	-7°	-7°	◀ Torque
0°	0°	+8°	+8°	+4°	+4°	+8°	+8°	0°	0°	◀ Ángulo
10.35.109	10.35.109	10.35.105	10.35.103	10.35.101	10.35.102	10.35.104	10.35.106	10.35.110	10.35.110	◀ Repuesto .022" (10 UN)



2.9mm	2.9mm	2.9mm	2.3mm	2.3mm	2.3mm	2.3mm	2.9mm	2.9mm	2.9mm	◀ Longitud
-17°	-12°	-6°	-6°	-6°	-6°	-6°	-6°	-12°	-17°	◀ Torque
+2°	+2°	+3°	0°	0°	0°	0°	+3°	+2°	+2°	◀ Ángulo
10.35.118	10.35.116	10.35.112	10.35.111	10.35.111	10.35.111	10.35.111	10.35.113	10.35.117	10.35.119	◀ Repuesto .022" (10 UN)



1 CASO • KIT 20 PIEZA


SLOT22 REF 10.35.906

10 CASOS • KIT 200 PIEZAS

SLOT22 REF 10.35.936

Morelli Max Prescripción M.B.T. with hooks in the canines and premolars

2.9mm	2.9mm	2.9mm	2.5mm	3.2mm	3.2mm	2.5mm	2.9mm	2.9mm	2.9mm	◀ Length
-7°	-7°	0°	+10°	+17°	+17°	+10°	0°	-7°	-7°	◀ Torque
0°	0°	+8°	+8°	+4°	+4°	+8°	+8°	0°	0°	◀ Angle
10.35.109	10.35.109	10.35.107	10.35.103	10.35.101	10.35.102	10.35.104	10.35.108	10.35.110	10.35.110	◀ Replacement .022" (10 UN)



2.9mm	2.9mm	2.9mm	2.9mm	2.3mm	2.3mm	2.3mm	2.9mm	2.9mm	2.9mm	◀ Length
-17°	-12°	0°	-6°	-6°	-6°	-6°	0°	-12°	-17°	◀ Torque
2°	+2°	+3°	0°	0°	0°	0°	+3°	+2°	+2°	◀ Angle
10.35.118	10.35.116	10.35.114	10.35.111	10.35.111	10.35.111	10.35.111	10.35.115	10.35.117	10.35.119	◀ Replacement .022" (10 UN)



1 CASE - 20 PIECE KIT

SLOT22 REF 10.35.905

10 CASES - 200 PIECE KIT


SLOT22 REF 10.35.935

PRESCRIPTION **EDGEWISE**

Morelli Slim Max Prescripción Edgewise

2.5mm	2.5mm	2.5mm	2.2mm	3.0mm	3.0mm	2.2mm	2.5mm	2.5mm	2.5mm	◀ Length
0°	0°	0°	0°	0°	0°	0°	0°	0°	0°	◀ Torque
0°	0°	0°	0°	0°	0°	0°	0°	0°	0°	◀ Angle
10.65.403	10.65.403	10.65.403	10.65.402	10.65.401	10.65.401	10.65.402	10.65.403	10.65.403	10.65.403	◀ Replacement .022" (10 UN)

10.65.303 10.65.303 10.65.303 10.65.302 10.65.301 10.65.301 10.65.302 10.65.303 10.65.303 10.65.303 Replacement .018" (10 UN)



2.5mm	2.5mm	2.5mm	2.2mm	2.2mm	2.2mm	2.2mm	2.5mm	2.5mm	2.5mm	◀ Length
0°	0°	0°	0°	0°	0°	0°	0°	0°	0°	◀ Torque
0°	0°	0°	0°	0°	0°	0°	0°	0°	0°	◀ Angle
10.65.403	10.65.403	10.65.403	10.65.402	10.65.402	10.65.402	10.65.402	10.65.403	10.65.403	10.65.403	◀ Replacement .022" (10 UN)

10.65.303 10.65.303 10.65.303 10.65.302 10.65.302 10.65.302 10.65.302 10.65.303 10.65.303 10.65.303 Replacement .018" (10 UN)

Slim Max



1 CASE - 20 PIECE KIT

SLOT22 REF 10.65.905

SLOT18 REF 10.65.904

10 CASES - 200 PIECE KIT

SLOT22 REF 10.65.935

SLOT18 REF 10.65.934

Slim Max



Morelli Slim Max Prescription Edgewise with hooks in the canines

Length	2.5mm	2.5mm	2.5mm	2.2mm	3.0mm	3.0mm	2.2mm	2.5mm	2.5mm	2.5mm
Torque	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
Angle	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
Replacement .022" (10 UN)	10.65.403	10.65.403	10.65.404	10.65.402	10.65.401	10.65.401	10.65.402	10.65.405	10.65.403	10.65.403
Replacement .018" (10 UN)	10.65.303	10.65.303	10.65.304	10.65.302	10.65.301	10.65.301	10.65.302	10.65.305	10.65.303	10.65.303

15	14	13	12	11	21	22	23	24	25
45	44	43	42	41	31	32	33	34	35

1 CASE - 20 PIECE KIT

SLOT22 REF 10.65.907

SLOT18 REF 10.65.906

10 CASES - 200 PIECE KIT

SLOT22 REF 10.65.937

SLOT18 REF 10.65.936

Length	2.5mm	2.5mm	2.5mm	2.2mm	2.2mm	2.2mm	2.2mm	2.5mm	2.5mm	2.5mm
Torque	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
Angle	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
Replacement .022" (10 UN)	10.65.403	10.65.403	10.65.405	10.65.402	10.65.402	10.65.402	10.65.402	10.65.404	10.65.403	10.65.403
Replacement .018" (10 UN)	10.65.303	10.65.303	10.65.305	10.65.302	10.65.302	10.65.302	10.65.302	10.65.304	10.65.303	10.65.303

PRESCRIPTION ROTH

nano



Morelli Nano Prescription Roth with hooks in the canines (9°) and premolars

Length	2.5mm	2.5mm	2.6mm	2.5mm	3.2mm	3.2mm	2.5mm	2.6mm	2.5mm	2.5mm
Torque	-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°
Angle	0°	0°	+9°	+9°	+5°	+5°	+9°	+9°	0°	0°
Replacement .022" (10 UN)	10.15.509	10.15.509	10.15.505	10.15.503	10.15.501	10.15.502	10.15.504	10.15.506	10.15.510	10.15.510

15	14	13	12	11	21	22	23	24	25
45	44	43	42	41	31	32	33	34	35

1 CASE - 20 PIECE KIT

SLOT22 REF 10.15.906

Length	2.5mm	2.5mm	2.6mm	2.2mm	2.2mm	2.2mm	2.2mm	2.6mm	2.5mm	2.5mm
Torque	-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°
Angle	0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°
Replacement .022" (10 UN)	10.15.516	10.15.514	10.15.512	10.15.511	10.15.511	10.15.511	10.15.511	10.15.513	10.15.515	10.15.517

Monobloc

Nickel Free



Morelli Monobloc (Nickel Free) Prescription Roth with hooks in the canines (13°) and premolars

Length	2.5mm	2.5mm	2.7mm	2.5mm	3.2mm	3.2mm	2.5mm	2.7mm	2.5mm	2.5mm
Torque	-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°
Angle	0°	0°	+13°	+9°	+5°	+5°	+9°	+13°	0°	0°
Replacement .022" (10 UN)	10.15.206	10.15.206	10.15.204	10.15.202	10.15.200	10.15.201	10.15.203	10.15.205	10.15.207	10.15.207

15	14	13	12	11	21	22	23	24	25
45	44	43	42	41	31	32	33	34	35

1 CASE - 20 PIECE KIT

SLOT22 REF 10.15.901

Length	2.5mm	2.5mm	2.6mm	2.3mm	2.3mm	2.3mm	2.3mm	2.6mm	2.5mm	2.5mm
Torque	-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°
Angle	0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°
Replacement .022" (10 UN)	10.15.214	10.15.212	10.15.210	10.15.208	10.15.208	10.15.208	10.15.208	10.15.211	10.15.213	10.15.215

Morelli Light Prescription Roth with canines (9°)

2.5mm	2.5mm	2.7mm	2.5mm	3.2mm	3.2mm	2.5mm	2.7mm	2.5mm	2.5mm	◀ Length
-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°	◀ Torque
0°	0°	+9°	+9°	+5°	+5°	+9°	+9°	0°	0°	◀ Angle
10.10.408	10.10.408	10.10.406	10.10.402	10.10.400	10.10.401	10.10.403	10.10.407	10.10.408	10.10.408	◀ Replacement .022" (10 UN)

15	14	13	12	11	21	22	23	24	25
45	44	43	42	41	31	32	33	34	35



1 CASE - 20 PIECE KIT
SLOT22 **REF 10.10.915**

2.7mm	2.7mm	2.7mm	2.2mm	2.2mm	2.2mm	2.2mm	2.7mm	2.7mm	2.7mm	◀ Length
-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°	◀ Torque
0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°	◀ Angle
10.10.416	10.10.414	10.10.412	10.10.410	10.10.410	10.10.410	10.10.410	10.10.413	10.10.414	10.10.416	◀ Replacement .022" (10 UN)

Morelli Light Prescription Roth with canines (13°)

2.5mm	2.5mm	2.7mm	2.5mm	3.2mm	3.2mm	2.5mm	2.7mm	2.5mm	2.5mm	◀ Length
-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°	◀ Torque
0°	0°	+13°	+9°	+5°	+5°	+9°	+13°	0°	0°	◀ Angle
10.10.408	10.10.408	10.10.404	10.10.402	10.10.400	10.10.401	10.10.403	10.10.405	10.10.408	10.10.408	◀ Replacement .022" (10 UN)
10.10.308	10.10.308	10.10.304	10.10.302	10.10.300	10.10.301	10.10.303	10.10.305	10.10.308	10.10.308	◀ Replacement .018" (10 UN)

15	14	13	12	11	21	22	23	24	25
45	44	43	42	41	31	32	33	34	35



1 CASE - 20 PIECE KIT
SLOT22 **REF 10.10.907**
SLOT18 **REF 10.10.906**

2.7mm	2.7mm	2.7mm	2.2mm	2.3mm	2.3mm	2.2mm	2.7mm	2.7mm	2.7mm	◀ Length
-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°	◀ Torque
0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°	◀ Angle
10.10.416	10.10.414	10.10.412	10.10.410	10.10.410	10.10.410	10.10.410	10.10.413	10.10.414	10.10.416	◀ Replacement .022" (10 UN)
10.10.316	10.10.314	10.10.312	10.10.310	10.10.310	10.10.310	10.10.310	10.10.313	10.10.314	10.10.316	◀ Replacement .018" (10 UN)

Morelli Light Prescription Roth with hooks in the canines (9°)

2.5mm	2.5mm	2.7mm	2.5mm	3.2mm	3.2mm	2.5mm	2.7mm	2.5mm	2.5mm	◀ Length
-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°	◀ Torque
0°	0°	+9°	+9°	+5°	+5°	+9°	+9°	0°	0°	◀ Angle
10.10.408	10.10.408	10.10.422	10.10.402	10.10.400	10.10.401	10.10.403	10.10.423	10.10.408	10.10.408	◀ Replacement .022" (10 UN)

15	14	13	12	11	21	22	23	24	25
45	44	43	42	41	31	32	33	34	35



1 CASE - 20 PIECE KIT
SLOT22 **REF 10.10.963**

2.9mm	2.7mm	2.7mm	2.2mm	2.2mm	2.2mm	2.2mm	2.7mm	2.7mm	2.9mm	◀ Length
-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°	◀ Torque
0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°	◀ Angle
10.10.416	10.10.414	10.10.426	10.10.410	10.10.410	10.10.410	10.10.410	10.10.427	10.10.414	10.10.416	◀ Replacement .022" (10 UN)

Light




1 CASE - 20 PIECE KIT

SLOT22 REF 10.10.961

SLOT18 REF 10.10.960

Morelli Light Prescription Roth with hooks in the canines (13°)

Length	2.5mm	2.5mm	2.7mm	2.5mm	3.2mm	3.2mm	2.5mm	2.7mm	2.5mm	2.5mm
Torque	-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°
Angle	0°	0°	+13°	+9°	+5°	+5°	+9°	+13°	0°	0°
Replacement .022" (10 UN)	10.10.408	10.10.408	10.10.420	10.10.402	10.10.400	10.10.401	10.10.403	10.10.421	10.10.408	10.10.408
Replacement .018" (10 UN)	10.10.308	10.10.308	10.10.320	10.10.302	10.10.300	10.10.301	10.10.303	10.10.321	10.10.308	10.10.308



Length	2.7mm	2.7mm	2.7mm	2.3mm	2.3mm	2.2mm	2.2mm	2.7mm	2.7mm	2.7mm
Torque	-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°
Angle	0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°
Replacement .022" (10 UN)	10.10.416	10.10.414	10.10.426	10.10.410	10.10.410	10.10.410	10.10.410	10.10.427	10.10.414	10.10.416
Replacement .018" (10 UN)	10.10.316	10.10.314	10.10.326	10.10.310	10.10.310	10.10.310	10.10.310	10.10.327	10.10.314	10.10.316

Light



1 CASE - 20 PIECE KIT


SLOT22 REF 10.10.981

10 CASES - 200 PIECE KIT

SLOT22 REF 10.10.997

Morelli Light Prescription Roth with hooks in the canines (9°) and premolars

Length	2.9mm	2.9mm	2.7mm	2.5mm	3.2mm	3.2mm	2.5mm	2.7mm	2.9mm	2.9mm
Torque	-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°
Angle	0°	0°	+9°	+9°	+5°	+5°	+9°	+9°	0°	0°
Replacement .022" (10 UN)	10.10.424	10.10.424	10.10.422	10.10.402	10.10.400	10.10.401	10.10.403	10.10.423	10.10.425	10.10.425



Length	2.9mm	2.9mm	2.7mm	2.2mm	2.2mm	2.2mm	2.2mm	2.7mm	2.9mm	2.9mm
Torque	-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°
Angle	0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°
Replacement .022" (10 UN)	10.10.430	10.10.428	10.10.426	10.10.410	10.10.410	10.10.410	10.10.410	10.10.427	10.10.429	10.10.431

Light



1 CASE - 20 PIECE KIT

SLOT22 REF 10.10.979


SLOT18 REF 10.10.978

10 CASES - 200 PIECE KIT

SLOT22 REF 10.10.996

Morelli Light Prescription Roth with hooks in the canines (13°) and premolars


Length	2.9mm	2.9mm	2.7mm	2.5mm	3.2mm	3.2mm	2.5mm	2.7mm	2.9mm	2.9mm
Torque	-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°
Angle	0°	0°	+13°	+9°	+5°	+5°	+9°	+13°	0°	0°
Replacement .022" (10 UN)	10.10.424	10.10.424	10.10.420	10.10.402	10.10.400	10.10.401	10.10.403	10.10.421	10.10.425	10.10.425
Replacement .018" (10 UN)	10.10.324	10.10.324	10.10.320	10.10.302	10.10.300	10.10.301	10.10.303	10.10.321	10.10.325	10.10.325




Length	2.9mm	2.9mm	2.7mm	2.2mm	2.2mm	2.2mm	2.2mm	2.7mm	2.9mm	2.9mm
Torque	-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°
Angle	0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°
Replacement .022" (10 UN)	10.10.430	10.10.428	10.10.426	10.10.410	10.10.410	10.10.410	10.10.410	10.10.427	10.10.429	10.10.431
Replacement .018" (10 UN)	10.10.330	10.10.328	10.10.326	10.10.310	10.10.310	10.10.310	10.10.310	10.10.327	10.10.329	10.10.331

Morelli Light Prescription Roth with hooks in the canines (13°) and premolars

2.9mm	2.9mm	2.7mm	2.5mm	3.2mm	3.2mm	2.5mm	2.7mm	2.9mm	2.9mm	◀ Length
-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°	◀ Torque
0°	0°	+13°	+9°	+5°	+5°	+9°	+13°	0°	0°	◀ Angle
10.10.424	10.10.424	10.10.420	10.10.436	10.10.434	10.10.435	10.10.437	10.10.421	10.10.425	10.10.425	◀ Replacement .022" (10 UN)



2.9mm	2.9mm	2.7mm	2.2mm	2.2mm	2.2mm	2.2mm	2.7mm	2.9mm	2.9mm	◀ Length
-22°	-17°	-11°	+10°	+10°	+10°	+10°	-11°	-17°	-22°	◀ Torque
0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°	◀ Angle
10.10.430	10.10.428	10.10.426	10.10.432	10.10.438	10.10.438	10.10.433	10.10.427	10.10.429	10.10.431	◀ Replacement .022" (10 UN)



Light




1 CASE - 20 PIECE KIT
SLOT22 REF 10.10.920


PRESCRIPTION **M.B.T.**

Morelli Light HP Two-dimensional Prescription M.B.T. with hooks in the canines

2.9mm	2.9mm	2.9mm	2.5mm	3.2mm	3.2mm	2.5mm	2.9mm	2.9mm	2.9mm	◀ Length
-7°	-7°	0°	+10°	+17	+17°	+10°	0°	-7°	-7°	◀ Torque
0°	0°	+8°	+8°	+4°	+4°	+8°	+8°	0°	0°	◀ Angle
10.35.412	10.35.412	10.35.417	10.35.403	10.35.401	10.35.402	10.35.404	10.35.418	10.35.412	10.35.412	◀ Replacement .022" (10 UN)



2.9mm	2.9mm	2.9mm	2.3mm	2.3mm	2.3mm	2.3mm	2.9mm	2.9mm	2.9mm	◀ Length
-17°	-12°	0°	-6°	-6°	-6°	-6°	0°	-12°	-17°	◀ Torque
0°	+2°	+3°	0°	0°	0°	0°	+3°	+2°	+2°	◀ Angle
10.35.415	10.35.413	10.35.420	10.35.419	10.35.419	10.35.419	10.35.419	10.35.421	10.35.414	10.35.416	◀ Replacement .022" (10 UN)




Light HP




1 CASE - 20 PIECE KIT
SLOT22 REF 10.35.924
SLOT18

Morelli Light Prescription M.B.T.

2.9mm	2.9mm	3.2mm	2.9mm	3.9mm	3.9mm	2.9mm	3.2mm	2.9mm	2.9mm	◀ Length
-7°	-7°	-7°	+10°	+17°	+17°	+10°	-7°	-7°	-7°	◀ Torque
0°	0°	+8°	+8°	+4°	+4°	+8°	+8°	0°	0°	◀ Angle
10.35.007	10.35.007	10.35.030	10.35.103	10.35.101	10.35.102	10.35.104	10.35.031	10.35.007	10.35.007	◀ Replacement .022" (10 UN)



2.9mm	3.2mm	3.2mm	2.2mm	2.2mm	2.2mm	2.2mm	3.2mm	3.2mm	2.9mm	◀ Length
-17°	-12°	-6°	-6°	-6°	-6°	-6°	-6°	-12°	-17°	◀ Torque
+2°	+2°	+3°	0°	0°	0°	0°	+3°	+2°	+2°	◀ Angle
10.35.015	10.35.013	10.35.032	10.35.111	10.35.111	10.35.111	10.35.111	10.35.033	10.35.014	10.35.016	◀ Replacement .022" (10 UN)



Light



1 CASE - 20 PIECE KIT
SLOT22 REF 10.35.904

PRESCRIPTION M.B.T

Light



1 CASE - 20 PIECE KIT

SLOT22 REF 10.35.903

Morelli Light Prescription M.B.T. with hooks in the canines and premolars

Length		2.9mm	2.9mm	2.7mm	2.5mm	3.2mm	3.2mm	2.5mm	2.7mm	2.9mm	2.9mm
Torque		-7°	-7°	0°	+10°	+17°	+17°	+10°	0°	-7°	-7°
Angle		0°	0°	+8°	+8°	+4°	+4°	+8°	+8°	0°	0°
Replacement .022" (10 UN)		10.35.034	10.35.034	10.35.025	10.35.103	10.35.101	10.35.102	10.35.104	10.35.026	10.35.035	10.35.035

15	14	13	12	11	21	22	23	24	25
45	44	43	42	41	31	32	33	34	35

Length		2.9mm	2.9mm	2.9mm	2.3mm	2.3mm	2.3mm	2.3mm	2.9mm	2.9mm	2.9mm
Torque		-17°	-12°	0°	-6°	-6°	-6°	-6°	0°	-12°	-17°
Angle		+2°	+2°	+3°	0°	0°	0°	0°	+3°	+2°	+2°
Replacement .022" (10 UN)		10.35.038	10.35.036	10.35.028	10.35.111	10.35.111	10.35.111	10.35.111	10.35.029	10.35.037	10.35.039

PRESCRIPTION ROTH

STANDARD



1 CASE - 20 PIECE KIT

SLOT22 REF 10.10.911

Morelli Standard Prescription Roth with canines (9°)

Length		2.5mm	2.5mm	3.2mm	2.9mm	3.9mm	3.9mm	2.9mm	3.2mm	2.5mm	2.5mm
Torque		-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°
Angle		0°	0°	+9°	+9°	+5°	+5°	+9°	+9°	0°	0°
Replacement .022" (10 UN)		10.10.208	10.10.208	10.10.206	10.10.202	10.10.200	10.10.201	10.10.203	10.10.207	10.10.208	10.10.208

15	14	13	12	11	21	22	23	24	25
45	44	43	42	41	31	32	33	34	35

Length		2.7mm	2.7mm	3.2mm	2.2mm	2.2mm	2.2mm	2.2mm	3.2mm	2.7mm	2.7mm
Torque		-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°
Angle		0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°
Replacement .022" (10 UN)		10.10.216	10.10.214	10.10.212	10.10.410	10.10.410	10.10.410	10.10.410	10.10.213	10.10.214	10.10.216

STANDARD



1 CASE - 20 PIECE KIT

SLOT22 REF 10.10.901

SLOT18 REF 10.10.900

Morelli Standard Prescription Roth with canines (13°)

Length		2.5mm	2.5mm	3.2mm	2.9mm	3.9mm	3.9mm	2.9mm	3.2mm	2.5mm	2.5mm
Torque		-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°
Angle		0°	0°	+13°	+9°	+5°	+5°	+9°	+13°	0°	0°
Replacement .022" (10 UN)		10.10.208	10.10.208	10.10.204	10.10.202	10.10.200	10.10.201	10.10.203	10.10.205	10.10.208	10.10.208
Replacement .018" (10 UN)		10.10.108	10.10.108	10.10.104	10.10.102	10.10.100	10.10.101	10.10.103	10.10.105	10.10.108	10.10.108

15	14	13	12	11	21	22	23	24	25
45	44	43	42	41	31	32	33	34	35

Length		2.7mm	2.7mm	3.2mm	2.2mm	2.2mm	2.2mm	2.2mm	3.2mm	2.7mm	2.7mm
Torque		-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°
Angle		0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°
Replacement .022" (10 UN)		10.10.216	10.10.214	10.10.212	10.10.410	10.10.410	10.10.410	10.10.410	10.10.213	10.10.214	10.10.216
Replacement .018" (10 UN)		10.10.116	10.10.114	10.10.112	10.10.310	10.10.310	10.10.310	10.10.310	10.10.113	10.10.114	10.10.116

Morelli Standard Prescription Roth with hooks in the canines (9°)

STANDARD

2.5mm	2.5mm	3.2mm	2.9mm	3.9mm	3.9mm	2.9mm	3.2mm	2.5mm	2.5mm	◀ Length
-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°	◀ Torque
0°	0°	+9°	+9°	+5°	+5°	+9°	+9°	0°	0°	◀ Angle
10.10.208	10.10.208	10.10.222	10.10.202	10.10.200	10.10.201	10.10.203	10.10.223	10.10.208	10.10.208	◀ Replacement .022" (10 UN)



1 CASE - 20 PIECE KIT
SLOT22 REF 10.10.953

15	14	13	12	11	21	22	23	24	25
45	44	43	42	41	31	32	33	34	35

2.7mm	2.7mm	3.2mm	2.2mm	2.2mm	2.2mm	2.2mm	3.2mm	2.7mm	2.7mm	◀ Length
-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°	◀ Torque
0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°	◀ Angle
10.10.216	10.10.214	10.10.226	10.10.410	10.10.410	10.10.410	10.10.410	10.10.227	10.10.214	10.10.216	◀ RReplacement .022" (10 UN)

MMorelli Standard Prescription Roth with hooks in the canines (13°)

STANDARD

2.5mm	2.5mm	3.2mm	2.9mm	3.9mm	3.9mm	2.9mm	3.2mm	2.5mm	2.5mm	◀ Length
-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°	◀ Torque
0°	0°	+13°	+9°	+5°	+5°	+9°	+13°	0°	0°	◀ Angle
10.10.208	10.10.208	10.10.220	10.10.202	10.10.200	10.10.201	10.10.203	10.10.221	10.10.208	10.10.208	◀ Replacement .022" (10 UN)
10.10.108	10.10.108	10.10.120	10.10.102	10.10.100	10.10.101	10.10.103	10.10.121	10.10.108	10.10.108	◀ Replacement .018" (10 UN)



1 CASE - 20 PIECE KIT

SLOT22 REF 10.10.951

SLOT18 REF 10.10.950

15	14	13	12	11	21	22	23	24	25
45	44	43	42	41	31	32	33	34	35

2.7mm	2.7mm	3.2mm	2.2mm	2.2mm	2.2mm	2.2mm	3.2mm	2.7mm	2.7mm	◀ Length
-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°	◀ Torque
0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°	◀ c
10.10.216	10.10.214	10.10.226	10.10.410	10.10.410	10.10.410	10.10.410	10.10.227	10.10.214	10.10.216	◀ Replacement .022" (10 UN)
10.10.116	10.10.114	10.10.126	10.10.310	10.10.310	10.10.310	10.10.310	10.10.127	10.10.114	10.10.116	◀ Replacement .018" (10 UN)

Morelli Standard Prescription Roth with hooks in the canines (9°) and premolars

STANDARD

2.9mm	2.9mm	3.2mm	2.9mm	3.9mm	3.9mm	2.9mm	3.2mm	2.9mm	2.9mm	◀ Length
-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°	◀ Torque
0°	0°	+9°	+9°	+5°	+5°	+9°	+9°	0°	0°	◀ Angle
10.10.224	10.10.224	10.10.222	10.10.202	10.10.200	10.10.201	10.10.203	10.10.223	10.10.225	10.10.225	◀ Replacement .022" (10 UN)



1 CASE - 20 PIECE KIT

SLOT22 REF 10.10.971

15	14	13	12	11	21	22	23	24	25
45	44	43	42	41	31	32	33	34	35

2.9mm	2.9mm	3.2mm	2.2mm	2.2mm	2.2mm	2.2mm	3.2mm	2.9mm	2.9mm	◀ Length
-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°	◀ Torque
0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°	◀ Angle
10.10.230	10.10.228	10.10.226	10.10.410	10.10.410	10.10.410	10.10.410	10.10.227	10.10.229	10.10.231	◀ Replacement .022" (10 UN)

STANDARD



1 CASE - 20 PIECE KIT
SLOT22 REF 10.10.969

Morelli Standard Prescription Roth with hooks in the canines (13°) and premolars

Length	3.2mm	3.2mm	3.2mm	2.9mm	3.9mm	3.9mm	2.9mm	3.2mm	3.2mm	3.2mm
Torque	-7°	-7°	-2°	+8°	+12°	+12°	+8°	-2°	-7°	-7°
Angle	0°	0°	+13°	+9°	+5°	+5°	+9°	+13°	0°	0°
Replacement .022" (10 UN)	10.10.224	10.10.224	10.10.220	10.10.202	10.10.200	10.10.201	10.10.203	10.10.221	10.10.225	10.10.225

15	14	13	12	11	21	22	23	24	25
45	44	43	42	41	31	32	33	34	35

Length	3.2mm	3.2mm	3.2mm	2.2mm	2.2mm	2.2mm	2.2mm	3.2mm	3.2mm	3.2mm
Torque	-22°	-17°	-11°	0°	0°	0°	0°	-11°	-17°	-22°
Angle	0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°
Replacement .022" (10 UN)	10.10.230	10.10.228	10.10.226	10.10.410	10.10.410	10.10.410	10.10.410	10.10.227	10.10.229	10.10.231

PRESCRIPTION M.B.T

STANDARD



1 CASE - 20 PIECE KIT
SLOT22 REF 10.35.901

Morelli Standard Prescription M.B.T.

Length	3.2mm	3.2mm	2.9mm	2.9mm	3.9mm	3.9mm	2.9mm	2.9mm	3.2mm	3.2mm
Torque	-7°	-7°	-7°	+10°	+17°	+17°	+10°	-7°	-7°	-7°
Angle	0°	0°	+8°	+8°	+4°	+4°	+8°	+8°	0°	0°
Replacement .022" (10 UN)	10.35.007	10.35.007	10.35.005	10.35.003	10.35.001	10.35.002	10.35.004	10.35.006	10.35.007	10.35.007

15	14	13	12	11	21	22	23	24	25
45	44	43	42	41	31	32	33	34	35

Longitud	2.9mm	2.9mm	2.9mm	2.9mm	2.9mm	2.9mm	2.9mm	2.9mm	2.9mm	2.9mm
Torque	-17°	-12°	-6°	-6°	-6°	-6°	-6°	-6°	-12°	-17°
Ángulo	+2°	+2°	+3°	0°	0°	0°	0°	+3°	+2°	+2°
Replacement .022" (10 UN)	10.35.015	10.35.013	10.35.011	10.35.011	10.35.111	10.35.111	10.35.011	10.35.012	10.35.014	10.35.016

STANDARD



1 CASE - 20 PIECE KIT
SLOT22 REF 10.35.902

Morelli Standard Prescription M.B.T. with hooks in the canines

Length	2.9mm	2.9mm	3.2mm	2.9mm	3.9mm	3.9mm	2.9mm	3.2mm	2.9mm	2.9mm
Torque	-7°	-7°	0°	+10°	+17°	+17°	+10°	0°	-7°	-7°
Angle	0°	0°	+8°	+8°	+4°	+4°	+8°	+8°	0°	0°
Replacement .022" (10 UN)	10.35.007	10.35.007	10.35.017	10.35.003	10.35.001	10.35.002	10.35.004	10.35.018	10.35.007	10.35.007

15	14	13	12	11	21	22	23	24	25
45	44	43	42	41	31	32	33	34	35

Length	2.9mm	2.9mm	3.2mm	2.3mm	2.3mm	2.3mm	2.3mm	3.2mm	2.9mm	2.9mm
Torque	-17°	-12°	0°	-6°	-6°	-6°	-6°	0°	-12°	-17°
Angle	+2°	+2°	+3°	0°	0°	0°	0°	+3°	+2°	+2°
Replacement .022" (10 UN)	10.35.015	10.35.013	10.35.119	10.35.111	10.35.111	10.35.111	10.35.111	10.35.020	10.35.014	10.35.016

Morelli Standard Prescription Edgewise

STANDARD

3.0mm	3.0mm	3.0mm	2.3mm	3.4mm	3.4mm	2.3mm	3.0mm	3.0mm	3.0mm	◀ Length
0°	0°	0°	0°	0°	0°	0°	0°	0°	0°	◀ Torque
0°	0°	0°	0°	0°	0°	0°	0°	0°	0°	◀ Angle
10.30.208	10.30.208	10.30.208	10.30.209	10.30.201	10.30.201	10.30.209	10.30.208	10.30.208	10.30.208	◀ Replacement .022" (10 UN)
10.30.106	10.30.106	10.30.106	10.30.107	10.30.101	10.30.101	10.30.107	10.30.106	10.30.106	10.30.106	◀ Replacement .018" (10 UN)
15	14	13	12	11	21	22	23	24	25	
45	44	43	42	41	31	32	33	34	35	



1 CASE - 20 PIECE KIT

SLOT22 **REF 10.30.901**

SLOT18 **REF 10.30.900**

3.0mm	3.0mm	3.0mm	2.3mm	2.3mm	2.3mm	2.3mm	3.0mm	3.0mm	3.0mm	◀ Longitud
0°	0°	0°	0°	0°	0°	0°	0°	0°	0°	◀ Torque
0°	0°	0°	0°	0°	0°	0°	0°	0°	0°	◀ Ángulo
10.30.208	10.30.208	10.30.208	10.30.209	10.30.209	10.30.209	10.30.209	10.30.208	10.30.208	10.30.208	◀ RReplacement .022" (10 UN)
10.30.106	10.30.106	10.30.106	10.30.107	10.30.107	10.30.107	10.30.107	10.30.106	10.30.106	10.30.106	◀ Replacement .018" (10 UN)

Morelli Slim Prescription Edgewise with hooks in the canines

Slim

2.5mm	2.5mm	2.5mm	2.2mm	3.0mm	3.0mm	2.2mm	2.5mm	2.5mm	2.5mm	◀ Length
0°	0°	0°	0°	0°	0°	0°	0°	0°	0°	◀ Torque
0°	0°	0°	0°	0°	0°	0°	0°	0°	0°	◀ Angle
10.65.206	10.65.206	10.65.204	10.65.203	10.65.201	10.65.201	10.65.203	10.65.205	10.65.206	10.65.206	◀ Replacement .022" (10 UN)
10.65.106	10.65.106	10.65.104	10.65.103	10.65.101	10.65.101	10.65.103	10.65.105	10.65.106	10.65.106	◀ Replacement .018" (10 UN)
15	14	13	12	11	21	22	23	24	25	
45	44	43	42	41	31	32	33	34	35	



1 CASE - 20 PIECE KIT

SLOT22 **REF 10.65.901**

SLOT18 **REF 10.65.900**

2.5mm	2.5mm	2.5mm	2.2mm	2.2mm	2.2mm	2.2mm	2.5mm	2.5mm	2.5mm	◀ Length
0°	0°	0°	0°	0°	0°	0°	0°	0°	0°	◀ Torque
0°	0°	0°	0°	0°	0°	0°	0°	0°	0°	◀ Angle
10.65.206	10.65.206	10.65.205	10.65.203	10.65.203	10.65.203	10.65.203	10.65.204	10.65.206	10.65.206	◀ Replacement .022" (10 UN)
10.65.106	10.65.106	10.65.105	10.65.103	10.65.103	10.65.103	10.65.103	10.65.105	10.65.106	10.65.106	◀ Replacement .018" (10 UN)

Morelli Slim Prescription Edgewise

Slim

2.5mm	2.5mm	2.5mm	2.2mm	3.0mm	3.0mm	2.2mm	2.5mm	2.5mm	2.5mm	◀ Length
0°	0°	0°	0°	0°	0°	0°	0°	0°	0°	◀ Torque
0°	0°	0°	0°	0°	0°	0°	0°	0°	0°	◀ Angle
10.65.206	10.65.206	10.65.206	10.65.203	10.65.201	10.65.201	10.65.203	10.65.206	10.65.206	10.65.206	◀ Replacement .022" (10 UN)
10.65.106	10.65.106	10.65.106	10.65.103	10.65.101	10.65.101	10.65.103	10.65.106	10.65.106	10.65.106	◀ Replacement .018" (10 UN)
15	14	13	12	11	21	22	23	24	25	
45	44	43	42	41	31	32	33	34	35	



1 CASE - 20 PIECE KIT

SLOT22 **REF 10.65.903**

SLOT18 **REF 10.65.902**

2.5mm	2.5mm	2.5mm	2.2mm	2.2mm	2.2mm	2.2mm	2.5mm	2.5mm	2.5mm	◀ Length
0°	0°	0°	0°	0°	0°	0°	0°	0°	0°	◀ Torque
0°	0°	0°	0°	0°	0°	0°	0°	0°	0°	◀ Angle
10.65.206	10.65.206	10.65.206	10.65.203	10.65.203	10.65.203	10.65.203	10.65.206	10.65.206	10.65.206	◀ Replacement .022" (10 UN)
10.65.106	10.65.106	10.65.106	10.65.103	10.65.103	10.65.103	10.65.103	10.65.106	10.65.106	10.65.106	◀ Replacement .018" (10 UN)

Slim




1 CASE ~ 20 PIECE KIT

SLOT18 REF 10.25.905

Morelli Slim Prescription Ricketts with hooks in the canines and premolars

Length	▶	2.5mm	2.5mm	2.5mm	2.2mm	3.0mm	3.0mm	2.2mm	2.5mm	2.5mm	2.5mm
Torque	▶	0°	0°	+7°	+14°	+22°	+22°	+14°	+7°	0°	0°
Angle	▶	0°	0°	+5°	+10°	+3°	+3°	+10°	+5°	0°	0°
Replacement .018" (10 UN)	▶	10.65.104	10.65.104	10.25.305	10.25.303	10.25.301	10.25.302	10.25.304	10.25.306	10.65.105	10.65.105



Length	▶	2.5mm	2.5mm	2.5mm	2.2mm	2.2mm	2.2mm	2.2mm	2.5mm	2.5mm	2.5mm
Torque	▶	0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°
Angle	▶	0°	0°	+5°	0°	0°	0°	0°	+5°	0°	0°
Replacement .018" (10 UN)	▶	10.65.105	10.65.105	10.25.307	10.65.103	10.65.103	10.65.103	10.65.103	10.25.308	10.65.104	10.65.104

Actual




1 CASE ~ 20 PIECE KIT

SLOT18 REF 10.25.901

Morelli Actual Prescription Ricketts with hooks in the canines and premolars

Length	▶	2.5mm	2.5mm	2.7mm	2.5mm	3.2mm	3.2mm	2.5mm	2.7mm	2.5mm	2.5mm
Torque	▶	0°	0°	+7°	+14°	+22°	+22°	+14°	+7°	0°	0°
Angle	▶	0°	0°	+8°	+10°	+3°	+3°	+10°	+8°	0°	0°
Replacement .018" (10 UN)	▶	10.25.109	10.25.107	10.25.105	10.25.103	10.25.101	10.25.102	10.25.104	10.25.106	10.25.108	10.25.110



Length	▶	2.5mm	2.5mm	2.7mm	2.2mm	2.2mm	2.2mm	2.2mm	2.7mm	2.5mm	2.5mm
Torque	▶	0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°
Angle	▶	0°	0°	+2°	0°	0°	0°	0°	+2°	0°	0°
Replacement .018" (10 UN)	▶	10.25.110	10.25.115	10.25.113	10.25.111	10.25.111	10.25.111	10.25.111	10.25.114	10.25.116	10.25.109

Actual Especial




1 CASE ~ 20 PIECE KIT

SLOT18 REF 10.25.900

Morelli Actual Special Prescription Ricketts with hooks in the canines and premolars

Length	▶	2.5mm	2.5mm	2.7mm	2.5mm	3.2mm	3.2mm	2.5mm	2.7mm	2.5mm	2.5mm
Torque	▶	0°	0°	+7°	+10°	+22°	+22°	+10°	+7°	0°	0°
Angle	▶	0°	0°	+8°	+10°	+3°	+3°	+10°	+8°	0°	0°
Replacement .018" (10 UN)	▶	10.25.125	10.25.107	10.25.119	10.25.123	10.25.101	10.25.102	10.25.124	10.25.120	10.25.108	10.25.126



Length	▶	2.5mm	2.5mm	2.7mm	2.2mm	2.2mm	2.2mm	2.2mm	2.7mm	2.5mm	2.5mm
Torque	▶	0°	0°	+7°	0°	0°	0°	0°	+7°	0°	0°
Angle	▶	0°	0°	+2°	0°	0°	0°	0°	+2°	0°	0°
Replacement .018" (10 UN)	▶	10.25.115	10.25.115	10.25.121	10.25.111	10.25.111	10.25.111	10.25.111	10.25.122	10.25.116	10.25.116


Bracket	Prescription	Use	Long 	Torque	Qty	Bonding		Welding	
						REF SLOTT 22	REF SLOTT 18	REF SLOTT 22	REF SLOTT 18
	Edgewise/ Rickets	$\frac{12 \ 11}{21 \ 22}$	3,0 mm	0°	10	10.30.202	10.30.102	-	-
	Edgewise	$\frac{42 \ 41}{31 \ 32}$	1,8 mm	0°	10	10.30.203	10.30.103	-	-
	Edgewise/ Rickets	$\frac{12 \ 22}{42 \ 41 \ 31 \ 32}$	2,3 mm	0°	10	10.30.205	10.30.104	-	-
	Edgewise	$\frac{15 \ 14 \ 13}{45 \ 44 \ 43} \mid \frac{23 \ 24 \ 25}{33 \ 34 \ 35}$	2,5 mm	0°	10	10.30.206	-	-	-
	Edgewise/ Rickets	$\frac{15 \ 14 \ 13}{45 \ 44 \ 43} \mid \frac{23 \ 24 \ 25}{33 \ 34 \ 35}$	2,3 mm	0°	10	10.30.207	10.30.105	-	-
	Edgewise/ Rickets	$\frac{12^* \ 22^*}{42 \ 41^{**} \ 31 \ 32^{**}}$ *Edgewise **Edgewise / Rickets	2,3 mm	0°	10	10.30.209	10.30.107	-	-
	Edgewise	$\frac{16}{46} \mid \frac{26}{36}$	3,4 mm	0°	10	10.30.210	10.30.108	-	-
	Edgewise/ Rickets	$\frac{15 \ 14 \ 13}{45 \ 44 \ 43} \mid \frac{33 \ 34 \ 35}{23 \ 24 \ 25}$	3,0 mm	0°	10	10.30.211	10.30.111	-	-
		$\frac{15 \ 14 \ 13}{45 \ 44 \ 43} \mid \frac{33 \ 34 \ 35}{23 \ 24 \ 25}$				10.30.212	10.30.112		
	Canal Vertical	-	1,8 mm	0°	10	10.40.201	-	-	-
	Edgewise	$\frac{11}{21}$	3,4 mm	0°	10	-	-	10.31.201	-
	Edgewise	$\frac{15 \ 14 \ 13}{45 \ 44 \ 43} \mid \frac{23 \ 24 \ 25}{33 \ 34 \ 35}$	2,6 mm	0°	10	-	-	10.31.206	-
	Edgewise	$\frac{15 \ 14 \ 13}{45 \ 44 \ 43} \mid \frac{23 \ 24 \ 25}{33 \ 34 \ 35}$	2,3 mm	0°	10	-	-	10.31.207	-

Discover the main differences between the **Max** and **Light** lines

INJECTED X MACHINED

The Max line is produced using the most modern MIM (Metal Injection Molding) technology, which results in a single-body piece, made of a single material, without joints or seams.

In the Light lines, the bracket bodies are produced using the machining process, followed by welding of the mesh base and hooks.



PART SIZE



The Max Line is an evolution of the Light Line, and therefore follows very similar dimensions.

The Max brackets are only larger than the Nano Line, which actually has reduced dimensions.

ACCESSORY ADAPTATION

When compared to the Light line, the hooks in the Mc line have a longer shaft, allowing better adaptation of intraoral elastics, metal ligatures and springs during anchoring mechanics.

BASE WITH MICRO PINS X SCREEN BASE




The bases of the Light line feature the traditional retaining screen welded to the bracket body, which is sandblasted to increase roughness, and a careful cleaning process completely removes sandblasting residue.

The bases of the entire Max line, with optimized geometry and surface treatment, are the most effective on the metal bracket market.

OPTIMIZED APPLICATION



Unlike the Light line, the Max line brackets also feature a long axis marking and a "V" notch between the claws, to facilitate the support of the clamp and correct placement during bonding.



Tubes



super kit TUBES

Simple tube with hook Prescription Roth

Prescription	Position	Length	Tube	Torque	Distal Offset	Qty	Bonding REF SLOT 22
Roth	17 16 26 27	3,5 mm	-	-14°	14°	200	20.11.921
				-30°	4°		
	47 46 36 37						

200

TUBES FOR SIMPLE GLUE

Slot 22

All Quadrants Prescription
Roth comm hook

4 Different
Tubes

20.11.211
Single Roth
Tube. 022"



20.11.212
Single Roth
Tube. 022"



20.11.221
Single Roth
Tube. 022"



20.11.222
Single Roth
Tube. 022"



Simple tube Prescription Edgewise

Prescription	Position	Length	Tube	Torque	Distal Offset	Qty	Bonding REF SLOT 22
Roth	17 16 26 27 47 46 36 37	5,2 mm	-	0°	0°	200	20.31.920

200

TUBES FOR SIMPLE GLUE

Slot 22

All Quadrants Edgewise
Prescription



SLTube



CONSECRATED QUALITY AUTOLIGADO MORELLI NOW IN TUBE FORMAT

Far beyond the self-ligating technique, the Morelli Self-ligating Tube combines the ease of inserting the wire into convertible tubes with low-friction mechanics due to its Passive Clip.

MAIN FEATURES OF SLTUBE

- Features Torque and Distal Offset at the base of the piece;
- Big foot base with micropins to ensure excellent retention during treatment;
- Easy insertion of wires, as it is not necessary to fit the wire into the slot entrance;
- Versatility: Can be used in Capellozza, MBT and Roth prescriptions;
- Nitinol Clip: resistance and elastic memory;
- Practical: does not require the application of ties;
- Economy: save chair time by opening and closing the Clip.

KIT

Description	Position		REF
SLTube .022" Glue Case - 4 Quadrants	17 16	26 27	20.13.900
	47 46	36 37	

Repositions on Page 38



PRESCRIÇÃO ROTH / Prescription Roth

Prescription	Position	Length	Tube	Torque	Distal	Qty	Bonding		Welding	
							REF	SLOT 22	REF	SLOT 18

Simple Self-ligating Tube - SLTube

*	17 16	3,5 mm	-	-10°	10°	05	20.13.211	-	-	-
*	26 27	3,5 mm	-	-10°	10°	05	20.13.212	-	-	-
*	36 37	3,5 mm	-	-25°	0°	05	20.13.221	-	-	-
*	47 46	3,5 mm	-	-25°	0°	05	20.13.222	-	-	-

*Compatível com as prescrições Roth, M.B.T. e Capellozza.

Simple tube with hook

Roth	17 16	3,5 mm	-	-14°	14°	10	20.11.211	-	-	-
Roth	26 27	3,5 mm	-	-14°	14°	10	20.11.212	-	-	-
Roth	36 37	3,5 mm	-	-30°	4°	10	20.11.222	-	-	-
Roth	47 46	3,5 mm	-	-30°	4°	10	20.11.221	-	-	-

Simple tube with hook

Roth	17 16	3,5 mm	-	-14°	14°	10	-	20.11.111	20.10.211	-
Roth	26 27	3,5 mm	-	-14°	14°	10	-	20.11.112	20.10.212	-
Roth	36 37	3,5 mm	-	-30°	4°	10	-	20.11.122	20.10.222	-
Roth	47 46	3,5 mm	-	-30°	4°	10	-	20.11.121	20.10.221	-

Simple mini tube with hook

Roth	17	2,5 mm	-	-14°	0°	10	20.11.213	20.11.113	-	-
Roth	27	2,5 mm	-	-14°	0°	10	20.11.214	20.11.114	-	-
Roth	37	2,5 mm	-	-25°	0°	10	20.11.224	20.11.124	-	-
Roth	47	2,5 mm	-	-25°	0°	10	20.11.223	20.11.123	-	-

Convertible Simple Bracket Tube

Roth M.B.T Capellozza	44 34	2,5 mm	-	-17°	0°	10	20.11.225	-	-	-
-----------------------------	-------	--------	---	------	----	----	-----------	---	---	---

Convertible double tube with hook

Roth	36	3,5 mm	-	-30°	4°	10	20.11.412	20.11.312	20.10.412	-
Roth	46	3,5 mm	-	-30°	4°	10	20.11.411	20.11.311	20.10.411	-

Prescription	Position	Length	Tubo	Torque	Distal	Qty	Bonding		Welding	
							REF	SLOT 22	REF	SLOT 18

Double convertible tube with hook

Roth	36	3,5 mm	cervical 1,15 mm .045"	-30°	4°	10	-	-	20.10.414	20.10.314
Roth	46	3,5 mm	cervical 1,15 mm .045"	-30°	4°	10	-	-	20.10.413	20.10.313



Double convertible tube with hook

Roth	16	3,5 mm	-	-14°	14°	10	20.11.417	-	20.10.417	-
Roth	26	3,5 mm	-	-14°	14°	10	20.11.418	-	20.10.418	-



Triple convertible tube with hook

Roth	16	3,5 mm	oclusal 1,15 mm .045"	-14°	14°	10	20.11.611	-	20.10.611	20.10.511
Roth	26	3,5 mm	oclusal 1,15 mm .045"	-14°	14°	10	20.11.612	-	20.10.612	20.10.512



Double convertible tube with hook

Especial	17 16	3,5 mm	-	0°	10°	10	20.11.419	-	-	-
Especial	26 27	3,5 mm	-	0°	10°	10	20.11.420	-	-	-



Compatible with prescription Roth

Simple convertible tube with hook

Capellozza Roth	17 16	3,5 mm	-	-10°	10°	10	20.19.200	-	20.19.100	-
Capellozza Roth	26 27	3,5 mm	-	-10°	10°	10	20.19.201	-	20.19.101	-
Capellozza Roth	36 37	3,5 mm	-	-25°	5°	10	20.19.203	-	20.19.103	-
Capellozza Roth	47 46	3,5 mm	-	-25°	5°	10	20.19.202	-	20.19.102	-



Double convertible tube with hook

Capellozza (Padrão II)	17 16 26 27	3,5 mm	-	-10°	0°	10	20.19.301	-	20.19.300	-
------------------------	---------------	--------	---	------	----	----	-----------	---	-----------	---



Triple convertible tube with hook

Capellozza (Padrão II)	16 26	3,5 mm	oclusal 1,15 mm .045"	-10°	0°	10	-	-	20.19.500	-
------------------------	---------	--------	-----------------------------	------	----	----	---	---	-----------	---



PRESCRIPTION EDGEWISE

Prescripción	Posición	Long.	Tubo	Torque	Distal	Ctd.	Encolado		Soldadura	
							REF SLOTT 22	REF SLOTT 18	REF SLOTT 22	REF SLOTT 18



Simple tube

Edgewise	$\frac{17}{47} \frac{16}{46}$ $\frac{26}{36} \frac{27}{37}$	5,2 mm	-	0°	0°	10	20.31.200	20.31.100	-	-
----------	---	--------	---	----	----	----	-----------	-----------	---	---



Edgewise	$\frac{17}{47} \frac{16}{46}$ $\frac{26}{36} \frac{27}{37}$	5,2 mm	-	0°	0°	10	-	-	20.30.200	20.30.100
----------	---	--------	---	----	----	----	---	---	-----------	-----------



Simple special tube with hook

Edgewise	$\frac{17}{47} \frac{16}{46}$ $\frac{26}{36} \frac{27}{37}$	3,3 mm	-	0°	0°	10	20.31.201	-	20.30.201	-
----------	---	--------	---	----	----	----	-----------	---	-----------	---



Simple tube with hook

Edgewise	$\frac{17}{47} \frac{16}{46}$ $\frac{26}{36} \frac{27}{37}$	5,2 mm	-	0°	0°	10	20.31.230	-	20.30.230	-
Edgewise	$\frac{17}{47} \frac{16}{46}$ $\frac{26}{36} \frac{27}{37}$	5,2 mm	-	0°	0°	10	20.31.260	-	20.30.260	-



Simple tube with hook

Edgewise	$\frac{17}{47} \frac{16}{46}$ $\frac{26}{36} \frac{27}{37}$	3,5 mm	-	0°	0°	10	20.31.241	20.31.141	20.30.241	-
----------	---	--------	---	----	----	----	-----------	-----------	-----------	---



Simple convertible tube with hook

Edgewise	$\frac{17}{47} \frac{16}{46}$ $\frac{26}{36} \frac{27}{37}$	3,5 mm	-	0°	0°	10	20.31.280	20.31.180	20.30.280	20.30.180
----------	---	--------	---	----	----	----	-----------	-----------	-----------	-----------



Double tube with hook

Edgewise	$\frac{17}{47} \frac{16}{46}$ $\frac{26}{36} \frac{27}{37}$	3,5 mm	-	0°	0°	10	20.41.263	-	20.40.263	-
Edgewise	$\frac{17}{47} \frac{16}{46}$ $\frac{26}{36} \frac{27}{37}$	3,5 mm	-	0°	0°	10	*20.41.363	-	*20.40.363	-

*Slot .018" for occlusal and .022" for cervical

Double tube with hook

Tubo doble con gancho

Ricketts	$\frac{17}{47} \frac{16}{46}$ $\frac{26}{36} \frac{27}{37}$	3,5 mm	-	0°	0°	10	-	20.41.163	-	20.40.163
----------	---	--------	---	----	----	----	---	-----------	---	-----------

Prescripción	Posición	Long.	Tubo	Torque	Distal	Ctd.	Encolado		Soldadura	
							REF	SLOT 22	REF	SLOT 18

Double convertible tube with hook

Edgewise Ricketts	16 46	26 36	3,5 mm	cervical 1,15 mm .045"	0°	0°	10	-	-	20.40.331	20.40.330
Edgewise Ricketts	16 46	26 36	3,5 mm	cervical 1,30 mm .051"	0°	0°	10	-	-	20.40.431	20.40.430



Double convertible tube

Edgewise Ricketts	16 46	26 36	3,5 mm	-	0°	0°	10	20.41.501	20.41.500	20.40.501	20.40.500
-------------------	----------	----------	--------	---	----	----	----	-----------	-----------	-----------	-----------



Triple tube with hook

Edgewise Ricketts	16 46	26 36	3,5 mm	central 1,15 mm .045"	0°	0°	10	-	-	20.50.181	20.50.180
-------------------	----------	----------	--------	-----------------------------	----	----	----	---	---	-----------	-----------



Simple tube with hook

M.B.T.	37 47	47	3,5 mm	-	-10°	0°	10	20.35.215	-	-	-
--------	----------	----	--------	---	------	----	----	-----------	---	---	---



Simple convertible tube

M.B.T.	17 26	16 27	3,5 mm	-	-14°	10°	10	20.35.113	-	20.35.103	-
M.B.T.	26 46	27 36	3,5 mm	-	-14°	10°	10	20.35.114	-	20.35.104	-



M.B.T.	46 46	36 36	3,5 mm	-	-20°	0°	10	20.35.212	-	20.35.202	-
--------	----------	----------	--------	---	------	----	----	-----------	---	-----------	---



Double convertible tube with hook

M.B.T.	46 46	36 36	3,5 mm	-	-20°	0°	10	-	-	20.35.401	-
--------	----------	----------	--------	---	------	----	----	---	---	-----------	---



PRESCRIPCIÓN

Prescription	Position	Length	Tube	Torque	Distal	Qty	Bonding		Welding	
							REF	SLOT 22	REF	SLOT 18



Double convertible tube with hook

M.B.T.	16	3,5 mm	cervical 1,15 mm .045"	-14°	10°	10	-	-	20.35.402	-
M.B.T.	26	3,5 mm	cervical 1,15 mm .045"	-14°	10°	10	-	-	20.35.403	-



Double convertible tube with hook for P.L.A.

M.B.T.	46 36	3,5 mm	cervical 1,15 mm .045"	-20°	0°	10	-	-	20.35.405	-
--------	---------	--------	------------------------------	------	----	----	---	---	-----------	---



Double convertible tube with hook

M.B.T.	46 36	3,5 mm		-20°	0°	10	20.35.416	-	-	-
--------	---------	--------	--	------	----	----	-----------	---	---	---



Triple convertible tube with hook

M.B.T.	16	3,5 mm	oclusal 1,15 mm .045"	-14°	10°	10	-	-	20.35.501	-
M.B.T.	26	3,5 mm	oclusal 1,15 mm .045"	-14°	10°	10	-	-	20.35.502	-

LINGUAL

Prescription	Position	Length	Torque	Distal	Qty	Welding 2 x 0,90mm 2 x .036"
--------------	----------	--------	--------	--------	-----	------------------------------------



Lingual tube with hook

Lingual	16 26 46 36	3,5 mm	0°	0°	10	20.61.206
---------	--------------------	--------	----	----	----	-----------

Accessories



Type	Qty	REF
Orthodontics accessory kit		
32 types of orthodontic accessories	240	30.90.900

super kit

Accessories
of orthodontics

32

Type Kit
of products
with 240 pieces.

The most orthodontic brand
complete with the main
accessories in one kit



Replacement Reference Codes

20.11.225  10 un. 01	20.61.206  10 un. 02	30.06.303  06 un. 03	30.10.205  30 un. 04	30.10.210  05 un. 05	30.20.101  10 un. 06	30.21.104  10 un. 07	30.21.107  05 un. 08
30.40.003  10 un. 09	30.40.004  10 un. 10	30.40.005  10 un. 11	30.40.008  10 un. 12	30.40.009  10 un. 13	30.40.010  10 un. 14	30.40.016  10 un. 15	30.40.021  20 un. 16
30.45.011  02 un. 17	30.45.012  02 un. 18	30.45.013  02 un. 19	30.45.014  02 un. 20	30.45.015  02 un. 21	30.45.016  02 un. 22	30.60.001  10 un. 23	30.60.002  10 un. 24
35.20.064  04 un. 25	35.20.065  04 un. 26	35.20.066  04 un. 27	35.20.067  04 un. 28	37.30.001  04 un. 29	37.30.002  04 un. 30	37.30.003  04 un. 31	37.30.004  04 un. 32

Number of pieces per cavity in this package



Comes with box organizer with brands numbered for all accessories.




Anchorage hook

Used in the orthodontic arch for anchoring chain elastic, springs, metal ligature wires, to assist treatments involving retraction, mesialization and anchorage. Fixation is better in rectangular arches because they are larger.


Type	Size	Qty	REF
------	------	-----	-----

Universal hook for TAD


slot.022"		05	30.40.010
-----------	---	----	-----------

Adjustable height hook, remove hooks by flexing to the desired size and finish


Ball curved tip hook

Right slot.022"		10	30.40.003
Left slot.022"			30.40.004

Ball straight tip hook

slot.022"		10	30.40.005
-----------	---	----	-----------

Ball universal tip hook

slot.022"		10	30.40.009
-----------	--	----	-----------

Type	Size	Qty	REF
------	------	-----	-----

Gurin with closed round slot

Universal	.022" ajustável adjustable	10	35.05.001
-----------	----------------------------	----	-----------

Gurin with closed round slot and hook

Right	.022" ajustável adjustable	10	35.05.002
Left			35.05.003

Gurin with closed round slot and right hook

Right	.022" ajustável adjustable	05	35.05.004
Left			35.05.005

Gurin with open rectangular round hook

Universal	.022" ajustável adjustable	10	35.05.011
-----------	----------------------------	----	-----------

Gurin with open rectangular slot and hook

Right	.022" ajustável adjustable	10	35.05.012
Left			35.05.013

Gurin


Son aparatos que, aplicados a las arcadas de ortodoncia, con la ayuda de la llave hexagonal*, permiten el anclaje de los aparatos directamente en el arco.



How to use Gurin


Hook for traction

The hook is bonded to the impacted tooth whose current remains anchored in the intraoral arch. The tooth extrusion is seen as the links appear, which, once they are counted, they identify the tooth movement and positioning.

Size	Qty	REF
Hook for traction of non-erupted teeth with chain		
	10	30.21.107




Hook for welding

	100	30.20.101
---	-----	-----------



Hook for traction of non-erupted teeth

	10	30.21.104
---	----	-----------





Button

Accessory used in dental movements such as: gyroversions, tractions, retractions or anchoring systems. They can retain elastic and ligature wires.



When to use the orthodontic button.


Type	Size	Qty	REF
Button			
For Welding		10	30.10.105
Convex for bonding		10	30.10.202
For Convex for bonding		10	30.10.203
For bonding		10	30.10.205
Composite for bonding		10	30.10.206
Ceramic for bonding		02	30.10.210





Size	Qty	REF
------	-----	-----

Stop for rectangular arches

 .022"	10	30.40.008
---	----	-----------

For Rectangular arches:
Effective to restrain or limit the movement of arches and orthodontic accessories, in addition to avoiding the discomfort caused by the displacement of the arch at the distal end




Stop for round arches

 .021" / Ø 0.54mm	25	30.40.020
--	----	-----------

For round arches:
Designed for midline demarcation or locking round orthodontic arches, mainly avoiding discomfort caused by displacement of the arch at the distal tube end. It limits lateral movement without causing deformation. It is used in all orthodontic techniques including self-ligating brackets.



Stop closed with applicator for wires up to .022"

 .021" / Ø 0.54mm	10	30.40.021
--	----	-----------

Check how to apply.



Size	Qty	REF
------	-----	-----

Stop and bite guide

 Base / Base 3.6mm x 4.0mm	10	30.60.001
--	----	-----------

The main function of the Stop and Guide is to reduce the impact of the upper teeth with the lower brackets, as it works by opening the bite, and at the same time provides reeducation and muscle rest.



Size	Qty	REF
------	-----	-----

Mini spur tongue educator

 Base / Base 2.8mm x 3.8mm	10	30.60.002
--	----	-----------

The Mini Spur induces the tongue to reduce the pressure on the anterior, upper and/or lower teeth, by assisting in digital suction treatments (finger sucking), correction of open bites in re-education of relapses.



Tube for activators with bar

Ø1,15mm Ø0.45"	10	30.05.145
Ø1,30mm Ø0.51"	10	30.05.151


Resorte

Los resortes Morelli proporcionan fuerzas ligeras y continuas, lo que permite intervalos de activación más prolongados. Son excelentes dispositivos para obtener fuerzas de tracción o compresión.

Los resortes de tracción cerrados cuentan con extremos preparados para adaptarse a los accesorios


Type	Length	Size	Qty	REF
Spring CrNi				
Open compression	1 meter	 ø interno / internal 0,80mm (.032")	01	35.20.060
Closed traction	1 meter	 ø interno / internal 0,80mm (.032")		35.20.061

Spring NiTi

Open compression	178mm	 ø interno / internal 0,76mm (.030")	03	35.20.062
------------------	-------	---	----	-----------

Spring NiTi open and closed

Open and closed compression	178mm	 ø interno / internal 1,00mm (.040")	03	35.20.072
-----------------------------	-------	---	----	-----------



To avoid the open spiral tip, it is recommended to cut the spring in the closed spiral

Type	Length	Size	Size	Qty	REF
------	--------	------	------	-----	-----

Resorte NiTi cerrado

Closed traction	7mm	10 a 15 mm	200gf	04	35.20.064
	9mm	15 a 25 mm			35.20.065
	12mm	25 a 35 mm			35.20.066
	15mm	35 a 45 mm			35.20.067



Spring NiTi closed for mini screw

Closed traction	7mm	10 a 15 mm	200gf	04	37.30.001
	9mm	15 a 55 mm			37.30.002
	12mm	25 a 35 mm			37.30.003
	15mm	35 a 45 mm			37.30.004



Containment

Orthodontic retention is used to maintain teeth in their new positions after orthodontic movement with the substantial function of providing stability of the results obtained once the orthodontic treatment is completed.

Type	Length	Size	Qty	REF
TwistForm Containment 3x3				
3 fios 	41,5 mm	 ø externo / external 0,50mm (.020")	10	30.06.303
3 fios 	41,5 mm	 ø externo / external 0,65mm (.026")	10	30.06.304



Check out the twistform application



Sliding Jig

Size	Type	Qty	REF
------	------	-----	-----

Sliding Jig for Mesialization and Distalization

12mm Ø1,30mm	Corto - Superior. Izquierdo / Inf. Derecho	05	30.45.011
12mm	Corto - Superior. Izquierdo / Inf. Derecho		30.45.012
18mm	Medio - Superior. Izquierdo / Inf. Derecho	05	30.45.013
18mm	Medio - Superior. Izquierdo / Inf. Derecho		30.45.014
32mm	Largo - Sup. Izquierdo / Inf. Derecho	05	30.45.015
32mm	Largo - Sup. Izquierdo / Inf. Derecho		30.45.016



Morelli Sliding Jig is an intraoral device for installation in the intraoral arch, which allows the anchoring springs or elastics in an adjustable position hook. The clinical force adjustment is performed by positioning the hook along the Ø 0.7mm sliding jig wire, using the Gurin*

*Gurin key page 103 REF: 75.02.011





Check the sliding
jig application.

Distalization



Mesialization



Type	Size	Qty	REF
------	------	-----	-----

Cross tube

-	.022"	10	30.40.006
---	-------	----	-----------

Cervical offset cross tube

-	.022"	10	30.40.016
---	-------	----	-----------

Cross Tube

Cross tube is fixed to the orthodontic arch to receive wires up to .021"x.025". Attached to the main arch, it allows the installation of segmented arches, cantilevers, handles and others. It assists in space closure, anterior segment intrusion and loss of anchorage.

Cervical offset cross tube
The tube cervical offset prevents the contact of the opposing tooth on the wire inserted in the hole.




Distalizer

This device is used to distalize maxillary molars and offers constant forces. The activation of the nitinol spring releases 60 to 80gf forces when subjected to 5 mm medium compressions.

Size (L)	Compression*	Force	Qty	REF
----------	--------------	-------	-----	-----

Molar distalizer type Jones universal

15 mm	3 mm	60gf	02	30.45.003
	4 mm	65gf		
	5 mm	70gf		
	6 mm	75gf		
	7 mm	80gf		



*Medium Values

Spring NiTi



Movable / occlusal eyelet for anchorage

Telescopic tube

Made of stainless steel, they are used in the composition of some orthopedic and orthodontic appliances having sliding movements on a wire or rod. Such as Mandibular Protraction Appliances (MPA), used to correct class II malocclusions

Indicated for wires up to	Qty	Length	REF
---------------------------	-----	--------	-----

Telescopic tube

diâmetro interno	0,80mm	05	330mm	30.40.030
	0,90mm			30.40.031
	1,40mm			30.40.032



Headed ball clasp

The Headed Ball Clasp is used in the manufacture of removable braces to create "retention" points through its ball which is attached to the teeth by pressure.

Dimensions	Qty	Length	REF
------------	-----	--------	-----

Headed ball clasp

Haste: ø 0,8mm; Bola: ø 1,5mm	50	38mm	30.40.041
----------------------------------	----	------	-----------



Palatal bar



The palatal bar is a device that connects to the lingual tubes. Used for different purposes, it can generate forces for gyroversion, intrusion or torque, as well as working as an anchorage reinforcement.

Size	Qty	REF - Mesial	REF - Distal
------	-----	--------------	--------------

Palatal bar with mesial and distal loop

*34mm	ø0,9mm	10	30.30.034	30.30.134
37mm			30.30.037	30.30.137
39mm			30.30.039	30.30.139
41mm			30.30.041	30.30.141
43mm			30.30.043	30.30.143
45mm			30.30.045	30.30.145
47mm			30.30.047	30.30.147
49mm			30.30.049	30.30.149
52mm			30.30.052	30.30.152
55mm			30.30.055	30.30.155

Fio Ø 0,9mm





Type	Size	Qty	REF
------	------	-----	-----

Orthodontic wire mesh

CrNi	Fio = $\varnothing 0,12$ mm 100 x 100mm	01	35.11.001
------	--	----	-----------

Orthodontic wire mesh

Used to make a base for bonding tubes and other parts. It uses stainless steel braided wires with 80 wires/inch.



Type	Size	Qty	REF
------	------	-----	-----

Silver solder

56% (Ag) Fusion = 650°C	Rolo / roll $\varnothing 0,05$ mm 5m	01	35.30.010
	Vareta / Stick $\varnothing 0,8$ mm 300mm	04	35.30.011
	Vareta / Stick $\varnothing 0,8$ mm 300mm	10	35.30.013

Silver solder

Silver solder wire, easy to apply, toxic material-free.



Size	Qty	REF
------	-----	-----

Solder flux

20 ml	01	35.30.020
-------	----	-----------

Solder flux

Material used as an indispensable flux to join metals through brazing, allowing uniform and bubble-free welds



Qty	REF
-----	-----

Orthodontic wax protection

10 cards	35.33.001
----------	-----------

Orthodontic wax

The Protection Wax developed by Morelli offers comfort during adaptation to the orthodontic appliance, in addition to avoiding emergency appointments and increasing patient collaboration at the beginning of the treatment. Flavored with a mild mint flavor

1 card	35.33.002
--------	-----------

Thermoforming Plates



Plates for Thermoforming

Clarís Plates are round in shape with Ø125mm. They are manufactured in PET-G laminated plate. They are used in the manufacture of orthodontic appliances through the thermoforming process.

They are resistant and safe, for formability by vacuum heating or pressure. In addition, they have excellent final shape memory, with high ductility, as they do not fracture under the forces applied during use.

Measures	Qty	REF
----------	-----	-----

Thermoforming Plates

Clarís - PET-G plate 0.50 mm - Ø125mm	10	66.01.050
Clarís - PET-G plate 0.76 mm - Ø125mm	10	66.01.076
Clarís - PET-G plate 1.00 mm - Ø125mm	10	66.01.100



Dental applications

- Aesthetic aligners
- Retainer plates
- Indirect bonding guides
- Bruxism plates
- Tins for whitening or fluoridation



Check how to use it

Absolute Anchorage





TITANIUM MINI SCREW WITH BUTTON HEAD

Ideal for inter-radicular application, it allows several possibilities of movement during orthodontic treatment, regardless of the patient cooperation.

TITANIUM MINI SCREW WITH BRACKET HEAD

In addition to anchoring springs, elastics or ligatures, it allows the use of cantilevers, and thus it offers new possibilities for orthodontic movement. The Bracket Head Mini Screw has two 0.021" x 0.025" slots perpendicular to each other in the screw head to accommodate rectangular wires precisely.

STEEL MINI SCREW WITH BUTTON HEAD

Ideal for extra-radicular application, the Steel Mini Screw, when compared to titanium, it presents more strength to bending and torsion.


MAIN CHARACTERISTICS OF MORELLI MINI SCREWS

- ▶ Active tip with self-drilling profile, specially designed to offer stability, secure fit in bone texture and strength for immediate orthodontic load applications;
- ▶ Transmucosal section of perfect adaptation and maintenance of peri-implant soft tissues;
- ▶ High biocompatibility, material tested and approved by ANVISA regarding its biological safety;
- ▶ 0.8mm hole (indicated for wire up to .021"x0.25") in the hex head, sized especially for the application of orthodontic accessories;
- ▶ Easy application, all Mini Screws are compatible with the same Morelli Tool kit;
- ▶ Availability, available in the main dental offices in Brazil with guaranteed supply at any time;
- ▶ Fair price, hoping to contribute to the evolution of Brazilian Orthodontics, we try to practice an affordable price.




Threading Length	Transmucosal	Diameter	Insertion		Qty	REF
------------------	--------------	----------	-----------	--	-----	-----

Titanium Orthodontic Mini Screw with Button Head




6mm	1mm	Ø1.5mm		Clockwise	2	37.10.101
8mm						37.10.102




6mm	2mm	Ø1.5mm		Clockwise	2	37.10.201
8mm						37.10.202
10mm						37.10.203



Titanium Orthodontic Mini Screw with Bracket Head




6mm	1mm	Ø1.5mm		Counter Clockwise	1	37.10.231
8mm						37.10.232



6mm	1mm	Ø1.5mm		Clockwise	1	37.10.241
8mm						37.10.242



Steel Orthodontic Mini Screw with Button Head



5mm	4mm	Ø2.0mm		Clockwise	1	37.10.405
7mm	4mm	Ø2.0mm		Clockwise	1	37.10.407
10mm	4mm	Ø2.0mm		Clockwise	1	37.10.410
12mm	4mm	Ø2.0mm		Clockwise	1	37.10.412
14mm	4mm	Ø2.0mm		Clockwise	1	37.10.414



All Morelli Mini Screws are sterilized.



Spear Tip

Size	Qty	REF
60mm	1	37.20.010

Wide Shank Hexagonal Wrench

41mm	1	37.20.011
------	---	-----------

Average Haste Hex Key

24mm	1	37.20.015
------	---	-----------

Short Hexagonal Haste Wrench

8mm	1	37.20.012
-----	---	-----------

Key Hexagonal Digital

25mm	1	37.20.013
------	---	-----------

Key Hand

92mm	1	37.20.014
------	---	-----------

Contra angle wrench

59mm	1	37.20.016
------	---	-----------

Tool Kit

100mm 60mm	1	37.20.001
---------------	---	-----------

Tool Kit Items: Spear Tip, Long Hex Wrench, Medium and Short, Digital Hexagonal Wrench and Hand Wrench.

Spear Tip

Used to facilitate entry into the gums and marking of the bone to be drilled by the Mini Screw.

Key Hexagonal

Long, medium and short to fit the Miniscrew. Then equipped with Nitinol retaining ring durable and safe.

Key Hexagonal Digital

Digitally, this wrench can be used in reduced spaces, intended for installing and removing Orthodontic Mini Screws.

Key Hand

Solid construction, fits comfortably in the palm of your hand and frees your fingers to turn the Mini Screw safely and firmly.

Contra angle wrench

Coupled to the dental contra-angle, the Key replaces the low-speed micromotor, transforming the contra-angle into an installation key angled.

Tool Kit


Complete set of tools whose organizing case is made of engineering polymer, resistant to deformation caused by high autoclave temperatures.

Bands



Universal band for molar

Universal bands have internal perimeters ranging from 28 to 42mm, corresponding to the teeth perimeters..



Perimeter mm	Qty	REF	
		47 46 36 37	17 16 26 27
28	10	40.01.280	40.02.280
28,5	10	40.01.285	40.02.285
29	10	40.01.290	40.02.290
29,5	10	40.01.295	40.02.295
30	10	40.01.300	40.02.300
30,5	10	40.01.305	40.02.305
31	10	40.01.310	40.02.310
31,5	10	40.01.315	40.02.315
32	10	40.01.320	40.02.320
32,5	10	40.01.325	40.02.325
33	10	40.01.330	40.02.330
33,5	10	40.01.335	40.02.335
34	10	40.01.340	40.02.340
34,5	10	40.01.345	40.02.345
35	10	40.01.350	40.02.350
35,5	10	40.01.355	40.02.355
36	10	40.01.360	40.02.360
36,5	10	40.01.365	40.02.365
37	10	40.01.370	40.02.370
37,5	10	40.01.375	40.02.375
38	10	40.01.380	40.02.380
38,5	10	40.01.385	40.02.385
39	10	40.01.390	40.02.390
39,5	10	40.01.395	40.02.395
40	10	40.01.400	40.02.400
40,5	10	40.01.405	40.02.405
41	10	40.01.410	40.02.410
41,5	10	40.01.415	40.02.415
42	10	40.01.420	40.02.420

Universal band for molar

Morelli anatomical bands reproduce the morphology of the molars, allowing optimal adaptation. Available for upper and lower molars, right and left sides.



Perimeter mm	Qty	REF	
		47 46 36 37	17 16 26 27
29	10	40.12.290	40.22.290
29,5	10	40.12.295	40.22.295
30	10	40.12.300	40.22.300
30,5	10	40.12.305	40.22.305
31	10	40.12.310	40.22.310
31,5	10	40.12.315	40.22.315
32	10	40.12.320	40.22.320
32,5	10	40.12.325	40.22.325
33	10	40.12.330	40.22.330
33,5	10	40.12.335	40.22.335
34	10	40.12.340	40.22.340
34,5	10	40.12.345	40.22.345
35	10	40.12.350	40.22.350
35,5	10	40.12.355	40.22.355
36	10	40.12.360	40.22.360
36,5	10	40.12.365	40.22.365
37	10	40.12.370	40.22.370
37,5	10	40.12.375	40.22.375
38	10	40.12.380	40.22.380
38,5	10	40.12.385	40.22.385
39	10	40.12.390	40.22.390
39,5	10	40.12.395	40.22.395
40	10	40.12.400	40.22.400
40,5	10	40.12.405	40.22.405
41	10	40.12.410	40.22.410
41,5	10	40.12.415	40.22.415
42	10	40.12.420	40.22.420



pieces RL + 5 pieces LL

5 pieces RU + 5 pieces LU



Band for upper premolar



upper

The premolar bands are anchoring devices used in maxillary expansion treatments associated with Haas and Hyrax, Jones Jig expander appliances, and Nance buttons. It can also be used in cases of recurrent bracket falls, where it is welded to the band.

Size	Qty	REF
		15 14 25 24
12	10	40.50.012
13	10	40.50.013
14	10	40.50.014
15	10	40.50.015
16	10	40.50.016
17	10	40.50.017
18	10	40.50.018
19	10	40.50.019
20	10	40.50.020
21	10	40.50.021
22	10	40.50.022
23	10	40.50.023
24	10	40.50.024
25	10	40.50.025
26	10	40.50.026
27	10	40.50.027
28	10	40.50.028
29	10	40.50.029
30	10	40.50.030
31	10	40.50.031
32	10	40.50.032



See bands seater on page 97

To order your Anatomical Band with Welded Tube, follow these steps:

Go to the page to place your order

[illegible]

Step 2

Select the Tube

Step 3


Identifique el código de intersección para realizar su pedido

Examples:

You will receive 1 Pair of the selected set (Right/Left).

Familia 47
Banda con Tubo

Family 48
Band with Tube



Code
47.10.007
2 un





Code
48.31.109
2 un

ANATOMICAL BAND with tube Lower				TUBES					
				SIMPLE	DOUBLE		DOUBLE + LINGUAL		
Prescription				ROTH	ROTH	ROTH *L.B.	M.B.T.	ROTH	LINGUAL
Torque				-30°	-30°	-30°	-20°	-30°	-
Offset				4°	4°	4°	0°	4°	-
Tube ø				-	-	.045"	-	-	-
Slot				.022"	.022"	.022"	.022"	.022"	2x .036"
									
Perímetro mm	Bandas	R	Ctd.	20.10.221	20.10.411	20.10.413	20.35.401	20.10.411	20.61.206
		L		20.10.222	20.10.412	20.10.414	-	20.10.412	20.61.206
32	40.12.320	R	2	47.10.007	47.20.007	-	-	-	-
32,5	40.12.325	R	2	47.10.008	47.20.008	-	-	-	-
33	40.12.330	R	2	47.10.009	47.20.009	-	-	-	-
33,5	40.12.335	R	2	47.10.010	47.20.010	-	-	-	-
34	40.12.340	R	2	47.10.011	47.20.011	-	-	-	-
34,5	40.12.345	R	2	47.10.012	47.20.012	47.20.112	47.20.312	48.20.112	-
35	40.12.350	R	2	47.10.013	47.20.013	47.20.113	47.20.313	48.20.113	-
35,5	40.12.355	R	2	47.10.014	47.20.014	47.20.114	47.20.314	48.20.114	-
36	40.12.360	R	2	47.10.015	47.20.015	47.20.115	47.20.315	48.20.115	-
36,5	40.12.365	R	2	47.10.016	47.20.016	47.20.116	47.20.316	48.20.116	-
37	40.12.370	R	2	47.10.017	47.20.017	47.20.117	47.20.317	48.20.117	-
37,5	40.12.375	R	2	47.10.018	47.20.018	47.20.118	47.20.318	48.20.118	-
38	40.12.380	R	2	47.10.019	47.20.019	47.20.119	47.20.319	48.20.119	-
38,5	40.12.385	R	2	47.10.020	47.20.020	47.20.120	47.20.320	48.20.120	-
39	40.12.390	R	2	47.10.021	47.20.021	-	-	-	-
39,5	40.12.395	R	2	47.10.022	47.20.022	-	-	-	-
40	40.12.400	R	2	47.10.023	47.20.023	-	-	-	-
40,5	40.12.405	R	2	47.10.024	47.20.024	-	-	-	-
41	40.12.410	R	2	47.10.025	47.20.025	-	-	-	-
41,5	40.12.415	R	2	47.10.026	47.20.026	-	-	-	-
42	40.12.420	R	2	47.10.027	47.20.027	-	-	-	-

ANATOMICAL BAND with tube Upper				TUBES					
				SIMPLE	DOUBLE		TRIPLE	TRIPLE + LINGUAL	
Prescription				ROTH	ROTH	M.B.T. *L.B.	ROTH	ROTH	LINGUAL
Tanque				-14°	-14°	-14°	-14°	-14°	-
Offset				14°	14°	10°	14°	14°	-
Tube ø				-	-	.045"	.045"	.045"	-
Slot				.022"	.022"	.022"	.022"	.022"	2x .036"
									
Perimeter mm	Bands	R	Qty	20.10.211	20.10.417	20.35.402	20.10.611	20.10.611	20.61.206
		L		20.10.212	20.10.418	20.35.403	20.10.612	20.10.612	-
32	40.22.320	R	2	47.11.007	-	-	-	-	-
32,5	40.22.325	R	2	47.11.008	-	-	-	-	-
33	40.22.330	R	2	47.11.009	-	-	47.31.009	48.31.109	-
33,5	40.22.335	R	2	47.11.010	-	-	47.31.010	48.31.110	-
34	40.22.340	R	2	47.11.011	47.21.011	47.21.111	47.31.011	48.31.111	-
34,5	40.22.345	R	2	47.11.012	47.21.012	47.21.112	47.31.012	48.31.112	-
35	40.22.350	R	2	47.11.013	47.21.013	47.21.113	47.31.013	48.31.113	-
35,5	40.22.355	R	2	47.11.014	47.21.014	47.21.114	47.31.014	48.31.114	-
36	40.22.360	R	2	47.11.015	47.21.015	47.21.115	47.31.015	48.31.115	-
36,5	40.22.365	R	2	47.11.016	47.21.016	47.21.116	47.31.016	48.31.116	-
37	40.22.370	R	2	47.11.017	47.21.017	47.21.117	47.31.017	48.31.117	-
37,5	40.22.375	R	2	47.11.018	47.21.018	47.21.118	47.31.018	48.31.118	-
38	40.22.380	R	2	47.11.019	47.21.019	47.21.119	47.31.019	48.31.119	-
38,5	40.22.385	R	2	47.11.020	47.21.020	47.21.120	47.31.020	48.31.120	-
39	40.22.390	R	2	47.11.021	47.21.021	47.21.121	47.31.021	48.31.121	-
39,5	40.22.395	R	2	47.11.022	-	-	47.31.022	48.31.122	-
40	40.22.400	R	2	47.11.023	-	-	47.31.023	48.31.123	-
40,5	40.22.405	R	2	47.11.024	-	-	47.31.024	48.31.124	-
41	40.22.410	R	2	47.11.025	-	-	47.31.025	48.31.125	-
41,5	40.22.415	R	2	47.11.026	-	-	47.31.026	48.31.126	-
42	40.22.420	R	2	47.11.027	-	-	47.31.027	48.31.127	-

NOTA: *L.B. = Tubo para LIP BUMPER (Placa Lábio Ativa)

Universal band kits

Size	Qty	Lower	REF	Uper	REF	Color
		7 6	6 7	7 6	6 7	
	260	-		40.02.900		
	260	40.01.900		-		
	130	-		40.02.901		
	130	40.01.901		-		

Kits con Bandas

Band kits present, in convenient packaging, all band sizes, in quantities sized according to their use. Standardized in red for upper band packaging and blue for lower bands. The boxes are also sold individually. (see page 112).



Kit composition

PERIMETER mm	28	28,5	29	29,5	30	30,5	31	31,5	32	32,5	33	33,5	34	34,5	35	30,5	36	36,5	37	37,5	38	38,5	39	39,5	40	40,5	41	41,5	42
Qty	4	4	4	4	4	4	4	8	8	8	8	8	10	10	10	14	16	16	20	20	20	18	10	8	4	4	4	4	4

Embalagem com 260 Unidades.

Kit composition

PERIMETER mm	28	28,5	29	29,5	30	30,5	31	31,5	32	32,5	33	33,5	34	34,5	35	30,5	36	36,5	37	37,5	38	38,5	39	39,5	40	40,5	41	41,5	42
Qty	2	2	2	2	2	2	2	4	4	4	4	4	5	5	5	7	8	8	10	10	9	5	4	2	2	2	2	2	2

Embalagem com 130 Unidades.



Blade Thickness x Height	Qty	REF
0,10 x 3,2mm .004" x .125"	8 metros	40.91.001
0,10 x 3,8mm .004" x .150"	8 metros	40.91.002
0,10 x 4,5mm .004" x .180"	8 metros	40.91.003
0,15 x 4,5mm .006" x .180"	8 metros	40.91.004

Banding Tapes in roll

Banding tapes are made of polished and malleable stainless steel, used for bandaging canines and other teeth.

Archs



Morelli intraoral archwire full-size diagram

Diagram CrNi

Cooper, Stainless Steel






Diagrama NiTi

Superelastic, Thermo-plus, Rhodium
y Reverse Curve

We offer three arch sizes to better adapt to different arch sizes and clinical objectives



Intraoral archwire CrNi

Identification		Qty			
Size			Upper	Lower	
					
	.012"	10	50.60.000	50.70.000	
	.014"	10	50.60.001	50.70.001	
	.016"	10	50.60.002	50.70.002	
	.018"	10	50.60.003	50.70.003	
	.020"	10	50.60.004	50.70.004	
	.016" x .016"	10	50.61.001	50.71.001	
	.016" x .022"	10	50.62.001	50.72.001	
	.017" x .025"	10	50.62.002	50.72.002	
	.018" x .025"	10	50.62.003	50.72.003	
	.019" x .025"	10	50.62.004	50.72.004	
	.021" x .025"	10	50.62.005	50.72.005	

Archwire CrNi

Morelli stainless steel orthodontic archwires feature three types of cross section: round, square and rectangular.



Identification		Qty				
Size			Upper	Lower	Small	
						
	.010"	10	-	50.70.010	-	
	.012"	10	50.60.011	50.70.011	50.80.011	
	.014"	10	50.60.012	50.70.012	50.80.012	
	.016"	10	50.60.013	50.70.013	50.80.013	
	.018"	10	50.60.014	50.70.014	50.80.014	
	.020"	10	50.60.015	50.70.015	50.80.015	
	.016" x .016"	10	50.61.002	50.71.002	50.81.002	
	.016" x .022"	10	50.62.011	50.72.011	50.82.011	
	.017" x .025"	10	50.62.012	50.72.012	50.82.012	
	.018" x .025"	10	50.62.013	50.72.013	50.82.013	
	.019" x .025"	10	50.62.014	50.72.014	50.82.014	
	.021" x .025"	10	50.62.015	50.72.015	50.82.015	

Archwire Superelastic NiTi

NiTi (nickel-titanium) alloys are relevant in orthodontic biomechanics, as they have super-elasticity and shape memory effects. They release smooth and constant forces compared to stainless steel. They do not undergo permanent deformation, that is, the wires cannot be bent.

Morelli archwires NiTi are presented in Small, Medium and Large models.



Archwire Super-elastic Reverse Curve NiTi

These archwires are useful in situations where the presence of curve in the archwire is important for the orthodontic biomechanics. They are available in different thicknesses and crosssections (round, square and rectangular)



Intraoral Archwire Super-elastic NiTi Reverse Curve – Spee

Identification				
			Large	Medium
Size		Qty		
	.012"	10	50.60.021	50.70.021
	.014"	10	50.60.022	50.70.022
	.016"	10	50.60.023	50.70.023
	.018"	10	50.60.024	50.70.024
	.020"	10	50.60.034	50.70.034
	.016" x .016"	10	50.61.003	50.71.003
	.016" x .022"	10	50.62.021	50.72.021
	.017" x .025"	10	50.62.022	50.72.022
	.019" x .025"	10	50.62.023	50.72.023

Archwire Super-elastic Aesthetic NiTi

Made of Nickel-Titanium (Nitinol) and covered by high technology polymer. Morelli Intraoral Archwires Aesthetic have a color close to dental enamel, resistant coating, flexible material and biocompatible, resistant to stains and color change.



Intraoral Archwire Super-elastic Aesthetic

Identificación				
			Large	Medium
Size		Qty		
	.012"	05	50.60.035	50.70.035
	.014"	05	50.60.036	50.70.036
	.016"	05	50.60.037	50.70.037
	.018"	05	50.60.038	50.70.038
	.017" x .025"	05	50.62.036	50.72.036
	.019" x .025"	05	50.62.037	50.72.037

*Aesthetic bows equal to and greater than 0.16" only apply to slot 22.

Archwire NiTi Thermo Plus

Thermoactivated NiTi Thermo-Plus alloy undergoes microstructural transformation with thermal change. When the Archwire is cooled, it becomes malleable and inactive. When it reaches the oral environment temperature (above 27°C), it is activated by acquiring shape memory, generating orthodontic force for tooth movements.



Intraoral Archwire Niti Thermo Plus

Identification					
Size		Qty	Large 	Medium 	Small
	.012"	10	50.60.224	50.70.224	50.80.224
	.014"	10	50.60.225	50.70.225	50.80.225
	.016"	10	50.60.226	50.70.226	50.80.226
	.018"	10	50.60.227	50.70.227	50.80.227
	.020"	10	50.60.228	50.70.228	50.80.228
	.016" x .016"	10	50.61.204	50.71.204	50.81.204
	.016" x .022"	10	50.62.224	50.72.224	50.82.224
	.017" x .025"	10	50.62.225	50.72.225	50.82.225
	.018" x .025"	10	50.62.226	50.72.226	50.82.226
	.019" x .025"	10	50.62.227	50.72.227	50.82.227
	.021" x .025"	10	50.62.228	50.72.228	50.82.228

Arco Intraoral Copper Niti

Identification				
Size		Qty	Upper 	Lower
	.014"	10	50.60.150	50.70.150
	.016"	10	50.60.151	50.70.151
	.018"	10	50.60.152	50.70.152
	.014" x .025"	10	50.62.150	50.72.150
	.016" x .022"	10	50.62.151	50.72.151
	.017" x .025"	10	50.62.152	50.72.152
	.018" x .025"	10	50.62.153	50.72.153
	.019" x .025"	10	50.62.154	50.72.154

Archwire Copper NiTi

Morelli Intraoral Archwires Copper NiTi have been specially designed to provide the best clinical performance to the self-ligating bracket system. The 6% copper edition in the alloy composition enhances the thermal properties of nickel-titanium, resulting in the delivery of ideal forces, with incomparable tooth movement efficiency.




Archwire Super-elastic NiTi Rhodium

Developed from archwires super-elastic NiTi, archwires with Rhodium coating provide excellent aesthetics combined with mechanical performance, keeping the release of light and constant clinical forces throughout the treatment. Morelli Archwires Rhodium are coated with a thin layer of rhodium, a precious metal with a silvery white color, which provides a whitish appearance of this metal, similar to the tooth color, ensuring an excellent visual result.



Intraoral Archwire Super-elastic NiTi Rhodium

Identification				
Size		Qty	Large	Medium
	.012"	01		
	.014"	01	50.60.135	50.70.135
	.016"	01	50.60.136	50.70.136
	.018"	01	50.60.137	50.70.137
	.018"	01	50.60.138	50.70.138
	.017" x .025"	01	50.62.136	50.72.136
	.019" x .025"	01	50.62.137	50.72.137

Archwire Kit for Self-ligating

The packaging itself is an explanatory booklet on the application of archwires and their respective forces. The Archwires and the suggestions for use presented do not exempt the careful assessment of the qualified professional, to adapt the treatment to the diagnosis and biological response



Intraoral Archwire Kit for Self-ligating

THE MOST COMPLETE ARCHWIRE SEQUENCE ON THE MARKET

Type	Qty	Size	REF
Super-elastic - NiTi	20 un*	Ø.010"	50.14.901
Thermo-Plus - NiTi		Ø.012" / Ø.014" / Ø.016" .016"x.016" / .016"x.022" .017"x.025" / .021"x.025"	
Chromium Nickel - CrNi		.019"x.025"	
Titanio - Molibdeno - TMA		.017"x.025"	

*10 different pairs

Each size contains 2 archwires, lower and upper

Sequence of Intraoral Archwires for Conventional Brackets
Archwire kit for Self ligating

Type	Qty	Size	REF
Superelastic - NiTi	18 un*	$\varnothing.012'' / \varnothing.014''$ $\varnothing.016'' / \varnothing.018''$ $.016'' \times .022''$ $.017'' \times .025''$ $.019'' \times .025''$	50.14.900
		.019'' x .025''	
Titânio - Molibdênio - TMA		.017'' x .025''	

*9 different pairs
Each size contains 2 archwires, lower and upper.

Sequence of Intraoral Archwires for Aesthetic Bracket
Archwire kit for Selfligating

Type	Qty	Size	REF
Superelastic Niti Esthetico	12 un*	$\varnothing.012'' / \varnothing.014''$ $\varnothing.016'' / \varnothing.018''$	50.14.902
		.017'' x .025'' .019'' x .025''	

*6 different pairs
Each size contains 2 archwires, lower and upper.

Sequence of Intraoral Archwires for Conventional Brackets

Sequence composed of NiTi, CrNi and TMA archwires, with archwire sizes to be used with conventional brackets.




Sequence of Intraoral Archwires for Aesthetic Brackets

Sequence composed of aesthetic superelastic NiTi archwires, in a sequence of archwire sizes to be used with Aesthetic Brackets.



Extraoral archwires

Morelli extraoral archwires are made of stainless steel, with high strength and malleability. Together with the several types of tray and extraoral appliances, these devices exhibit versatility in the angle of application, allowing the achievement of orthodontic or orthopedic forces. Accompanies exclusive gray bag for archwires.



Extraoral archwires

Type	Size	Qty	REF
Standard	 Ø1,15mm (.045")	1	50.10.001
	 Ø1,30mm (.051")	1	50.11.001


Type	Size	Qty	REF
Medium Hard	 Ø1,15mm (.045")	1	50.10.002
	 Ø1,30mm (.051")	1	50.11.002
Medium Hard Small	 Ø1,15mm (.045")	1	50.10.003
Medium Mild	 Ø1,15mm (.045")	1	50.10.005



Accompanies exclusive grey bag.

Lip-bumper

Morelli lip-bumpers are made of stainless steel and injected plastic, with excellent polishing and surface smoothness. They can be used to reinforce anchorage, move molars and treat lip hypotonia.




REF 1: Lip-bumper without hook


REF 2: Lip-bumper without hook

Lip-bumper

Colors		Size	QTY	REF 1	REF 2
	Blue	Ø 1,15mm – .045"	1	50.40.002	50.50.002
	Pearl	Ø 1,15mm – .045"	1	50.40.006	50.50.006
	Pink	Ø 1,15mm – .045"	1	50.40.007	50.50.007

Wires





Size	Qty	code	Qty	REF
------	-----	------	-----	-----

Morest CrNi wire in stick

	.012" - 0,30mm	10	55.01.012	-	-
	.014" - 0,35mm	10	55.01.014	-	-
	.016" - 0,40mm	10	55.01.016	-	-
	.018" - 0,45mm	10	55.01.018	-	-
	.020" - 0,50mm	10	55.01.020	-	-
	.024" - 0,60mm	10	55.01.024	20	55.06.060
	.028" - 0,70mm	10	55.01.028	20	55.06.070
	.032" - 0,80mm	10	55.01.032	20	55.06.080
	.036" - 0,90mm	10	55.01.036	20	55.06.090
	.040" - 1,00mm	10	55.06.100	-	-
	.047" - 1,20mm	10	55.06.120	-	-
	.016" x .016" 0.40 x 0,40mm	10	55.02.002	-	-
	.016" x .022" 0.40 x 0,55mm	10	55.03.011	-	-
	.017" x .022" 0.43 x 0,55mm	10	55.03.016	-	-
	.017" x .025" 0.43 x 0,63mm	10	55.03.012	-	-
	.018" x .025" 0.45 x 0,63mm	10	55.03.013	-	-
	.019" x .025" 0.48 x 0,63mm	10	55.03.014	-	-
	.021" x .025" 0.53 x 0,63mm	10	55.03.015	-	-

Morest CrNi wire in stick

Stainless steel wire in stick, available in the necessary diameters for all stages of orthodontic treatment

Size	Qty	REF
------	-----	-----

Twist-flex CrNi wire in stick

3 wires		.015" - 0,38mm	10	55.04.315
		.018" - 0,45mm	10	55.04.318
		.020" - 0,50mm	10	55.04.320
4 wires		.015" - 0,38mm	10	55.04.415

Twist-flex CrNi wire in stick

Twist-flex wire is indicated for the first stages of the orthodontic treatment, as well as for making lingual retaining bars.

Size	Qty	REF
------	-----	-----

Twist-flex CrNi wire in stick for lingual bar

3 wires		.032" - 0,80mm	10	55.04.332
---------	---	----------------	----	-----------

#TIP

Also check out the Twistform thread for containment on the page 49

Blue tip morealloy CoCr wire in stick

Cobalt-chromium wires, with strength levels suitable for the bioprogressive technique.

TMA wire in stick

The special titanium molybdenum alloy of Morelli TMA wires brings numerous advantages in the completion of orthodontic treatments due to the characteristics of bending capacity, welding and shape memory, in addition to the wire great elasticity.

TMA brand is sole ownership of Dental Morelli Ltda. in Brazil.

Size	Qty	REF
------	-----	-----


Blue tip morealloy CoCr wire in stick

 .016" x .016" 0.40 x 0,40mm	10	55.02.001
 .016" x .022" 0.40 x 0,55mm	10	55.03.001
 .017" x .025" 0.43 x 0,63mm	10	55.03.002

Size	Qty.	REF
------	------	-----

Alambre varilla TMA

 .012" - 0,30mm	10	56.01.001
 .014" - 0,35mm	10	56.01.002
 .016" - 0,40mm	10	56.01.003
 .018" - 0,45mm	10	56.01.004
 .020" - 0,50mm	10	56.01.005
 .032" - 0,80mm	10	56.01.006
 .036" - 0,90mm	10	56.01.007
 .016" x .016" 0.40 x 0,40mm	10	56.02.001
 .016" x .022" 0.40 x 0,55mm	10	56.03.001
 .017" x .025" 0.43 x 0,63mm	10	56.03.003
 .018" x .025" 0.45 x 0,63mm	10	56.03.004
 .019" x .025" 0.48 x 0,63mm	10	56.03.005
 .021" x .025" 0.53 x 0,63mm	10	56.03.006



Wires in Rolls

Stainless steel wires, with a round cross-section, are used for different purposes in orthodontics, such as in the manufacture of removable appliances. Available in 50g or 500g rolls.

Size	Approx. Size	REF
------	--------------	-----

Hard elastic wire CrNi 50g

 .012" - 0,30mm	90m	55.01.030
 .014" - 0,35mm	66m	55.01.035
 .016" - 0,40mm	50m	55.01.040
 .018" - 0,45mm	40m	55.01.045
 .020" - 0,50mm	32m	55.01.050
 .024" - 0,60mm	22m	55.01.060
 .028" - 0,70mm	16m	55.01.070
 .032" - 0,80mm	12m	55.01.080
 .036" - 0,90mm	10m	55.01.090
 .040" - 1,00mm	08m	55.01.100
 .047" - 1,20mm	05m	55.01.120



Wire in Roll CrNi

Stainless steel wires in roll, available in necessary diameters for all stages of the orthodontic treatment.



Hard elastic wire CrNi 500g

Size	Approx. Size	REF
○ .028" x 0,70mm	165m	55.01.570
○ .032" x 0,80mm	126m	55.01.580
○ .036" x 0,90mm	100m	55.01.590
○ .040" x 1,00mm	81m	55.01.591



Size	Approx. Size	REF
------	--------------	-----

Hard elastic wire CrNi 500g

⌒ ø 1,30 x 0,65mm	5m	55.01.130
⌒ ø 1,50 x 0,75mm	5m	55.01.150
⌒ ø 1,75 x 0,90mm	5m	55.01.175
⌒ ø 2,00 x 1,00mm	5m	55.01.200




Ligature wire CrNi 50g

○ .008" - 0,20mm	201m	55.01.208
○ .010" - 0,25mm	130m	55.01.210
○ .012" - 0,30mm	90m	55.01.212




Aesthetic wire for ligature CrNi in roll

(D1) Internal Wire .010" - 0,25mm	10m	55.01.220
(D2) External Wire .014" - 0,35mm		



Elastics



ELASTICS UNIONS

WITH AN TIM ICRO BIAL AD D ITIVE

Bandage elastics

The new line of joints presents super modern and exclusive colors for modular and canes. And better: all with antimicrobial additives, an innovation that provides more durability and conservation to the elastics.



Bastón



Modular



Cód 75.01.028



Request your color sample here

COR	BASTÓN	Ctd. Anillos	Modular	Ctd. Anillos
-----	--------	--------------	---------	--------------

Elastics for ligation

	Amethyst Purple	60.03.166	1000	60.06.166	500
	Bubblegum Rose	60.03.168	1000	60.06.168	500
	Pitaya Rose	60.03.165	1000	60.06.165	500
	Coral Orange	60.03.167	1000	60.06.167	500
	Lime Green	60.03.161	1000	60.06.161	500
	Caribbean Blue	60.03.162	1000	60.06.162	500
	Sky Blue	60.03.163	1000	60.06.163	500
	Petroleum Blue	60.03.164	1000	60.06.164	500
	Graphite gray	60.03.160	1000	60.06.160	500



ANTIMICROBIAL TECHNOLOGY

The incorporation of the antimicrobial additive to the polymer matrix prevents the development and fixation of fungi and bacteria. This innovation reduces unwanted stains and odors, in addition to preserving the mechanical properties of the elastics.



EXCLUSIVE COLORS

The unions have carefully selected and superior quality pigments, which allow the creation of a more stable, modern and vibrant color palette. There are 9 exclusive options to suit different tastes and styles



SHAPE MEMORY AND MECHANIAL PROPERTIES

Using excellent raw materials with antimicrobial additives, we produce elastics with more suitable mechanical properties. In addition, we heat treat the material, ensuring greater shape memory of the elastics.

PLIC CLAMPS FOR MODULARS

- Tip with perfect fit
- Easy application of the elastic to the bracket
- Also great for applying stops, bows, springs and hooks

Printed colors are approximate and may vary.

Color	Qty Rings	Elastic Ligature in Canes	Modular
-------	--------------	---------------------------------	---------

Elastic ligatures

	Grey	5000	60.03.701	–
	Crystal	5000	60.03.700	–

	Yellow	1000	60.03.305	60.06.104
	Blue	1000	60.03.308	60.06.107
	Baby Blue	1000	60.03.304	60.06.103
	Crystal Blue	1000	60.03.354	60.06.154
	Navy Blue	1000	60.03.316	60.06.115
	White Lace	1000	60.03.327	60.06.127
	Cristal	1000	60.03.300	60.06.100
	Grey	1000	60.03.301	60.06.101
	Crystal Orange	1000	60.03.352	60.06.152
	Lilac	1000	60.03.328	60.06.128
	Ivory	1000	60.03.320	60.06.120
	Crystal Pink	1000	60.03.355	60.06.155
	Silver	1000	60.03.317	60.06.117
	Black	1000	60.03.311	60.06.110
	Pink	1000	60.03.312	60.06.111
	Baby Pink	1000	60.03.319	60.06.119
	Rosa Perla	1000	60.03.326	60.06.126
	Purple	1000	60.03.313	60.06.112
	Miscellaneous	1000	60.03.302	60.06.118
	Crystal Mix	1000	60.03.358	60.06.158
	Crystal Grape	1000	60.03.359	60.06.159
	Crystal Green	1000	60.03.353	60.06.153
	Crystal Lime Green	1000	60.03.350	60.06.150
	Moss Green	1000	60.03.307	60.06.106
	Sea Green	1000	60.03.315	60.06.114
	Pearls Green	1000	60.03.324	60.06.124
	Red	1000	60.03.310	60.06.109
	Crystal Red	1000	60.03.351	60.06.151
	Wine	1000	60.03.322	60.06.122

The printed colors are approximate, subject to change

Elastic ligatures

Morelli elastic ligatures use polyurethane controlled in its own manufacturing process and non-toxic pigmentation. Featuring different colors in canes and modular versions. These elastics keep efficient levels of memory and low rates of deformation in the oral environment.

Elastic Ligature in Canes






Modular



Request your color
sample here

Elastic chain

Elastic chains are useful in space closure mechanics, anterior retraction, and tooth traction. They are made from high quality material, offering elasticity and memory. Available in short, medium and long links versions, colored and in 4.5m packages (gray and crystal) and 1.5m packages (other colors).



	Color	Length	Short	Medium	Long
Elastic chain					
	Grey	4,5m	60.05.501	60.05.511	60.05.521
	Crystal	4,5m	60.05.500	60.05.510	60.05.520
	Yellow	1,5m	60.05.105	60.05.205	60.05.305
	Blue	1,5m	60.05.108	60.05.208	60.05.308
	Baby Blue	1,5m	60.05.104	60.05.204	60.05.304
	Crystal Blue	1,5m	60.05.154	60.05.254	60.05.354
	Navy Blue	1,5m	60.05.116	60.05.216	60.05.316
	White Lace	1,5m	60.05.127	60.05.227	60.05.327
	Crystal Orange	1,5m	60.05.152	60.05.252	60.05.352
	Lilac	1,5m	60.05.128	60.05.228	60.05.328
	Ivory	1,5m	60.05.120	60.05.220	60.05.320
	Crystal Pink	1,5m	60.05.155	60.05.255	60.05.355
	Silver	1,5m	60.05.117	60.05.217	60.05.317
	Black	1,5m	60.05.111	60.05.211	60.05.311
	Pink	1,5m	60.05.112	60.05.212	60.05.312
	Baby Pink	1,5m	60.05.119	60.05.219	60.05.319
	Pearl Pink	1,5m	60.05.126	60.05.226	60.05.326
	Purple	1,5m	60.05.113	60.05.213	60.05.313
	Crystal Grape	1,5m	60.05.159	60.05.259	60.05.359
	Crystal Green	1,5m	60.05.153	60.05.253	60.05.353
	Crystal Lime	1,5m	60.05.150	60.05.250	60.05.350
	Moss Green	1,5m	60.05.107	60.05.207	60.05.307
	Sea Green	1,5m	60.05.115	60.05.215	60.05.315
	Pearls Green	1,5m	60.05.124	60.05.224	60.05.324
	Red	1,5m	60.05.110	60.05.210	60.05.310
	Crystal Red	1,5m	60.05.151	60.05.251	60.05.351
	Wine	1,5m	60.05.122	60.05.222	60.05.322

The printed colors are approximate, subject to change

ELÁSTICOS INTRAORAIS

Black LATEX

The Morelli Intraoral Elastic Line now has the **Black Latex** version. With an antimicrobial additive, this is an innovation that brings even more aesthetics and durability to elastics.



ANTIMICROBIAL TECHNOLOGY

The incorporation of the antimicrobial additive into the elastomer matrix prevents the development and attachment of fungi and bacteria, preserving the mechanical properties of the elastics, in addition to ensuring the maintenance of the **Black Latex** color.



MORE DISCRETION

The lighter shade of natural latex stands out in the patient's oral cavity. Black latex gives more discretion when opening the mouth, which can help with acceptance during



Black
LATEX



LATEX



WITH MINT
SCENT

	CODE	QTY
● Medium 1/8" = 3,20mm	60.02.410	1000
● Medium 3/16" = 4,80mm	60.02.411	1000
○ Medium 1/4" = 6,35mm	60.02.412	1000
○ Medium 5/16" = 7,94mm	60.02.413	1000



Color	Ctd.	Modular
-------	------	---------

Elastic ligatures

Pearl Blue	500	60.06.010
------------	-----	-----------

Elastic Pearl Blue

Morelli Pearl Blue Elastics They are ideal for all orthodontic treatments, and especially indicated for treatment with aesthetic braces. Due to its high resistance to stains and bluish color, it has an excellent aesthetic result!



Color	Length	Short	Medium	Long
-------	--------	-------	--------	------

Elástico en cadena

Pearl Blue	1,5m	60.05.101	60.05.201	60.05.301
------------	------	-----------	-----------	-----------

Force	Color	Size	Qty	REF
-------	-------	------	-----	-----

Latex intraoral elastics 3/4"

Ligth	Natural	3/4" = 19,05mm	100	60.01.330
Medium		3/4" = 19,05mm	100	60.01.340

Latex intraoral elastics

The natural latex elastic properties are widely used by Orthodontics through dental intercuspation of several levels. They have a mint scent

Force	Color	Size	Qty	REF
-------	-------	------	-----	-----

Latex intraoral elastics

Ligth	Natural	1/8" = 3,20mm	1000	60.01.200
		3/16" = 4,80mm	1000	60.01.201
		1/4" = 6,35mm	1000	60.01.302
		5/16" = 7,94mm	1000	60.01.303
		3/8" = 9,53mm	1000	60.01.304
Medium	Natural	1/8" = 3,20mm	1000	60.01.310
		3/16" = 4,80mm	1000	60.01.311
		1/4" = 6,35mm	1000	60.01.203
		5/16" = 7,94mm	1000	60.01.205
		3/8" = 9,53mm	1000	60.01.314
Heavy	Natural	1/8" = 3,20mm	1000	60.01.320
		3/16" = 4,80mm	1000	60.01.321
		1/4" = 6,35mm	1000	60.01.322
		5/16" = 7,94mm	1000	60.01.323
		3/8" = 9,53mm	1000	60.01.208



Clinical protocols and tips to use intraoral elastics.



Force	Color	Size	Qty	REF
-------	-------	------	-----	-----

Latex intraoral elastics

Ligth		 1/2" = 12,7mm	100	60.01.210
Medium	Natural	 1/2" = 12,7mm	100	60.01.211
Heavy		 1/2" = 12,7mm	100	60.01.212



Silicone intraoral elastics

Morelli Silicone intraoral elastics were developed to meet the mechanical needs required by brackets and accessories with an aesthetic, biocompatible and anti-allergic appeal. They have a mint scent.

Force	Color	Size	Qty	REF
-------	-------	------	-----	-----

Silicone intraoral elastics

Medium	Translucent	 1/8" = 3,20mm	1000	60.08.110
		 3/16" = 4,80mm	1000	60.08.111
		 1/4" = 6,35mm	1000	60.08.112
		 5/16" = 7,94mm	1000	60.08.113

The nominal force will stabilize only after the first stretch



Separator elastic

Morelli separator elastics meet the elasticity and resistance specifications, promoting an efficient opening of interproximal space.

Force	Color	Size	Qty	REF
-------	-------	------	-----	-----

Separator elastics

In bulk	 Blue	 5/32" = 4,0mm	1000	60.04.200
	 Green	 3/16" = 4,80mm	1000	60.04.400



Modular	 Blue	 5/32" = 4,0mm	500	60.04.201
	 Green	 3/16" = 4,80mm	500	60.04.401





Color	Size	Diam.	REF
-------	------	-------	-----

Tube for protection

	Crystal		0.75mm	60.05.400
	Grey		0.75mm	60.05.401
	Blue		0.75mm	60.05.403
	Crystal Pink		0.75mm	60.05.405
	Crystal		0.95mm	60.05.410
	Grey		0.95mm	60.05.411
	Blue		0.95mm	60.05.402
	Crystal Pink		0.95mm	60.05.404



*Diameters: 0.75mm for .012" up to .016"x.022" wires
0.95mm for .017" up to .021"x.021" wires

Tube for protection

This tube function is to protect the mucosa and the wires themselves, providing reinforcement and comfort to the patient. In addition to protection, the color tubes provide a decorative effect and the possibility of combining colors with other ligatures used in orthodontic treatment.

Promoting the aesthetic appeal the patient desires in a safe way.



Color	Qty	REF
-------	-----	-----

Rotator elastics

	Cristal	100	60.07.100
	Cinza	100	60.07.101



The rotator elastics are positioned on the brackets fins, in order to promote the rotational movement of the tooth.



Universal rotator elastics

	Cristal	100	60.07.200
	Cristal	100	60.07.201





Consulta
cómo usarlo.

Functional Orthopedics






Size	Opening	Qty.	REF
------	---------	------	-----

Expansor standard

	 Expansion screw 5 Opening 7 + 0.8mm	7mm	10	65.05.001
		7mm	1	65.05.101
	 Expansion screw 9 Opening 11 + 0.8mm	11mm	10	65.05.009
		11mm	1	65.05.109

Los expansores de Morelli se utilizan en aparatos fijos o móviles para la expansión o disyunción de las arcadas dentarias. Están disponibles en varios tamaños, modelos y aperturas.

Expansor universal

	 Expansion screw 5 Opening 7 + 0.8mm	7mm	10	65.05.002
		7mm	1	65.05.102
	 Expansion screw 7 Opening 9 + 0.8mm	9mm	10	65.05.003
		9mm	1	65.05.103
	 Expansion screw 9 Opening 11 + 0.8mm	11mm	10	65.05.008
		11mm	1	65.05.108

*Expansion
All expansion screws come with an activation key code 75.02.026 - page 103



Size	Opening	Qty	REF
------	---------	-----	-----

Mini expansion screw

	6,5mm	10	65.05.006
	6,5mm	1	65.05.106



	9mm	10	65.05.007
	9mm	1	65.05.107



Unilateral mini expansion screw

	5mm	1	65.05.156
--	-----	---	-----------



Unilateral expansion screw

	6.5mm	1	65.05.151
--	-------	---	-----------



Fan type expansion screw

	8mm	1	65.05.160
--	-----	---	-----------




*Expansion
Opening

Size	Opening	Qty	REF
------	---------	-----	-----

Expansor Hyrax


Palatal expansion screw

	7mm	1	65.05.010
	9mm	1	65.05.011
	11mm	1	65.05.012
	13mm	1	65.05.013

Hyrax Type
Hyrax type expansion screws are indicated to correct poor maxillomandibular relations. The most common clinical situations to apply are:

- Mixed/early permanent dentition cases;
- Posterosuperior bilateral crossbite, as a result of deficiency of the maxilla apical base;
- Anterior or posterior crossbites that are not too severe.

Rapid fan expansion screw

	10mm	1	65.05.014
--	------	---	-----------

All Expansion Screw are accompanied by an activation key code 75.01.038 – Page 103


Tipo Leque

Hyrax Type The rapid fast fan expansion screw is intended for the anterior expansion of the upper arch, mounted in the anterior palatal region, commonly used in the treatment of patients with craniofacial anomaly (cleft palate).

Size	Qty	REF
------	-----	-----


Escudo de Bimler

Bimler shield

	10	65.01.010
---	----	-----------

Aparato que ensambla y une los arcos inferiores del modelador elástico Bimler.
El escudo de Bimler se utiliza en casos de distocclusión con o sin apiñamiento y falta de espacio para el canino

Equiplan planas

	10	65.01.040
---	----	-----------

ortho protector

Multipurpose Lift Plate



The Ortho Protector is a bite lift plate that acts as a muscle relaxant.

It makes it possible to **temporarily** treat occlusal stress problems, such as bruxism and temporomandibular disorders. In a practical and comfortable way, even in patients with orthodontic braces.

In addition to assisting in the protection of laminated facets and rehabilitation after implant surgery.



No moldings and with adaptation instant on any patient

Resistant Material

Made of polyurethane, which provides durability, even with occlusal tightening.

Occlusal Splints

They easily adjust to all types and sizes of arches, including children aged 5 years and older.



Size Adjustments

There are 4 positions available for better adaptation and comfort during use.



Front Handle

Fits anatomically on the buccal side of the lower incisors.


Ortho Protector

Ortho Protector is an accessory developed for daily use to guarantee the patient, even with bruxism, the night protection of dental crowns and joints.



Size	Opening	Qty	REF
 Adjustable 	Blue	1	65.03.001
	Rose	1	65.03.002

Comes with box for storage and transportation; and instructions for use for the patient



It is possible to use occlusal plates in 4 position for better adjustment and comfort.



Watch the video with the full explanation on the Ortho Protector use. Point your cell phone camera to scan the QR code on the side.



Headgear For Safety Breakaway devices





Size	Qty	REF
------	-----	-----

Headgear for Safety Breakaway

Small	1	70.31.001
Medium	1	70.31.002
Large	1	70.31.003

Headgear

They are used for extraoral anchorage, allowing the application of forces with



Size	Qty	REF
------	-----	-----

Headgear - Adjustable traction

one size	1	70.35.001
----------	---	-----------



Size	Qty	REF
------	-----	-----

Chin cup complete - high traction

one size	1	70.36.001
----------	---	-----------

Mentonera

Para contener el crecimiento mandibular, aplicando fuerzas en dirección parietal. Se proporciona completa o solamente con el apoyo del mento



Chin cup single

one size	1	70.36.002
----------	---	-----------


Cervical flap

For application of extra-oral force with cervical traction, using breakaway.



Size	Qty	REF
------	-----	-----

Universal cervical flap for low traction

	1	70.20.032
---	---	-----------

Petit Facial Mask

Anchorage device used for reverse traction of the maxilla with posteroanterior advancement. Morelli petit facial mask promotes at the same time, a forward tilting action down and back of the jaw.

Color	Qty	REF
-------	-----	-----

Petit facial mask

Blue	1	70.50.001
Rose	1	70.50.002
Gris	1	70.50.003



The petit Morelli face mask can be used in two ways: fixed chin guard or articulated chin guard.



Articulated Chin Cup



Fixed Chin Cup



Qty	Light Force 350gf	Medium Force 500gf	Heavy Force 750gf
-----	----------------------	-----------------------	----------------------

Breakaway

It unlocks the clip spontaneously if the force applied exceeds the limit.

Breakaway with traction device

2	70.42.102	70.42.202	70.42.302
---	-----------	-----------	-----------

Confira como
utilizar



Size	Qty	REF
------	-----	-----

Parietal elastic double clasp - low traction

One size	1	70.20.011
----------	---	-----------



Instruments





Applicator Tip

Extractor Tip

Size	Qty	REF
------	-----	-----

Elastic ligature applicator/extractor

	1	75.01.001
--	---	-----------

	1	75.01.002
--	---	-----------

Ligature applicator

This instrument is used to remove and place elastic ligatures.

They have special stainless steel double active tips.



Check how to open Morelli Self-ligating clip.



Elastic Ligature Applicator

	1	75.01.003
--	---	-----------

	1	75.01.004
--	---	-----------

Bracket positioner

Made of stainless steel, Morelli bracket positioner features an exclusive articulated tip, which facilitates positioning in posterior teeth.

This instrument allows the correct positioning of the bracket slot in relation to the occlusal or incisal tooth.







Bracket positioner

	3,5mm	1	75.01.006
	4,0mm		
	4,5mm		
	5,0mm		
	2,5mm		
	3,0mm	1	75.01.039
	5,5mm		
	6,0mm		

Band seater

Indicated to accommodate the band on the tooth.

Tips made of stainless steel with formats that favor the adaptation of orthodontic bands. Produced in gray nylon, they resist the bite pressure and the sterilization temperatures in an autoclave.

Size	Qty	Format	REF
Band seater			
	1	Triangular tips	75.01.011
	1	Invertible triangular tip	75.01.012
	1	Round tip	75.01.013
	1	Punta cuadrada	75.01.014



Lip retractor

Morelli lip retractors are useful for optimizing isolation during direct bonding of accessories, as well as taking intraoral pictures.

Size	Colors	Qty	REF
Retractor labial			
Small	Amber	2	75.01.033
Large	Amber	2	75.01.043
Small	Amber	2	75.01.036
Large	Amber	2	75.01.046
Small	White	1	75.01.133
Large	White	1	75.01.143
Small	White	1	75.01.136
Large	White	1	75.01.146





Material	Colors	Size	REF	Qty
----------	--------	------	-----	-----

Impression trays

Impression trays

Plastic	Grey	01	75.01.051	1 pair
	Yellow	02	75.01.052	
	Pink	03	75.01.053	
	Orange	04	75.01.054	
	Green	05	75.01.055	
	Red	06	75.01.056	
	Blue	07	75.01.057	
	Purple	08	75.01.058	

Morelli impression trays have excellent alginate retention due to the design of the grooves on the base, in addition to allowing a complete impression of the dental arch, due to the larger lateral height.

They can be autoclaved or chemically sterilized.



Tower to form arches

The tower is an instrument used to form orthodontic arches from square or rectangular wire rods of several thicknesses. Morelli Tower is made of stainless steel with an engineering polymer handle featuring excellent finish and resistance to autoclaves.

Size	Qty	REF
.016"	1	75.02.005
.017"		
.018"		
.019"		
.021"		



Detacher

It removes the cap from single, double and triple convertible tubes. Its fitting also allows the use for scalpel blades. The detacher blades are available for replacement.

Size	Qty	REF
130mm	1	75.01.024

Detacher blade

33mm	5	75.01.025
------	---	-----------




Tensiometer

Mechanical tensiometer used to measure forces in intraoral and extraoral elastics. They have two tips for easy positioning.

COLORS	Force	Qty	REF
	100gf to 1600gf	1	75.02.006
	50gf to 500gf (cN)	1	75.02.008
	25gf to 250gf (cN)	1	75.02.009





Size	Qty	REF
------	-----	-----

Orthodontic tweezer



1

75.01.022


Orthodontic tweezer

The orthodontic tweezer was designed to allow the safe and perfect gripping of brackets and accessories.

The handle end has a unique finish to facilitate the bracket placement and positioning.



Short tweezer for brackets position



1

75.01.064



Pinza para tuboPinza para tubo




1


75.01.027

Pinza para tubo

La pinza de posicionamiento de tubos se desarrolló para ayudar en el posicionamiento durante la unión de tubos molares y accesorios en lugares de difícil acceso.



Pinza "Plic"



1

75.01.028

Pinza Plic

La pinza "Plic", a través de la punta empotrada, ofrece probadas ventajas en la colocación de elásticos, resortes, topes redondos, etc.



Sugerencia de aplicación.


Bands marker

This instrument is used to mark accessory heights on orthodontic bands.

It has a swivel mechanism with sizes from 3.5mm to 5.0mm.

Size	Qty	REF
------	-----	-----

This instrument is used to mark

	3,5mm	1	75.01.007
	4,0mm		
	4,5mm		
	5,0mm		





Bands marker

Made of stainless steel, this instrument has heat-treated tips. It is used to mark bands and position brackets and accessories.

Size	Qty	REF
------	-----	-----

Boone gauge

	3,5mm	1	75.01.023
	4,0mm		
	4,5mm		
	5,0mm		
	2,0mm	1	75.01.029
	2,5mm		
	3,0mm		
	3,5mm		
	4,0mm	1	75.01.030
	4,5mm		
	5,0mm		
	5,5mm		




Bender

Morelli bender is a convenient device to facilitate bending in thick wires (Hyrax type Expansion Rod, etc.). It has openings and slots for fitting wires, allowing the use of wires up to Ø1.5mm thick.

Size	Qty	REF
------	-----	-----

Wire and rods bender

	1	75.02.001
---	---	-----------



Distal arm bending

Used to bend and finish the orthodontic wire distally, avoiding discomfort to the patient. Available in versions with a hole (closed) or with a slot (open) in addition to the appropriate curvature to facilitate its application.

Size	Qty	REF
------	-----	-----

Distal arm bending - closed slot

.022"	1	75.02.022
-------	---	-----------



Distal arm bending - opened slot

.022"	1	75.02.023
-------	---	-----------

Tienen puntas rectas y curvas para un mejor ajuste y aplicación.



Curve

Size	Qty	REF
------	-----	-----

Arm bending for wire ligature - curve

.022"	1	75.02.074
-------	---	-----------

Arm bending for wire ligature

The arm bending for ligature are useful for bending and accommodating metallic and elastic ligatures in the canals of brackets and other orthodontic appliances.



Straight

Arm bending for wire ligature - straight/curve tip

.022"	1	75.02.075
-------	---	-----------



Hexagonal key

Hexagonal key for fixing Gurin nut (page 46).

Size	Qty	REF
------	-----	-----

Hexagonal key for gurin

	1	75.02.011
---	---	-----------



Elastic positioner

Positioners (applicators) facilitate theplacement of elastics made by thepatient, according to the professional guidance.

Color	Qty	REF
-------	-----	-----

Posicionador de elásticos

Black	20	75.01.015
-------	----	-----------

See retainer box for elastic n°3 on page 103



Template

Indicated for making orthodontic cephalometric tracings. It allows the professional to get a better view of the patient problems, which provides a closer to ideal diagnosis

Size	Qty	REF
------	-----	-----

Benvenega cephalometric

190 x 200mm	1	75.01.020
-------------	---	-----------




Expansion screw key

It allows activating the expansion screw safely and conveniently.

Long rod	5	75.01.038
----------	---	-----------

Short rod	10	75.02.026
-----------	----	-----------



**Diameter****Longitud****Qty****REF**

Zirconia bur

Zirconia low rotation bur

Ø 3 mm

23 mm

5

75.03.001

Indicated for removing orthodontic resin safely, without damaging or scratching the dental enamel.

Use in low rotation contraangle, removing on average up to 5 arcades.

Efficiency

Morelli bur was developed specially to remove the remaining resin after the appliance removal.

Reusable

High strength to autoclave sterilization cycles and easy cleaning of the grooves, as they are straight.

Great finish

The bur quality causes low vibration, which generates less heat and perfectly smooth surfaces.

Differentiated material

The only removal bur on the market injected entirely in zirconia.

Durability

The Zirconia Bur has excellent cutting durability, strength, does not damage or contaminate the teeth. It is equally effective in any rotation direction, which prevents errors and increases its life.

Safety

Its multi-laminated active tip does not hurt the gum, which provides a controlled resin removal, eliminating the risk of cracking, scratching the enamel or damaging the soft tissue.



Check the step by step, from removal to polishing.

Disposable point for finishing

Ø 3 mm

23 mm

20

75.03.002

Disposable point

Excellent finishing

The point for finishing provides a great final finish of the enamel.

Safety

Its rounded format is safer for the gum and allows access to areas close to it, eliminating any remaining resin and ensuring an excellent finishing.

Material quality


Its polymer composition with fiberglass has adequate abrasiveness for the removal of adhesive residues, ensuring an effective and safe finishing against damage to dental enamel.

Convenience

Morelli tip can be applied dry and in any rotation direction. As it is disposable, it does not require sterilization.

Ideal to recover the brightness and smoothness of the dental enamel after using Morelli Zirconia Bur. Use in low rotation contra-angle, removing on average up to 1 arcade.

ORAL HYGIENE



Escova Dental BioFlex®

Bioflex® toothbrushes have arrived bringing a different, comfortable and safe way to take care of your oral health. Developed with a different design from conventional toothbrushes, Bioflex® toothbrushes have 393 "rubberized" bristles that provide excellent removal of bacterial plaque, protecting gums and mucous membranes from injuries. Because they have soft bristles, they also preserve restorations, orthodontic appliances, implants and other oral structures.

Sequence composed of NiTi, CrNi and TMA archwires, with archwire sizes to be used with conventional brackets.

High technology applied in production

Its great differential is due to the injected "monobloc" bristles which, unlike traditional "tufted" bristles, have fewer recesses and therefore favor the hygiene of the brush, preventing the accumulation of residue and the proliferation of bacteria.



Healthy and complete mouth cleaning

On the back of the bristles, the brush also has a "tongue scraper", which allows you to remove the tongue "coating".

With innovative design and differentiated features

The Bioflex® line was developed to meet the needs of all patients.

With innovative design and differentiated features

The shape of the handle provides a good fit in the hand, making brushing easier.

Bioflex Toothbrush

Bioflex® Brushes came bringing a different, comfortable and safe way to take care of your mouth health. Developed with a different design from conventional toothbrushes, Bioflex® brushes have 393 “rubbery” bristles that provide excellent bacterial plaque removal, protecting gums and mucous membranes from injuries. Because they have soft bristles, they also preserve restorations, orthodontic appliances, implants and other oral structures.

Color	Qty	REF
Bioflex Toothbrush		
White/Pink	01 unit	78.01.200
Blue/Lemon	01 unit	78.01.201
Pink/Purple	01 unit	78.01.202
White/Blue	01 unit	78.01.203



Display

Con 10 pinceles surtidos, el Display es ideal para exhibirlo en la recepción de la oficina.

Color	Qty	REF
Bioflex Toothbrush Soft Assorted on Display		
Assorted	1 kit 10 brushes	78.01.210

Kit composed of:
2 brushes (White/Pink), 3 brushes, (Blue/Lemon), 3 brushes, (Pink/Purple), 2 brushes, (White/Blue).





Color	Qty	REF
-------	-----	-----

Interdental brush

White	1 Kit	78.01.001
Blue	1 Kit	78.01.002
Yellow	1 Kit	78.01.003
Pink	1 Kit	78.01.004

Size	Qty	REF
------	-----	-----

Refills of interdental brush

Cylindrical P	1 blister 6 units	78.01.101
Cylindrical M	1 blister 6 units	78.01.102
Cylindrical G	1 blister 6 units	78.01.103



Color	Qty	REF
-------	-----	-----

Tongue cleaner

Blue and Pink	1 pair 2 units	78.02.001
Blue and Yellow	1 pair 2 units	78.02.002

Interdental brush

Interdental brushes are used for cleaning between teeth (interproximal), preventing tartar and periodontal diseases. They are very efficient for cleaning orthodontic appliances and reach places with difficult access for other cleaning instruments. Morelli interdental brushes help prophylaxis, favoring the evolution of orthodontic treatment, preventing stains and cavities in the teeth as well as gingival infections.



white coated wire: better bristle fixation and patient comfort

Kit composed of:
3 brushes, 1 case and 1 cable.

Tongue cleaner

Morelli tongue cleaner has anatomical and functional design, suitable for removing tongue coating (responsible for bad breath and bacterial plaque formation), with maximum efficiency and care, without damaging the tongue papillae.

Compartments



Box

In addition to being beautiful, Morelli Pocket Boxes are super resistant, convenient and hygienic. They have a vent on the top of the cap, which ensures air circulation, keeping both the box and its contents always dry. The Pocket Box is indicated for transport and storage of several items, such as: Retaining plate, bruxism plate, whitening plate, mobile appliances and aligners.

Colors	Qty	REF
Pocket box for mobile appliances and elastics		
Blue	10	80.04.211
Pink	10	80.04.212
Off white	05	80.04.214
Black Piano	05	80.04.215
Turquoise Blue	05	80.04.216
Navy Blue	05	80.04.217
Cherry Red	05	80.04.218
Magenta Rose	05	80.04.219
Tangerine Orange	05	80.04.220
Lime Green	05	80.04.221



Colors	Qty	Box	REF
--------	-----	-----	-----

 Assorted	10	Removable Appliances	80.03.004
--	----	----------------------	-----------



 Assorted	10	Removable Appliances	80.04.001
--	----	----------------------	-----------



 Assorted	10	For Elastics	80.04.002
--	----	--------------	-----------


Accompany Device for place and remove elastics!

Elastic positioner single code 75.01.015 - page 92



Size	Qty	Code
------	-----	------

Zip bag for elastics

	50	80.05.001
---	----	-----------





Size	Use	Qty	REF
------	-----	-----	-----

Box for upper bands

	Universal	1	80.02.011
--	-----------	---	-----------

Organizer box

With a rounded bottom and numbered bands and brackets, Morelli boxes are ideal for storing accessories, providing convenience and organization to the office.

They are made of high impact polystyrene.



Box for lower bands


	Universal	1	80.02.021
--	-----------	---	-----------



Size	Qty	Code
------	-----	------

Box for brackets

	1	80.02.030
--	---	-----------



Multiuser organizer box

	1	80.06.020
--	---	-----------

Accompanied by 64 adhesive labels, which can be glued in the Box to identify the materials..

Adhesivos



ORTHOBONDPLUS[®] + ORTHOPRIMER

Color Change Adhesive

Made by those who **know most** about brackets

for those who **understand more** about orthodontics

The **ORTHOBOND** it presents **HIGH PERFORMANCE** for bonding metallic, ceramic and composites, tubes and other accessory.

Yields more than **15 CASES COMPLETE**

The **ORTHOPRIMER** is a fluid adhesive used as a base for activating photopolymerizable orthodontic adhesives. It has high fluidity for micro penetration into enamel prisms, enhancing bonding in brackets and accessories.



Cód. 85.01.017



Cód. 85.01.016

The use of **ORTHOBONDPLUS** in conjunction with **ORTHOPRIMER** it provides a considerable increase in **EFFICIENCY OF COLLAGES**.

Yields more than **50 CASES COMPLETE**

SYSTEM COLOR CHANGE

When carrying out light curing, the adhesive loses its pink tone and has a color similar to that of the tooth.




Adhesive
Orthobond Plus

Size	QTY.	REF
------	------	-----

Orthobond Plus Color Change Adhesive

4 g	1	85.01.017
-----	---	-----------



Adhesivo
OrthoPrimer

Size	QTY.	REF
------	------	-----

OrthoPrimer

5 g	1	85.01.016
-----	---	-----------



Consult essential union advice in our
Orthodontic Union Manual.
Request through our website:
www.morelli.com.br/loja/catalogobrindes

Point your cell phone camera
to watch our collage talk



OPENBITE

Bite Lift Resin

**More agility in disocclusion,
and safety in removal**



Cód. 85.01.018

IDENTIFICATION SAFE
OCCLUSAL LIFT

Bite Lift Resin Open Bite

Size	QTY.	REF
Bite Lift Resin - Open Bite		
4 g	1	85.01.018



IDENTIFICATION SAFE

Formula with permanent blue pigment, which allows easy visualization from application to removal of occlusal lifts. Avoiding unnecessary wear on tooth enamel.

RESISTANCE

Its high adhesion and durability provides excellent mechanical and adhesive resistance. Avoiding wear on chewing.

LIGHT-CURING

Maximum resin curing control, to facilitate modeling and positioning of the lift.

SELF-ADHESIVE

This recipe has good adhesion and can some cases, providing more agility to the clinical stage.

VISCOSITY

Open Bite has the ideal consistency, whats makes it versatile to build the bite lifts.

MANEQUIM MORELLI

Morelli offers you, at cost price, quality material to help you in your daily clinical routine!



Cód: 91.99.001

The Morelli Mannequin comes with an exclusive box with 3D, high-resolution images of each bracket, making it easier to demonstrate the differences between each piece to your patient.

- ¼ Mouth Self-Ligating SLI
- ¼ PolySapphire mouthpiece
- ¼ Ceramic Mouth
- ¼ Max Mouth
- 01 Mini screw
- 02 NiTi .010" Archwires
- 04 Single Roth Tubes

Your patient will feel more secure knowing up close the pieces they will use.

In addition to being an excellent sales tool!

Olá Doutor!

QUEREMOS TE AJUDAR

ATRAIR MAIS PACIENTES

A Morelli sempre prezou pela saúde e pela importância do sorriso.

O **Encontre um Dentista** é uma iniciativa da Morelli para aproximar você, dentista, do seu paciente. Você pode participar com seu perfil pessoal, ou com o perfil de sua clínica, sem qualquer custo.



- ✓ **100% Gratuito**
- ✓ **Disponível para todo o Brasil**
- ✓ **Versão mobile e desktop**
- ✓ **Seja encontrado por mais pacientes próximos a você!**

QUER SER ENCONTRADO?



Aponte a câmera do seu celular e Faça seu cadastro

CONHEÇA A FERRAMENTA



Aponte a câmera do seu celular ou acesse: www.encontreumdentista.com.br

Para realizar o cadastro esteja logado na loja, acesse Meus Dados > Encontre um Dentista > Dados Cadastrais.

#JUNTOSPELOBRASIL
#JUNTOSPELAORTODONTIA

Uma iniciativa: **MORELLI**



GESTÃO 360° para seu consultório



**TOTALMENTE ONLINE
E GRATUITO**



Gestão Financeira

Produtividade gera necessidade de controle financeiro, assim você possui todo o controle de suas receitas e despesas.



Acesso remoto

Acesso de qualquer lugar a qualquer momento e com qualquer dispositivo ligado a internet.



Aplicativo BlueDental

O App está disponível para Android e Iphone para consulta de agendamentos e leitura de código de barras, para entrada e saída de materiais.

**Um Gerenciador de Clínicas completo e dedicado
para odontologia. Agora disponível para você!**

EVOLUA E DIGITALIZE com a gente



Gestão para sua Clínica Odontológica de maneira Simples e Acessível

Confira as vantagens do Bluedental e entenda
como podemos gerir seu consultório



Agendamentos

Faça agendamentos nos horários que você atende, mantendo encaixes em uma visualização complementar.



Gestão de Estoque

Atualize rapidamente seu estoque com a leitura por código de barras pelo celular e o upload de listas.



Envio de mensagens

Através do e-mail ou WhatsApp é possível confirmar, consultar, parabenizar o paciente pelo seu aniversário e agendar retorno.



Odontograma

Crie odontogramas personalizados para cada paciente. Inclua procedimentos por dente, faces e arcadas.



Geração de orçamento

Direto no odontograma é possível calcular o valor do tratamento e imprimir, salvar ou enviar orçamentos.



Automação de documentos

Recibos, contratos, odontogramas, fluxos de pagamentos, prontuários, tudo isso gerado em poucos cliques.

Seus Dados Seguros



Além de termos uma rotina para backup dos seus dados, eles são criptografados, isso permite ter a certeza que ninguém mais terá acesso aos seus dados.



Com foco na qualidade e eficiência contínua, a Morelli tem orgulho em avançar com a incorporação da BlueDental, um software para gestão de clínicas, que oferece uma tecnologia acessível a todos os dentistas.

ACESSE E UTILIZE GRATUITAMENTE



FACTORY VISIT



- **+ 17 years of program**
- **+ 500 visits made**
- **+ 12,000 visitors**



**AN AMAZING E
XPERIENCE**



**IN THE LARGEST
ORTHODONTICS CENTER
IN SOUTH AMERICA**



**SPECIALIZED IN
ORTHODONTIC
PRODUCTS**

**COME AND MEET
MORELLI TOO**



SCHEDULE YOUR VISIT
VISIT: WWW.MORELLI.COM.BR
INFORMATION: 0800-0121455

WHERE TO BUY OUR PRODUCTS

DISTRIBUTORS

We have a network of more than 200 distributors throughout the country.

Use the QR CODE to find the one closest to you.



ONLINE STORE

Through the Morelli website, we keep our portfolio available, with more than 1,700 items, in addition to other amenities:



**FREE SHIPPING ON
ORDERS OVER
R\$199.00**



**DISCOUNT FOR
PAYMENT IN CASH**



**INSTALLME
NTS IN UP**

OUR SERVICE CHANNELS



(15) 3238 8200

WHATSAPP



0800 012 1455

0800 703 1455



export@morelli.com.br

Opening hours: 8:00 a.m. to 5:48 p.m., Monday through Friday (Brasília time),
excluding holidays

We have qualified technical support to resolve all your questions.



☎ (15) 3238.8200 WHATSAPP
☎ 0800 0121455 • 0800 7031455
✉ VENDAS@MORELLI.COM.BR
🖱 WWW.MORELLI.COM.BR

📷 MORELLIORTODONTIA
f MORELLI.ORTODONTIA
▶ DENTALMORELLI
in DENTAL-MORELLI-LTDA