



ORTHODONTIC INSTRUMENTS

Pliers & Instruments Catalog

2026 EDITION

Precision orthodontic pliers, cutters and instruments — with live availability: **In Stock** and **Pre-Order**.



★ CERTIFIED QUALITY

✓ ISO 9001:2008

✓ ISO 13485:2016

✓ CE CA006210

orthodentalusa.com

Online store & VIP pricing

+1 (954) 655-8487

Call or WhatsApp

sales@orthodentalusa.com

1500 Weston Rd Suite 204, Weston, FL 33326, USA

WELCOME

Instruments built for the orthodontic chairside.

Orthodontalusa Corp supplies orthodontic practitioners, clinics and laboratories with professional-grade supplies from our base in Weston, Florida. This catalog presents our full range of orthodontic pliers, cutters and chairside instruments — and tells you at a glance what is **in stock today** and what you can **pre-order**.

★ CERTIFIED QUALITY

✓ ISO 9001:2008

✓ ISO 13485:2016

✓ CE CA006210

426

models in this catalog

237

ready to ship from stock

189

available to pre-order

Craftsmanship you can feel

Medical-grade stainless steel

01

Pliers and instruments are made from AISI 420 medical-grade stainless steel for strength and long service life.

Hardened cutting tips

02

Extra-hard, durable Immunity-HS steel tip inserts and tungsten-carbide inserts on distal end cutters, ligature cutters and hard-wire cutters, joined with special welding for extreme stability.

Precision joints

03

High-precision joint assembly keeps tips aligned and gives frictionless, smooth action. Rivet, box-joint and hex-screw designs with diamond-honed cutting edges.

Ergonomic & safe

04

All sharp edges are rounded for safer, easier handling. Every plier is ultrasonically cleaned, lubricated and tested with its recommended wire sizes.

Autoclave-ready

05

Ideal for steam-autoclave sterilization. See the care & sterilization guide at the back of this catalog.

Local support

06

Real people in Weston, FL — call, WhatsApp or email for availability, quotes and pre-orders.

HOW IT WORKS

In Stock or Pre-Order — always clear.

IN STOCK

Ready to ship

The model is in our Weston, FL inventory today and is listed on our online store. Order at orthodentalusa.com, by phone or WhatsApp.

PRE-ORDER

Reserve yours

The model is part of the range but is not currently in our inventory. Place a pre-order and our team will confirm availability and the estimated delivery time.

Reading a product row

- 1 **Reference** — Orthodentalusa part number to quote when ordering.
- 2 **Model & size** — instrument name, length and tip / joint type.
- 3 **Cutting capacity** — maximum wire the cutter is designed for.
- 4 **Link** — opens the matching product on orthodentalusa.com when it is listed there.

How to order

ONLINE
orthodentalusa.com
PHONE / WHATSAPP

+1 (954) 655-8487

EMAIL
sales@orthodentalusa.com

Product photographs are for identification; finishes and handle plating may vary slightly by production lot. Specifications are subject to change without notice. Pliers are for use by dental professionals only.

CONTENTS

What is inside

01	Distal End Cutters	19 models · 7 in stock	5
02	Flush Cutters	10 models · 7 in stock	9
03	Ligature & Pin Cutters	25 models · 12 in stock	12
04	Hard Wire & Heavy Duty Cutters	13 models · 9 in stock	18
05	Bracket & Band Removing Pliers	12 models · 8 in stock	21
06	Weingart Pliers	11 models · 5 in stock	24
07	Bird Beak & Light Wire Pliers	15 models · 7 in stock	27
08	Arch Forming & Nance Pliers	16 models · 7 in stock	31
09	Loop Forming, Three-Jaw & Tweed Pliers	20 models · 13 in stock	35
10	Jaraback, Hollow Chop & Utility Pliers	20 models · 11 in stock	40
11	Universal, Adams & Specialty Pliers	23 models · 11 in stock	45
12	Cinch-Back, Step & Torquing Pliers	27 models · 11 in stock	50
13	Lingual Series Pliers	7 models · 6 in stock	56
14	Aligner Pliers	19 models · 4 in stock	58
15	Mathieu Needle Holders & Mosquito Forceps	29 models · 22 in stock	62
16	Tweezers & Cotton Pliers	16 models · 14 in stock	68
17	Scalers & Orthodontic Instruments	30 models · 20 in stock	72
18	Gauges & Mapping Tools	33 models · 11 in stock	79
19	Spatulas, Knives & Scissors	17 models · 12 in stock	86
20	Photography & Intraoral Accessories	23 models · 3 in stock	90
21	Sterilization, Trays & Kits	15 models · 12 in stock	95
22	Additional In-Stock Instruments	26 models · 25 in stock	99
	Quick Reference Index (all models)	availability at a glance	105
	Materials, Warranty & Care	sterilization guide	109

01

SECTION Distal End Cutters

19

MODELS

7

IN STOCK

12

PRE-ORDER



N21150317

IN STOCK

Distal End Cutter, Universal Cut 12cm

Max: .021" (.54mm) x .025" (.64mm) • Min: .012" (.31mm)

Fitted with Immunity-HS Steel tips insert, with safety hold ability. Rivet joint. Stainless steel.

[View on orthodontalusa.com →](#)


N21150319

IN STOCK

Distal End Cutter, Universal Cut. Long Neck 12cm

Max: .021" (.54mm) x .025" (.64mm) • Min: .012" (.31mm)

Fitted with Immunity-HS Steel tips insert, with safety hold ability. Rivet joint. Stainless steel.

[View on orthodontalusa.com →](#)


N21150320

IN STOCK

Distal End Cutter, Universal Cut. Small Neck 11.5cm

Max: .021" (.54mm) x .025" (.64mm) • Min: .012" (.31mm)

Fitted with Immunity-HS steel tips insert, with safety hold ability. Rivet joint. Stainless steel.

[View on orthodontalusa.com →](#)


N21150321

IN STOCK

Distal End Cutter, Universal Cut, Small Neck, Maxi Look 12cm

Max: .021" (.54mm) x .025" (.64mm) • Min: .012" (.31mm)

Fitted with Immunity-HS Steel tips insert, with safety hold ability. Rivet joint. Stainless steel.

[View on orthodontalusa.com →](#)

★ CERTIFIED QUALITY

✓ ISO 9001:2008

✓ ISO 13485:2016

✓ CE CA006210



DEN-002

● PRE-ORDER

Distal End Cutter, Universal Cut Box Joint 11.5cm

Max: .021" (.54mm) x .025" (.64mm) • Min: .012" (.31mm)

Fitted with Immunity-HS Steel tips insert, with safety hold ability. Box joint. Stainless steel.

[Request pre-order →](#)


N21150326

● PRE-ORDER

Distal End Cutter, Universal Cut. Long Neck 12cm

Max: .021" (.54mm) x .025" (.64mm) • Min: .012" (.31mm)

Fitted with Immunity-HS Steel tips insert, with safety hold ability. Rivet joint. Stainless steel.

[Request pre-order →](#)


N21150327

● PRE-ORDER

Distal End Cutter, Universal Cut

Max: .021" (.54mm) x .025" (.64mm) • Min: .012" (.31mm)

Box Joint 12cm Fitted with immunity-HS Steel tips insert, with safety hold ability. Rivet joint. Stainless steel.

[Request pre-order →](#)


N21150328

● IN STOCK

Distal End Cutter, Universal Cut. Slim Pattern, Box Joint 12cm

Max: .021" (.54mm) x .025" (.64mm) • Min: .012" (.31mm)

Fitted with Immunity-HS Steel tips insert, with safety hold ability. Rivet joint. Stainless steel.

[View on orthodontalusa.com →](#)


N21150571

● PRE-ORDER

Distal End Cutter 12 cm

Max: .021" (.54mm) x .025" (.64mm) • Min: .012" (.31mm)

Fitted with Immunity-HS Steel tips insert, with safety hold ability. Rivet Joint 1/3 handle Gold Platted. Stainless steel.

[Request pre-order →](#)



N21150565

● PRE-ORDER

Distal End Cutter Universal 12.5cm

Max: .021" (.54mm) x .025" (.64mm) • Min: .012" (.31mm)

Fitted with Immunity-HS Steel tips insert, with safety hold ability. Box joint. Stainless steel.

[Request pre-order →](#)


N21150606

● PRE-ORDER

Distal End Cutter Universal 12.5cm

Max: .021" (.54mm) x .025" (.64mm) • Min: .012" (.31mm)

Fitted with Immunity-HS Steel tips insert, with safety hold ability. Rivet joint. Stainless steel.

[Request pre-order →](#)


N21150956

● PRE-ORDER

Distal End Cutter Universal 11.5cm

Max: .021" (.54mm) x .025" (.64mm) • Min: .012" (.31mm)

Fitted with Immunity-HS Steel tips insert, with safety hold ability. Rivet joint. Stainless steel.

[Request pre-order →](#)


N21150334

● PRE-ORDER

Distal End Cutter, Mini Head 13cm

Max: .021" (.54mm) x .025" (.64mm) • Min: .012" (.31mm)

Fitted with Immunity-HS Steel tips insert, Mini head allow access to narrow area with safety hold ability. Rivet Joint Stainless Steel.

[Request pre-order →](#)


N21150335

● PRE-ORDER

Distal End Cutter, Box Joint 12.5cm

Max: .021" (.54mm) x .025" (.64mm) • Min: .012" (.31mm) NICKEL TITANIUM

Fitted with Immunity-HS steel tips insert, with safety hold ability. Box joint 1/3 handle Gold plated. Stainless Steel.

[Request pre-order →](#)



N21150317L

IN STOCK

Distal End Cutter, Universal Cut. Long Handles 14cm

Max: .021" (.54mm) x .025" (.64mm) • Min: .012" (.31mm)

Fitted with Immunity-HS Steel tips insert, with safety hold ability. Rivet joint.
Stainless steel.

[View on orthodontalusa.com →](#)


N21151000

PRE-ORDER

Distal End Cutter

Max: .021" (.54mm) x .025" (.64mm) • Min: .012" (.31mm)

Fitted with Immunity-HS Steel tips insert, with safety hold ability. Rivet joint.
Stainless steel.

[Request pre-order →](#)


N21150606L

PRE-ORDER

Distal End Cutter, Universal 15cm

Max: .021" (.54mm) x .025" (.64mm) • Min: .012" (.31mm)

Fitted with Immunity-HS Steel tips insert, with safety hold ability. Rivet joint.
Stainless steel.

[Request pre-order →](#)


DEN-001

PRE-ORDER

Distal End Cutter, Universal Cut. Long Handles 14cm

Max: .021" (.54mm) x .025" (.64mm) • Min: .012" (.31mm)

Fitted with Immunity-HS Steel tips insert, with safety hold ability. Rivet joint.
Stainless steel.

[Request pre-order →](#)


N21150318

IN STOCK

Distal End Cutter, Universal Cut. Long Slim pattern 13cm

Max: .021" (.54mm) x .025" (.64mm) • Min: .012" (.31mm)

Fitted with Immunity-HS Steel tips insert, with safety hold ability. Rivet joint.
Stainless steel.

[View on orthodontalusa.com →](#)

02

SECTION

Flush Cutters

10

MODELS

7

IN STOCK

3

PRE-ORDER



N21150325

IN STOCK

Distal End Cutter, Universal Cut. Long Handles 14.5cm

Max: .021" (.54mm) x .025" (.64mm) • Min: .012" (.31mm)

Fitted with Immunity-HS Steel tips insert, with safety hold ability. Rivet joint. Stainless steel.

[View on orthodontalusa.com →](#)


N20S001

IN STOCK

Distal End Cutter, Universal Cut. Heavy pattern, Box joint, long Handles 15cm

Max: .021" (.54mm) x .028" (.71mm) • Min: .012" (.31mm)

Fitted with Immunity-HS Steel tips insert, with safety hold ability. Stainless steel. 1/3 handle Gold Plated.

[View on orthodontalusa.com →](#)


N21150322

IN STOCK

Flush Cutter without Hold 45 degree angled 12cm

Max: .016" (.41mm) Min: .012" (.31mm)

Fitted With Immunity-HS Steel Tips Insert, Without Safety hold ability. Rivet Joint. Stainless Steel.

[View on orthodontalusa.com →](#)


N21150323

IN STOCK

Flush Cutter with Hold 12cm

Max: .021" (.54mm) .028" (.71mm) • Min: .012" (.31mm)

This is a special designed flush cutter which holds the cut part of wires with the help of 1.40mm wire welded on the other side. Fitted With Immunity-HS Steel tips insert, with safety hold ability. Rivet joint. Stainless steel.

[View on orthodontalusa.com →](#)

★ CERTIFIED QUALITY

✓ ISO 9001:2008

✓ ISO 13485:2016

✓ CE CA006210



N21150323H

IN STOCK

Flush Cutter with Hold 13.5cm

Max: .021" (.54mm) x .025" (.64mm) • Min: .012" (.31mm)

This is a special designed flush cutter which holds the cut part of wires with the help of 0.45mm SS Sheet Spring, screwed on the inner side of the beak. Fitted With Immunity-HS Steel Tips Insert, With Safety hold ability.

Order: sales@orthodentalusa.com


N21150323F

IN STOCK

Flush Cutter with Hold, Heavy pattern, Long Handle 14cm

Max: .021" (.54mm) x .025" (.64mm) • Min: .014" (.36mm) NICKEL TITANIUM

This is a special designed flush cutter which holds the cut part of wires with the help of 1.20mm wire welded on the inner side. Fitted With Immunity-HS Steel tips insert, with safety hold ability. Rivet Joint. Stainless Steel.

Order: sales@orthodentalusa.com


N21150606F

PRE-ORDER

Flush Cutter with Hold 12.5 cm

Max: .021" (.54mm) x .025" (.64mm) • Min: .014" (.36mm) NICKEL TITANIUM

This is a special designed flush cutter which holds the cut part of wires with the help of 1.40mm wire welded on the other side. Fitted With Immunity-HS Steel tips insert, with safety hold ability. Rivet Joint. Stainless Steel.

Request pre-order →



N21150324

IN STOCK

Slim Flush Cutter with Hold 12.5cm

Max: .021" (.54mm) x .025" (.64mm) • Min: .014" (.36mm) NICKEL TITANIUM

This is a special designed flush cutter which holds the cut part of wires with the help of 1.40mm wire welded on the other side. Fitted With Immunity-HS Steel tips insert, with safety hold ability. Rivet joint. Stainless steel.

View on orthodentalusa.com →


N21150575

PRE-ORDER

Flush Cutter without Hold 12.5cm

Max: .021" (0.54mm) x .025" (.64mm) Min: .016" (0.41mm)

This is a special designed flush cut feature of this instrument allows the archwire to be cut flush with the end of the buccal tube, eliminating the possibility of patient discomfort resulting from archwire protrusion.Box Joint.

Request pre-order →



N21150345

● PRE-ORDER

Flush Cutter with Hold 13 cm

Cutting Capacity: Max: .021" (.54mm) x .025" (.64mm) • Min: .012" (.31mm)

This is a special designed flush cutter which holds the cut part of wires with the help of 1.40mm wire welded on the other side. Fitted With Immunity- HS Steel tips insert, with safety hold ability. Box Joint. Stainless Steel.

[Request pre-order →](#)

03

SECTION

Ligature & Pin Cutters

25

MODELS

12

IN STOCK

13

PRE-ORDER



N21150330

IN STOCK

Ligature & Pin Cutter straight 12cm

Max: .012" (.31mm) Soft Wires

Fitted With Immunity-HS Steel Tips Insert, Beak Length 12mm, Diamond honed edges, Rivet Joint. Stainless Steel.

[View on orthodontalusa.com →](#)


N21150330L

PRE-ORDER

Ligature & Pin Cutter straight Long Handles 14cm

Max: .012" (.31mm) Soft Wires

Fitted With Immunity-HS Steel Tips Insert, Beak Length 12mm, Diamond honed edges, Rivet Joint. Stainless Steel.

[Request pre-order →](#)


N21150342

IN STOCK

Ligature & Pin Cutter 15 degree angled 12cm

Max: .012" (.31mm) Soft Wires

Fitted With Immunity-HS Steel Tips Insert, Beak Length 12mm, Diamond honed edges, Rivet Joint. Stainless Steel.

[View on orthodontalusa.com →](#)


N21150862

PRE-ORDER

Ligature & Pin Cutter 15 degree Regular Beak T/C 12cm

Max: .012" (.31mm) to .016" (.41mm)

Fitted With Tungsten Carbide Tips Insert, Beak Length 8mm, Diamond honed edges, Rivet Joint. Stainless Steel. Soft Wires

[Request pre-order →](#)

★ CERTIFIED QUALITY

✓ ISO 9001:2008

✓ ISO 13485:2016

✓ CE CA006210



N21150572

● PRE-ORDER

Ligature & Pin cutter 12 cm

Max: .012" (.31mm) Soft Wires

Fitted With Immunity-HS Steel Tips Insert, Beak Length 8mm, Diamond honed edges, Rivet Joint. 1/3 Handle gold plated Stainless Steel.

[Request pre-order →](#)


N21150605

● PRE-ORDER

Ligature and Pin Cutter 12 cm

Max: .012" (.31mm) Soft Wires

Fitted With Immunity-HS Steel Tips Insert, Beak Length 12mm, Diamond honed edges, Rivet Joint. Stainless Steel.

[Request pre-order →](#)


N21150815

● PRE-ORDER

Ligature Cutter Long Handles 7 degree 14cm

Max: .012" (.31mm) Soft Wires

Fitted With Immunity-HS Steel Tips Insert, Beak Length 12mm, Diamond honed edges, Rivet Joint. Stainless Steel.

[Request pre-order →](#)


N21150816

● PRE-ORDER

Ligature Cutter 7 degree 12.5cm

Max: .012" (.31mm) Soft Wires

Fitted With Immunity-HS Steel Tips Insert, Beak Length 12mm, Diamond honed edges, Rivet Joint. Stainless Steel.

[Request pre-order →](#)


N21150817

● PRE-ORDER

Ligature Cutter 15 degree 12.5cm

Max: .012" (.31mm) Soft Wires

Fitted With Immunity-HS Steel Tips Insert, Beak Length 12mm, Diamond honed edges, Rivet Joint. Stainless Steel.

[Request pre-order →](#)



N21151002

● PRE-ORDER

Ligature & Pin Cutter wide beak 15 degree angled, 11cm

Max: .012" (.31mm) Soft Wires

Fitted With Immunity-HS Steel Tips Insert, Beak Length 8mm, Diamond honed edges, Rivet Joint. Stainless Steel.

[Request pre-order →](#)


N21150332

● IN STOCK

Mini Ligature & Pin Cutter 15 degree

Max: .012" (.31mm) Soft Wires

angled, Slim Beak 12cm Fitted With Immunity-HS Steel Tips Insert, Beak Length 9mm, Diamond honed edges, Rivet Joint. Stainless Steel.

[View on orthodontalusa.com →](#)


N21150332L

● PRE-ORDER

Mini Ligature & Pin Cutter 15 degree

Max: .012" (.31mm) Soft Wires

angled, Slim Beak Long Handles 14cm Fitted With Immunity-HS Steel Tips Insert, Beak Length 9mm, Diamond honed edges, Rivet Joint. Stainless Steel.

[Request pre-order →](#)


N21150333

● IN STOCK

Ligature & Pin Cutter wide beak 15 degree angled, 11cm

Max: .012" (.31mm) Soft Wires

Fitted With Immunity-HS Steel Tips Insert, Beak Length 8mm, Diamond honed edges, Rivet Joint. Stainless Steel.

[View on orthodontalusa.com →](#)


N21150352

● IN STOCK

Micro Ligature & Pin Cutter T/C 15 degree angled 11cm

Max: .012" (.31mm) to .016" (.41mm)

Fitted With Immunity-HS Steel Tips Insert, Beak Length 9mm, Diamond honed edges. Rivet Joint. Stainless Steel. Wires

[View on orthodontalusa.com →](#)



N21150351

IN STOCK

Micro Ligature & Pin Cutter 15 degree light model 11.5 cm

Max: .012" (.31mm) Soft Wires

Fitted With Immunity-HS Steel Tips Insert, Beak Length 9.5mm, Diamond honed edges, Rivet Joint. Stainless Steel.

Order: sales@orthodentalusa.com


N21150818

PRE-ORDER

Micro Ligature Cutter 7 degree 12cm

Max: .012" (.31mm) Soft Wires

Fitted With Immunity-HS Steel Tips Insert, Beak Length 12mm, Diamond honed edges, Rivet Joint. Stainless Steel.

Request pre-order →



N21150354

IN STOCK

Micro Ligature & Pin Cutter 15 degree angle 11 cm

Max: .012" (.31mm) Soft Wires

Fitted With Immunity-HS Steel Tips Insert, Beak Length 8mm, Diamond honed edges, Rivet Joint. Stainless Steel.

Order: sales@orthodentalusa.com


N21150603

PRE-ORDER

Mini Ligature and pin Cutter 12cm

Max: .012" (.31mm) Soft Wires

Fitted With Immunity-HS Steel Tips Insert, Beak Length 8mm, Diamond honed edges, Rivet Joint. 1/3 Handle Gold Platted. Stainless Steel.

Request pre-order →



N21150329

IN STOCK

Ligature Cutter 45 degree angled wide handles very fine tips 13 cm

Max: .012" (.31mm) Min: .008" (.20mm)

Fitted With Immunity-HS Steel Tips Insert, Tips length 4.5mm. Rivet Joint. Stainless Steel Soft Wires

View on orthodentalusa.com →



N21150336

● PRE-ORDER

Ligature & Pin cutter, Box Joint 12.5 cm

Max: .012" (.31mm) Soft Wires

Fitted With Immunity-HS Steel Tips Insert, Beak Length 9.5mm, Diamond honed edges, Box Joint. 1/3 Handle Gold Platted. Stainless Steel.

[Request pre-order →](#)


N21150S002

● IN STOCK

Micro Ligature & Pin Cutter 7 degree angle 12.5cm

Max: .012" (.31mm) Soft Wires

Fitted With Immunity-HS Steel Tips Insert, Beak Length 7.5mm, Diamond honed edges, Box Joint. 1/3 Handle Gold Platted. Stainless Steel.

[View on orthodontalusa.com →](#)


N21150S003

● IN STOCK

Ligature & Pin Cutter straight 13cm

Max: .012" (.31mm) Soft Wires

Fitted With Immunity-HS Steel Tips Insert, Beak Length 9.5mm, Diamond honed edges, Box Joint. 1/3 Handle Gold Platted. Stainless Steel.

[View on orthodontalusa.com →](#)


N21150331

● IN STOCK

Mini Ligature & Pin Cutter Str, Slim

Max: .012" (.31mm) Soft Wires

Beak 12cm Fitted With Immunity-HS Steel Tips Insert, Beak Length 8mm, Diamond honed edges, Box Joint. 1/3 Handle Gold Platted. Stainless Steel.

[View on orthodontalusa.com →](#)


N21150355

● IN STOCK

Ligature & Pin Cutter str, Light Slim Model 7 degree 12.5cm

Max: .012" (.31mm) · Min: .008" (.20mm)

Fitted With Immunity-HS Steel Tips Insert, Beak Length 10mm, Diamond honed edges, Box Joint. Stainless Steel. Soft Wires

[Order: sales@orthodontalusa.com](mailto:sales@orthodontalusa.com)



N21150566

● PRE-ORDER

Ligature and pin Cutter 12cm

Max: .012" (.31mm) Soft Wires

Fitted With Immunity-HS Steel Tips Insert, Beak Length 8mm, Diamond honed edges, Box Joint. 1/3 Handle Gold Platted. Stainless Steel.

[Request pre-order →](#)

04

SECTION

Hard Wire & Heavy Duty Cutters

13

MODELS

9

IN STOCK

4

PRE-ORDER



N21150353

IN STOCK

Hard Wire Cutter Straight 11.5cm

Max: .021" (.54mm) to .028" (.71mm)

Fitted With Tungsten Carbide Tips Insert, Beak Length 8mm, Diamond honed edges, Rivet Joint. Stainless Steel. Hard Wires

[View on orthodontalusa.com →](#)


N21150353A

IN STOCK

Hard Wire Cutter 15 degree 11.5cm

Max: .021" (.54mm) to .028" (.71mm)

Fitted With Tungsten Carbide Tips Insert, Beak Length 8mm, Diamond honed edges, Rivet Joint. Stainless Steel. Hard Wires

[View on orthodontalusa.com →](#)


N21150353L

IN STOCK

Hard Wire Cutter Str Long Beak 12cm

Max: .021" (.54mm) to .028" (.71mm)

Fitted With Tungsten Carbide Tips Insert, Beak Length 9mm, Diamond honed edges, Rivet Joint. Stainless Steel. Hard Wires

[View on orthodontalusa.com →](#)


N21150353H

IN STOCK

Hard Wire Cutter 15 degree Long

Max: .021" (.54mm) to .028" (.71mm),

Beak 13.5cm. Fitted With Tungsten Carbide Tips Insert, Beak Length 11mm, Diamond honed edges, Rivet Joint. Stainless Steel. Hard Wires

[View on orthodontalusa.com →](#)

★ CERTIFIED QUALITY

✓ ISO 9001:2008

✓ ISO 13485:2016

✓ CE CA006210



N21150353B

IN STOCK

Hard Wire Cutter Black, 15 degree

Max: .021" (.54mm) to .025" (.64mm)

12.5cm. Fitted With Tungsten Carbide Tips Insert, Beak Length 8.5mm, Diamond honed edges, Hex Screw Joint. Back Chromed Stainless Steel. Hard Wires

[View on orthodontalusa.com →](#)


N21150560

PRE-ORDER

Hard Wire Cutters 13cm

Max: .021" (.54mm) to .028" (.71mm)

Fitted With Tungsten Carbide Tips Insert, Beak Length 8mm, Diamond honed edges, Box Joint. Stainless Steel. Hard Wires

[Request pre-order →](#)


N21150819

PRE-ORDER

Hard Wire Cutter str 12.5cm

Max: .021" (.54mm) to .028" (.71mm)

Fitted With Tungsten Carbide Tips Insert, Beak Length 8mm, Diamond honed edges, Rivet Joint. Stainless Steel. Hard Wires

[Request pre-order →](#)


N21150820

PRE-ORDER

Hard Wire Cutter 15 degree 12.5cm

Max: .021" (.54mm) to .028" (.71mm)

Fitted With Tungsten Carbide Tips Insert, Beak Length 8mm, Diamond honed edges, Rivet Joint. Stainless Steel. Hard Wires

[Request pre-order →](#)


N21150356

IN STOCK

Heavy Duty Cutter, Multi-action

Max: .035" (.89mm) Hard Wires

Stainless Steel. 14cm Ideal to cut heavy wires easily with the help of spring action cutting mechanism without putting extra force. Stainless Steel.

[View on orthodontalusa.com →](#)



N21150357

IN STOCK

Heavy Duty Cutter, Multi-action. Carbon Steel. 14.5cm

Max: .035" (.89mm) Hard Wires

Ideal to cut heavy wires easily with the help of spring action cutting mechanism without putting extra force. Stainless Steel beak black coated.

[View on orthodontalusa.com →](#)


N21150358

PRE-ORDER

Heavy Duty Cutter, Multi-action Mechanism 9" T/C 23cm

Max: .040" (1.00mm) Hard Wires

Fitted with T/C Tips, Ideal to cut heavy wires easily with the help of 4 hinge mechanism without putting extra force. Stainless Steel. 1/3 handles gold plated.

[Request pre-order →](#)


N21150359

IN STOCK

Heavy Duty Cutter, Multi-action Mechanism 7" T/C 18cm

Max: .035" (.89mm) Hard Wires

Fitted with T/C Tips, Ideal to cut heavy wires easily with the help of 4 hinge mechanism without putting extra force. Stainless Steel. 1/3 handles gold plated.

[View on orthodontalusa.com →](#)


N21150360

IN STOCK

Heavy Wire Cutter. 15.5cm

Max: .030" (.76mm) Hard Wires

Good for cutter heavy wires upto .030", without spring action. Slim shaped beak. Hardened Stainless Steel.

[View on orthodontalusa.com →](#)

05

SECTION

Bracket & Band Removing Pliers

12

MODELS

8

IN STOCK

4

PRE-ORDER



N21150380

IN STOCK

Bracket Removing Plier 12.5cm angled

Sharp beak wedge can easily go under bracket base. Lifts off bracket from tooth surface with no discomfort to patient due to Firm Round Tips Grip assures, easy in horizontal and vertical use. Rivet Joint. Stainless Steel.

[View on orthodontalusa.com →](#)


N21150381

IN STOCK

Bracket Removing Plier 13cm Straight

Sharp beak wedge can easily go under bracket base. Lifts off bracket from tooth surface with no discomfort to patient due to Firm Round Tips Grip assures, easy in horizontal and vertical use. Rivet Joint. Stainless Steel.

[View on orthodontalusa.com →](#)


N21150383

IN STOCK

Bracket Removing Plier 12.5cm angled T/C Tips

Sharp beak wedge can easily go under bracket base. Lifts off bracket from tooth surface with no discomfort to patient due to Firm Sloped Tips Grip assures, easy in horizontal and vertical use. Box Joint. Stainless Steel.

[View on orthodontalusa.com →](#)


N21150577

PRE-ORDER

Bracket Removing Pliers Angled Box Joint 13cm

Sharp beak wedge can easily go under bracket base. Lifts off bracket from tooth surface with no discomfort to patient due to Firm Sloped Tips Grip assures, easy in horizontal and vertical use. Box Joint. Stainless Steel.

[Request pre-order →](#)

★ CERTIFIED QUALITY

✓ ISO 9001:2008

✓ ISO 13485:2016

✓ CE CA006210



N21150371

IN STOCK

Band Removing Pliers long neck 13cm

Ideal for posterior band removing. Replaceable high heat Pads, work upward and downward to remove band. Pads are high Heat resistant autoclavable. Rivet Joint. Stainless Steel. PADS are available in 4mm, 5mm, 6mm, dia.

[View on orthodontalusa.com →](#)


N21150373

IN STOCK

Band removing pliers short neck 13cm

Ideal for posterior band removing. Replaceable high heat Pads, work upward and downward to remove band. Pads are high Heat resistant autoclavable. Rivet Joint. Stainless Steel. PADS are available in 4mm, 5mm, 6mm dia.

[View on orthodontalusa.com →](#)


N21150374

IN STOCK

Band Removing Pliers 13.5cm

Ideal for posterior band removing. Replaceable high heat Pads, work upward and downward to remove band. Pads are high Heat resistant autoclavable. Rivet Joint. Stainless Steel. PADS are available in 4mm, 5mm, 6mm dia.

[View on orthodontalusa.com →](#)


N21150812

PRE-ORDER

Band Removing Pliers Long Neck 13.5cm

Ideal for posterior band removing. Replaceable high heat PADS, work upward and downward to remove band. Pads are high Heat resistant autoclavable. Rivet Joint. Stainless Steel. PADS are available in 4mm, 5mm, 6mm, dia.

[Request pre-order →](#)


N21150824

PRE-ORDER

Band Removing Pliers Long beak with T/C Tip 14.5cm

Ideal for posterior band removing. Replaceable high heat PADS, work upward and downward to remove band. Pads are high Heat resistant autoclavable. Box Joint. Stainless Steel. PADS are available in 4mm, 5mm, 6mm, dia.

[Request pre-order →](#)



N21150578

● PRE-ORDER

Band Removing Pliers Box Joint 14cm

Ideal for posterior band removing. Replaceable high heat Pads, work upward and downward to remove band. Pads are high Heat resistant autoclavable. Box Joint. Stainless Steel. PADS are available in 4mm, 5mm, 6mm, dia.

[Request pre-order →](#)


TEF-201

● IN STOCK

Teflon Studs / Pads

These are band remover spare pads non sterilizable These pads are available in 4mm, 5mm, 6mm Dia tip size These pads comes in set of 5 pads.

[View on orthodontalusa.com →](#)


PHI-105

● IN STOCK

High Heat Band Removing Pads

These are high heat sterilization band remover spare pads These can resist heat up to 400° F These pads are available in 4mm, 5mm, 6mm Dia tip size These pads comes in set of 5 pads.

[View on orthodontalusa.com →](#)

06

SECTION Weingart Pliers

11

MODELS

5

IN STOCK

6

PRE-ORDER



N21150378

IN STOCK

Weingart Pliers 13cm

Serrated jaws for firm grip, angled shape gives smooth working with easy access to the specified job. Rivet Joint. Stainless Steel.

[View on orthodontalusa.com →](#)


N21150377

IN STOCK

Weingart Pliers 13cm

Serrated jaws for firm grip, angled shape gives smooth working with easy access to the specified job. Rivet Joint. Stainless Steel.

[View on orthodontalusa.com →](#)


N21150379

IN STOCK

Weingart Pliers 13.5cm

Serrated jaws for firm grip, angled shape gives smooth working with easy access to the specified job. Box Joint. Stainless Steel.

[View on orthodontalusa.com →](#)


N21150395

PRE-ORDER

Weingart Pliers Box Joint 14 cm

Serrated jaws for firm grip, angled shape gives smooth working with easy access to the specified job. Box Joint. Stainless Steel.

[Request pre-order →](#)

★ CERTIFIED QUALITY

✓ ISO 9001:2008

✓ ISO 13485:2016

✓ CE CA006210



N21150396

● PRE-ORDER

Slim Weingart Pliers Box Joint 14.5 cm

Slim Serrated jaws for firm grip, angled Slim shape gives smooth working with easy access to the specified job. Box Joint. Stainless Steel.

[Request pre-order →](#)


N21150375

● IN STOCK

Slim Weingart Pliers 13.5cm

Serrated jaws for firm grip, angled shape gives smooth working with easy access to the specified job. Rivet Joint. Stainless Steel.

[View on orthodontalusa.com →](#)


N21150376

● IN STOCK

Slim Weingart Pliers 13.5cm

Serrated jaws for firm grip, angled slim and micro shape gives smooth working with easy access to the specified job. Rivet Joint. Stainless Steel.

[View on orthodontalusa.com →](#)


N21150638

● PRE-ORDER

Slim Weingart Pliers Narrow

Mini Beak 60 degree 12.5cm.. Serrated jaws for firm grip, angled slim and micro 60° shape gives more precise working with easy access to the specified job. Rivet Joint. Stainless Steel.

[Request pre-order →](#)


N21150863

● PRE-ORDER

Weingart Pliers (MINI) 13cm

Serrated jaws for firm grip, angled shape gives smooth working with easy access to the specified job. Rivet Joint. Stainless Steel.

[Request pre-order →](#)



N21150397

● PRE-ORDER

Slim Weingart 90 degree 13.5cm

Serrated jaws for firm grip, angled slim shape gives smooth working with easy access to the specified job. Specially Lingual Area. Rivet Joint. Stainless Steel.

[Request pre-order →](#)


N21150957

● PRE-ORDER

Slim Weingart Pliers 13cm

Serrated jaws for firm grip, angled shape gives smooth working with easy access to the specified job. Rivet Joint. Stainless Steel.

[Request pre-order →](#)

07

SECTION Bird Beak & Light Wire Pliers

15

MODELS

7

IN STOCK

8

PRE-ORDER



N21150390

IN STOCK

Bird Beak Pliers 12cm

Max: 0.030" (.76 mm) Wires

Utility Pliers popular for working on round & rectangular wires. Tips are hardened. Rivet Joint. Stainless Steel.

[View on orthodontalusa.com →](#)


N21150958

PRE-ORDER

Bird Beak Pliers Fine 11.5cm

Max: 0.030" (.76 mm) Wires

Utility Pliers popular for working on round & rectangular wires. Tips are hardened. Fine beaks gives more precision for working on round & rectangular wires Rivet Joint. Stainless Steel.

[Request pre-order →](#)


N21150576

PRE-ORDER

Bird Beak Pliers Box Joint 12.5cm

Max: .030" (.76 mm) Wires

Utility Pliers popular for working on round & rectangular wires. Tips are hardened. Rivet Joint. Stainless Steel.

[Request pre-order →](#)


N21150813

PRE-ORDER

Slim Sloped Bird Beak 12.5cm

Max: .019" (.48 mm) x .025" (.64 mm)

Utility pliers popular for working on round & rectangular wire. Tips are hardened for Making Smooth Loops. Box Joint. Stainless Steel. Wires

[Request pre-order →](#)

★ CERTIFIED QUALITY

✓ ISO 9001:2008

✓ ISO 13485:2016

✓ CE CA006210



N21150637

● PRE-ORDER

Slim Sloped Bird Beak 12.5cm

Max: .019" (.48 mm) x .025" (.64 mm)

Utility pliers popular for working on round & rectangular wire. Tips are hardened for Making Smooth Loops. Rivet Joint. Stainless Steel. Wires

[Request pre-order →](#)


N21150388

● IN STOCK

Slim Bird Beak Pliers Tapered Beak 12.5cm

Max: .019" (.48 mm) x .025" (.64 mm)

Utility pliers popular for working on round & rectangular wire. Tips are hardened for Making Smooth Loops. Box Joint. Stainless Steel. Wires

[View on orthodontalusa.com →](#)


N21150390H

● IN STOCK

Bird Beak Pliers Heavy and serrated jaw 12.5cm

Max: .036" (0.91mm) Hard wires

Heavy duty bird beak is ideal for functional appliance fabrication. Square tip is serrated to grip and secure lab size wires while bending. Rivet Joint. Stainless Steel

[View on orthodontalusa.com →](#)


N21150822

● PRE-ORDER

Lab Bird Beak Pliers with Serrated Tips 14cm

Max: .036" (0.91mm) Hard wires

Heavy duty bird beak is ideal for functional appliance fabrication. Square tip is serrated to grip and secure lab size wires while bending. Box Joint. Stainless Steel

[Request pre-order →](#)


N21150823

● PRE-ORDER

Bird Beak Long Beak for Heavy Wires 12.5cm

Max: .036" (0.91mm) Hard wires

Fitted With Immunity-HS Steel Tips Insert, this Strong Wire Banding Pliers is ideal for functional appliance fabrication for heavy Wires, Rivet Joint, Stainless Steel.

[Request pre-order →](#)



N21150631

PRE-ORDER

Light Wire Bending Pliers with grooves tips 12.5cm

Fitted with immunity-HS Steel tips inserts, Utility Pliers popular for working on round wires. with the addition of 2 grooves at the tips. the first groove is 1mm & the second groove is 2mm from the tip on the square side.

[Request pre-order →](#)


N21150392SB

PRE-ORDER

Light Wire Bending Plier Small Beak 12.4cm

Max: .021" (.54 mm) x .025" (.64 mm)

Utility Pliers popular for working on round & Rectangular wires. Tips are hardened. Rivet Joint. Stainless Steel. Wires

[Request pre-order →](#)


N21150392

IN STOCK

Light Wire Bending Pliers Long 12.5cm

Max: .021" (.54 mm) x .025" (.64 mm)

Utility Pliers popular for working on round & Rectangular wires. Tips are hardened. Rivet Joint. Stainless Steel. Wires

[View on orthodontalusa.com →](#)


N21150390W/C

IN STOCK

Bird Beak Pliers with cutters 12cm

Max: .021" (.54 mm) x .025" (.64 mm)

Utility Pliers popular for working with round & rectangular wires. Cutter formation makes this Plier popular among doctors for working and cutting at the same time. Tips are hardened. Rivet Joint. Stainless Steel. Wires

[View on orthodontalusa.com →](#)


N21150389

IN STOCK

Bird Beak pliers with cutters 12.5cm Tapered Beak

Utility pliers popular for working with round & rectangular wires Cutter formation makes this plier popular among doctors for working and cutting at the same time. Tips are hardened. Box Joint. Stainless Steel.

[View on orthodontalusa.com →](#)



N21150391

IN STOCK

Light wire bending Pliers with cutters 13 cm

Max: .021" (.54 mm) x .025" (.64 mm) Wires

Utility Pliers popular for working with round & rectangular wires. Cutter formation makes this Plier popular among doctors for working and cutting at the same time. Long and Slim beak gives precise. Tips are hardened.

[View on orthodentalusa.com](https://orthodentalusa.com) →

08

SECTION

Arch Forming & Nance Pliers

16

MODELS

7

IN STOCK

9

PRE-ORDER



N21150413

IN STOCK

TWEED Arch Forming Pliers 12cm

Max: .022" (.56mm) x .028 (.71mm)

The Pliers is designed for double back or triple back bending. Rivet Joint. Stainless Steel. Wires /

[View on orthodontalusa.com →](#)


N21150417

IN STOCK

TWEED Arch Forming Pliers 12.5cm

Max: .022" (.56mm) x .028 (.71mm)

The Pliers is designed for double back or triple back bending. Wide Handles for extra grip. Rivet Joint. Stainless Steel. Wires /

[View on orthodontalusa.com →](#)


N21150303

PRE-ORDER

Tweed Arch Bending Pliers

Max: .022" (.56mm) x .028 (.71mm)

Box Joint 12.5 cm. The Pliers is designed for double back or triple back bending. Box Joint. Stainless Steel. Wires

[Request pre-order →](#)


N21150412

IN STOCK

Slim Arch Forming Pliers 12.5cm

Max: .022" (.56mm) x .028 (.71mm)

The pliers is designed for double back or triple back bending. Rivet Joint. Stainless Steel. Wires

[Order: sales@orthodontalusa.com](mailto:sales@orthodontalusa.com)

★ CERTIFIED QUALITY

✓ ISO 9001:2008

✓ ISO 13485:2016

✓ CE CA006210



N21150411

IN STOCK

Rectangular Arch Forming Pliers 12.5cm

Max: .022" (.56mm) x .028 (.71mm)

The pliers is designed for double back or triple back bending. River Joint. Stainless Steel. Wires

[View on orthodontalusa.com →](#)


N21150568

PRE-ORDER

Arch forming pliers box joint 13cm

Max: .022" (.56mm) x .028 (.71mm)

The pliers is designed for double back or triple back bending. Box Joint. Stainless Steel. Wires

[Request pre-order →](#)


N21150410

IN STOCK

Lingual Arch Forming Pliers 12cm

Max: .036" (.91mm) Wires

The pliers is designed for double back or triple back bending can accommodate .030, .036 wires. Rivet Joint. Stainless Steel. 4 Hole. Size:- .80mm, .90mm, 1.25mm, .90mm

[View on orthodontalusa.com →](#)


N21150405

PRE-ORDER

Lingual Arch Forming Pliers 12cm

Max: .036" (.91mm) Wires

The pliers is designed for double back or triple back bending can accommodate .030, .036 wires. Rivet Joint. Stainless Steel. 4 Hole. Size:- .90mm, 1.20mm, 1.50mm, 1.40mm

[Request pre-order →](#)


N21150407

PRE-ORDER

Arch Forming Pliers 12cm

Max: .036" (.91mm) Wires

For individual construction and activation of lingual / palated bone. using these pliers, the bent tips of the wire will fit into all lingual / palated tubes. Rivet Joint. Stainless Steel. 4 Hole. Size:- .90mm, .90mm, 1.00mm,

[Request pre-order →](#)



N21150382

IN STOCK

Nance Pliers 12cm

For Precise forming of loops in four lengths 3, 4, 5, & 6 mm. Rounded edges prevent scoring wire. Rivet Joint. Stainless Steel.

[View on orthodontalusa.com →](#)


N21150864

PRE-ORDER

Nance Pliers 12.5cm

Max: .028" (.71 mm) Wires

For precise forming of loops in four lengths 3, 4, 5 & 6mm. Rounded edges prevent scoring wire. Rivet Joint. Stainless Steel.

[Request pre-order →](#)


N21150580

PRE-ORDER

Nance Pliers Box Joint 12.5 cm

Max: .028" (.71 mm) Wires

For precise forming of loops in four lengths 3, 4, 5 & 6mm. Rounded edges prevent scoring wire. Box Joint. Stainless Steel.

[Request pre-order →](#)


N21150639

PRE-ORDER

Nance Pliers Double Sided Loop Forming 12cm

Max: .028" (.71 mm) Wires

For Precise forming of loops in four lengths 2.50mm, 4.00mm, 6.50mm, & 7.50mm. Rounded edges prevent scoring wire. Double sided steps formation can be used performing loop either sides. Rivet Joint. Stainless Steel.

[Request pre-order →](#)


N21150761

PRE-ORDER

Nance Double Sided Pliers 12cm

Max: .028" (.71 mm) Wires

For Precise forming of loops Rounded edges prevent scoring wire. Rivet Joint. Stainless Steel. Set of 2 Pieces

[Request pre-order →](#)



N21150304

PRE-ORDER

Optical Loop Forming Pliers Wide Handles 12.5cm

Fine finished round & concave beaks, Excellent for bending of round and edge-wise, affording firm grip, without nicking wire. Can be used for loop tip back and gable bends.

[Request pre-order →](#)


N21150911

IN STOCK

Optical Pliers 12.5cm

Finely finished round and concave beaks, both of T/C Tips inserts. Excellent for bending of round and edge-wise, affording firm grip, without nicking wire. Can be used for loop tip back and gable bends. Box Joint. Stainless Steel

[View on orthodontalusa.com →](#)

09

SECTION

Loop Forming, Three-Jaw & Tweed Pliers

20

MODELS

13

IN STOCK

7

PRE-ORDER



N21150915

● IN STOCK

Loop Forming Pliers, Small Beak 12cm

Max: .028" (0.71mm) Hard wires

Designed for bending small springs and loops. Rivet Joint. Stainless Steel

[View on orthodontalusa.com →](#)



N21150916

● IN STOCK

Loop Forming Pliers, Long Beak 12.5cm

Max: .028" (0.71mm) Hard wires

Designed for bending small springs and loops. Long Beak provides access area for working. Rivet Joint. Stainless Steel

[View on orthodontalusa.com →](#)



N21150955

● IN STOCK

Ruhland Loop Forming Pliers 12.5cm

Max: .043" (1.10mm) Hard wires

Designed for bending coffin springs, labial arches, loops etc. Rivet Joint. Stainless Steel

Order: sales@orthodontalusa.com



N21150385

● IN STOCK

Three Jaw Pliers 12.5cm

Max: .028" (.71 mm) Wires

Perfect aligned tips, tapered for work on delicate bending. Rivet Joint. Stainless Steel.

[View on orthodontalusa.com →](#)

★ CERTIFIED QUALITY

✓ ISO 9001:2008

✓ ISO 13485:2016

✓ CE CA006210



N21150579

● PRE-ORDER

Three Prong Pliers Box Joint 13 cm

Max: .030" (.76 mm) Wires

Perfect aligned tips, Rounded for work on delicate bending. Box Joint. Stainless Steel.

[Request pre-order →](#)



N21150386

● IN STOCK

Aderer Pliers 13cm

Max: .030" (.76 mm) Wires

Perfect aligned tips, tapered for work on delicate bending. Box Joint. Stainless Steel.

[View on orthodontalusa.com →](#)



N21150305

● PRE-ORDER

Three Jaw Pliers Rounded 12cm

Max: .030" (.76 mm) Wires

Perfect aligned tips, Rounded for work on delicate bending. Box Joint. Stainless Steel.

[Request pre-order →](#)



N21150370

● PRE-ORDER

NiTi Wire Bender (3 Jaw) Slim

Max: .028" (.71 mm) Wires

Box Joint 13cm. Perfect aligned tips, Rounded for work on delicate bending. Box Joint. Stainless Steel.

[Request pre-order →](#)



N21150387

● IN STOCK

Mini Aderer Pliers 11.5cm

Max: .020" (.51 mm) Wires

Perfect aligned tips, tapered for work on delicate: bending. Box Joint. Stainless Steel.

[View on orthodontalusa.com →](#)



N21153023

● PRE-ORDER

Lab Heavy Three Prong / Heavy Aderer 14.5cm

Max: .36" (.92 mm) Wires

Perfect aligned tips, tapered for work on Heavy Lab Wire bending. Box Joint. Stainless Steel.

[Request pre-order →](#)



N21150422

● IN STOCK

Tweed Pliers 12.5cm

Max: .022" (.56mm) x .028" (.71mm)

The pliers gives perfect rings to the wire forming 1.25,, 2.00mm and 2.50mm radius. Rivet Joint. Stainless Steel. Wires

[View on orthodontalusa.com →](#)



N21150420

● IN STOCK

Tweed Pliers 12.5cm

Max: .022" (.56mm) x .028" (.71mm)

The pliers gives perfect rings to the wire forming 1.25mm, 1.75mm and 2.25mm radius. Rivet Joint. Stainless Steel. Wires

[View on orthodontalusa.com →](#)



N21150423

● IN STOCK

Tweed forming Pliers 13.5 cm

Max: .022" (.56mm) x .028" (.71mm)

Have a fixed, one-piece tip with three forming sections and a convex/concave beak for precisely forming Omega loops. The tips and body are formed from one piece of steel for strength and longevity.

[View on orthodontalusa.com →](#)



N21150424

● IN STOCK

Omega Loop Forming Pliers, Tweed Type 12.5cm

Max: .022" (.56mm) x .025" (.64mm)

The pliers gives perfect rings to the wire forming 1.25mm, 1.75mm and 2.00mm radius Replaceable loop forming tip ability. Rivet Joint. Stainless Steel. Extra loop forming tip 115-0424T Extra Hex Screw tip 115-0424S Wires

[View on orthodontalusa.com →](#)



N21150425

● IN STOCK

Young Loop Bending Pliers 12.5cm

Max: .028" (0.71mm) Hard wires

Universal pliers for different size loops. Ideal for labial arches. Rivet Joint. Stainless Steel

[View on orthodontalusa.com →](#)



N21150339

● PRE-ORDER

Young Pliers Box Joint 13.5 cm

Max: .028" (0.71mm) Hard wires

Universal pliers for different size loops. Ideal for labial arches. Box Joint. Stainless Steel

[Request pre-order →](#)



N21150435

● PRE-ORDER

Auxiliary Spring Forming Pliers 13cm

Ideal pliers for forming auxiliary spring components on .012", .014" and .016" round wires. Wide handles give firm gripping while working. Box Joint. Stainless Steel.

[Request pre-order →](#)



N21150421

● IN STOCK

Loop Forming Pliers with cutters. Tweed Type 12.5cm

Max: .016" (.41mm) x .022" (.56mm) • Max: .020" (.51mm) Round Wire

Also known as Kim Pliers. The plier gives perfect rings to the wire forming .060", makes this plier popular and gives additional working with time saving. Rivet Joint. Stainless Steel Rectangular Wire

[View on orthodontalusa.com →](#)



N21150426

● PRE-ORDER

Kim Pliers

Max: .016" (.41mm) x .022" (.56mm) • Max: .020" (.51mm) Round Wire

The plier gives perfect rings to the wire forming .060", makes this plier popular and gives additional working with time saving. Rivet Joint. Stainless Steel Rectangular Wire

[Request pre-order →](#)



N21150393

IN STOCK

Jaraback Pliers 13cm

Max: .025" (.64mm) Wires

Ideal for precise bending and forming loops in light wires up to .020". Set of 3 grooves assures a firm grip. Rivet Joint. Stainless Steel.

[View on orthodentalusa.com](https://orthodentalusa.com) →

10

SECTION

Jaraback, Hollow Chop & Utility Pliers

20

MODELS

11

IN STOCK

9

PRE-ORDER



N21150338

● PRE-ORDER

Jaraback Pliers Box joint (Slim) 14 cm

Max: .025" (.64mm)

Ideal for precise bending and forming loops in light wire. Set of 3 grooves assures a firm grip. Box Joint. Stainless Steel.

[Request pre-order →](#)


N21150393W/C

● IN STOCK

Jaraback Pliers with cutters 13.5cm

Max: .016" (.41mm) .022" (.56mm)

Ideal for precise bending and forming loops in light wire. Set of 3 grooves makes it more versatile. Box Joint. Stainless Steel.

[View on orthodontalusa.com →](#)


N21150394

● IN STOCK

Jaraback Pliers with T/C cutters 13.5cm

Max: .022" (.56mm) x .025" (.64mm) • Min: .012" (.31mm) x .016" (.41mm)

Ideal for precise bending and forming loops in light wires Set of 3 grooves assures a firm grip. T/C Tips Cutter formation makes it more versatile and cut hard wire easily. Box Joint. Stainless Steel.

[View on orthodontalusa.com →](#)


N21150408

● IN STOCK

Arch Forming Pliers 12cm

The smooth concave and convex surface allows forming of arch wires. Spring action. Rivet Joint. Stainless Steel.

[View on orthodontalusa.com →](#)

★ CERTIFIED QUALITY

✓ ISO 9001:2008

✓ ISO 13485:2016

✓ CE CA006210



N21150414

IN STOCK

Hollow Chop Arch Forming Pliers 12cm

The smooth concave and convex plain surface allows forming of arch wires easily. Rivet Joint. Stainless Steel.

[View on orthodontalusa.com →](#)


N21150409

IN STOCK

Hollow Chop Arch Forming pliers with 3 Grooves 12cm

The smooth concave and convex surface with grooves allows forming of arch wires easily. Sizes of grooves .016" .018" .022". (.41, .46, .56mm) Rivet Joint. Stainless Steel.

[View on orthodontalusa.com →](#)


N21150419

PRE-ORDER

Arch Forming Pliers 13.5cm

Ideal Plier to create a smooth arch on both round wires up to .040" and rectangular wires up to .021" x .028". Beveled edges eliminate the nicking of the wire during operation. Rivet Joint. Stainless Steel

[Request pre-order →](#)


N21150436

PRE-ORDER

Wire Bending Pliers wide curve 12.5cm

Ideal for making Gradual "accentuated curve" on .019x.025" SS wire. Wider curved and smooth edges gives fine and scoring free wire bending. Box Joint. Stainless Steel.

[Request pre-order →](#)


N21150570

PRE-ORDER

Hollow Chop Pliers with Grooves

Box Joint 12.5cm The smooth concave and convex surface with grooves allows forming of arch wires easily. Sizes of grooves .016" .018" .022". (.41, .46, .56mm) Rivet Joint. Stainless Steel.

[Request pre-order →](#)



N21150559

● PRE-ORDER

Hollow Chop Pliers Plain Box Joint 12.5cm

The smooth concave and convex plain surface allows forming of arch wires easily. Box Joint. Stainless Steel.

[Request pre-order →](#)


N21150432

● IN STOCK

How Utility Pliers Straight 13.5cm

Fine Serrated 3mm dia beak gives firm grip during use. The thin beak provides excellent visibility during the placement and removal of arch wires and general use. Rivet Joint. Stainless Steel

[View on orthodontalusa.com →](#)


N21150430

● IN STOCK

How Utility Pliers angled 13.5cm

Fine Serrated 3mm dia beak gives firm grip during use. The thin beak provides excellent visibility during the placement and removal of arch wires and general use angled beak gives more utility for easy access. Rivet Joint.

[View on orthodontalusa.com →](#)


N21150431

● IN STOCK

How Utility Pliers angled 13cm

Fine Serrated 4mm dia beak gives firm grip during use. The thin beak provides excellent visibility during the placement and removal of arch wires and general use angled beak gives more utility for easy access. Rivet Joint.

[View on orthodontalusa.com →](#)


N21150427

● IN STOCK

Adhesive Removing Pliers 13.5cm

Tips design allows access to remove adhesive in practically any area of the mouth. Double ended carbide rod tip is replaceable for repeated uses for long time. Replaceable Opposing pad tip with stand high heat sterilization.

[View on orthodontalusa.com →](#)



N21150959

● PRE-ORDER

Adhesive Removing Pliers 13cm

Tips design allows access to remove adhesive in practically any area of the mouth. Double ended stainless steel rod tip is replaceable for repeated uses for long time.

[Request pre-order →](#)


N21150415

● IN STOCK

Band Crimping Pliers 12.5cm

This plier formation makes the work of operator to crimp and reform bands with ease and accuracy. Rivet Joint. Stainless Steel

[View on orthodontalusa.com →](#)


N21150418

● IN STOCK

Crown Crimping Pliers 12cm

Ideal Plier for contouring gingival surface of preformed bands to better fit tooth anatomy. Rivet Joint. Stainless Steel

[View on orthodontalusa.com →](#)


N21150814

● PRE-ORDER

Tube Crimping Plier (BB finish 1/3 Polish)

Max: .022" (.56mm) x .025" (.64mm) Wires

The pliers is designed for double back or triple back bending. River Joint. Stainless Steel.

[Request pre-order →](#)


N21150402

● PRE-ORDER

Wire Bending Plier with 1 groove 12.5cm

Ideal for precise bending and forming loops in wires. Box Joint. Stainless Steel

[Request pre-order →](#)



N21150401

● PRE-ORDER

Wire Bending Plier with 1 groove 12.5cm with T/C inserts Tips

Ideal for precise bending and forming loops in wires. Box Joint. Stainless Steel

[Request pre-order →](#)

11

SECTION

Universal, Adams & Specialty Pliers

23

MODELS

11

IN STOCK

12

PRE-ORDER



N21150399

IN STOCK

Universal Pliers Flat Pattern 15.5cm

Universal Pliers for wire and pin roughening, Cutting and bending. Stainless Steel
Bends Max: .038" (.97mm) Cuts Max: .028" (.71mm)

Order: sales@orthodentalusa.com


N21150400

IN STOCK

Universal Pliers 16cm

Universal Pliers for wire and pin roughening, Cutting and bending. Stainless Steel
Bends Max: .038" (.97mm) Cuts Max: .028" (.71mm)

View on orthodentalusa.com →


N21150404

IN STOCK

Wald Sachs Universal Pliers 15cm

This versatile plier may be used as Pin and Ligature, Clasp Bending and Cutting Pliers. Box Joint. Stainless Steel

View on orthodentalusa.com →


N21150064

IN STOCK

Adams Pliers No. 64. 14.5cm

Max: .036" (.91mm)

Utility Pliers for wire bending of different wire Sizes. Box Joint. Stainless Steel.

View on orthodentalusa.com →

★ CERTIFIED QUALITY

✓ ISO 9001:2008

✓ ISO 13485:2016

✓ CE CA006210



N21150065

● PRE-ORDER

Adams Pliers No. 65. 15cm

Max: .036" (.91 mm)

Utility Pliers for wire bending of different wire sizes. Box Joint. Stainless Steel

[Request pre-order →](#)


N21150625

● PRE-ORDER

Adams Pliers 12.5cm

Utility Pliers for wire bending of different Sizes. Rivet Joint. Stainless Steel.

[Request pre-order →](#)


N21150626

● PRE-ORDER

Universal Pliers Box Joint, 14cm

Universal pliers for making different sizes loops.and wire bending. Box Joint. Stainless Steel

[Request pre-order →](#)


N21150801

● PRE-ORDER

Ligature Forming Pliers 13cm

Easily bends ligature wire into a performed ligature, Box Joint. Stainless Steel.

[Request pre-order →](#)


N21150804

● IN STOCK

Ligature Forming Pliers Improved Design 13cm

Easily bends ligature wire into a performed ligature, Formed Ligature Wire Height is 6mm Box Joint. Stainless Steel.

[View on orthodontalusa.com →](#)



N21150632

● PRE-ORDER

Four Step Loop Forming Pliers 12.5cm

Four Step Pliers designed for automatically forming Loops on round wire up to .025" Wire one easy motion of the Plier. Designed to form accurate to consistent Loops 2mm, 3mm, 4mm & 5mm in Height, Box Joint. Stainless Steel.

[Request pre-order →](#)


N21150906

● PRE-ORDER

Band Slitter Pliers 12cm

This cutter has diamond shaped groove to safely trap band during removal. Diamond honed cutting surfaces. Rivet Joint. Stainless Steel.

[Request pre-order →](#)


N21150950

● PRE-ORDER

Renold Pliers 12.5cm

Ideal plier curved beak is useful when forming and contouring wires. Rivet Joint. Stainless Steel

[Request pre-order →](#)


N21150917

● IN STOCK

Pesso Pliers 13cm

Ideal plier long beak one side flat and one side round for forming collars. Rivet Joint. Stainless Steel

[Order: sales@orthodentalusa.com](mailto:sales@orthodentalusa.com)


N21150406

● IN STOCK

Johnson Contouring Pliers 13cm

With Curvature of the beaks only on one side. For contouring and adapting bands and crowns. Rivet Joint. Stainless Steel

[View on orthodentalusa.com →](https://orthodentalusa.com)



N21150912

● PRE-ORDER

Hook Crimping Pliers

Rivet Joint 12cm. Designed to crimp stops, hooks, and posts to archwires. The slot fits over the base of the hook while holding it tight while crimping. Rivet Joint. Stainless Steel. For all type of Hooks

[Request pre-order →](#)


N21150573

● PRE-ORDER

Hook Crimping Pliers 11.5 cm

Designed to crimp stops, hooks, and posts to archwires. The slot fits over the base of the hook while holding it tight while crimping. Rivet Joint. Stainless Steel. For all type of Hooks

[Request pre-order →](#)


N21150901

● IN STOCK

Hook Crimping Pliers 12.5cm

Designed to crimp stops, hooks, and posts to archwires. The slot fits over the base of the hook while holding it tight while crimping. Box Joint. Stainless Steel. For all type of Hooks


[View on orthodontalusa.com →](#)


N21150569

● PRE-ORDER

Hook Crimping Pliers Box joint 13cm

Designed to crimp stops, hooks, and posts to archwires. The slot fits over the base of the hook while holding it tight while crimping Special finish BB & 1/3 Handles Polish. Box Joint. Stainless Steel. For all type of Hooks

[Request pre-order →](#)


N21150904

● IN STOCK

Hook Crimping Pliers Cvd 12cm

Designed to crimp stops, hooks, and posts to archwires. The slot fits over the base of the hook while holding it tight while crimping. Angled beak allows for precise placement of hook on archwire. Rivet Joint. Stainless Steel

[View on orthodontalusa.com →](#)



N21150907

IN STOCK

Band Contouring Pliers Spring Action 12.5cm

Special design plier formation makes the work of operator to crimp and reform bands with ease and accuracy. Concave and Convex beak side gives accurate and firm gripping. Spring formation gives easy handling. Rivet Joint.

Order: sales@orthodentalusa.com


N21150910

IN STOCK

Stop "V" Bend Pliers 12.5cm

Max: .019" (.48mm) x .025" (.64mm) Wire

Special design to make 1mm "V" bends to shorten arch wires or provide a positive stop. Excellent for placing stops in nickel titanium wires. Rivet Joint. Stainless Steel.

View on orthodentalusa.com →


N21150890

PRE-ORDER

Stop V bend Pliers 12.5cm

Max: .019" (.48mm) x .025" (.64mm) Wire

Special design to make 1mm "V" bends to shorten arch wires or provide a positive stop. Excellent for placing stops in nickel titanium wires. Rivet Joint. Stainless Steel.

Request pre-order →



N21150574

PRE-ORDER

NITI Cinch Back Pliers Fine Beak 12.5cm

Max: .021" (.54 mm) x .25" (.64mm)

This plier formation makes the work of cinching back the end of the wire easy from a difficult access lingual area. Slim Beak makes the working more comfortable. Rivet Joint. Stainless Steel NiTi Wires

Request pre-order →

12

SECTION

Cinch-Back, Step & Torquing Pliers

27

MODELS

11

IN STOCK

16

PRE-ORDER



N21150802

IN STOCK

Niti Cinch Back Pliers. Slim Beak 14cm

Max: .012" (.31 mm) x .018" (.46mm)

This plier formation makes the work of cinching back the end of the wire easy from a difficult access lingual area. Slim Beak makes the working more comfortable. Rivet Joint. Stainless Steel NiTi Wires

[View on orthodontalusa.com →](#)


N21150810

PRE-ORDER

Niti Cinch Back Pliers 13.5cm. Short Handle

Max: .012" (.31 mm) x .018" (.46mm)

This plier formation makes the work of cinching back the end of the wire easy from a difficult access lingual area. Slim Beak makes the working more comfortable. Fine serration for firm gripping of wire. Rivet Joint.

[Request pre-order →](#)


N21150809

PRE-ORDER

Niti Cinch Back Pliers Short Head Long Handles 15.5cm

Max: .012" (.31 mm) x .018" (.46mm)

This plier formation makes the work of cinching back the end of the wire easy from a difficult access lingual area. Slim Beak makes the working more comfortable. Fine serration for firm gripping of wire. Rivet Joint.

[Request pre-order →](#)


N21150437

PRE-ORDER

Ceramic Bracket Remover 45 degree angled 12cm

Ideal to remove delicate ceramic bracket with its sharp edges and 45 degree angled. Rivet Joint. Stainless Steel.

[Request pre-order →](#)

★ CERTIFIED QUALITY

✓ ISO 9001:2008

✓ ISO 13485:2016

✓ CE CA006210



N21150875

IN STOCK

Ceramic Bracket remover with groove STR 13cm

Max: .022" (.56mm) x .028" (.71mm)

Sharp beak wedge can easily go under ceramic bracket base while the arch wire is still attached lifts off bracket without breaking and no discomfort to patient due to firm round grip assures,

[View on orthodontalusa.com →](#)


N21150876

IN STOCK

Ceramic Bracket Removing plier Angled 13cm

Sharp beak wedge can easily go under ceramic bracket base while the arch wire is still attached lifts of bracket without bracking and no discomfort to patient due to firm round wire assures,

[View on orthodontalusa.com →](#)


N21150930

IN STOCK

Cap Removing Pliers 12.5cm

Designed to remove caps of the convertible tubes easily without damaging bracket, fine and firm fit of blade beak provides easy working for the specific work Rivet Joint. Stainless Steel

[View on orthodontalusa.com →](#)


N21150930B

IN STOCK

REPLACEMENT BLADE

Order: sales@orthodontalusa.com


N21150633

PRE-ORDER

Double Step Bending Pliers "Z" Bending short 12.5 cm short bend 4.5 mm

Designed to put a Z bend root tip back bend (with one easy squeeze) in the Archwire without removing the Archwire. Also know as 2nd order step bending short bend. 4.5mm. Box Joint. Stainless Steel.

[Request pre-order →](#)



N21150634

● PRE-ORDER

Double Step Bending Pliers "Z" Bending long 12.5 cm long bend 7 mm

Designed to put a Z bend root tip back bend (with one easy squeeze) in the Archwire without removing the Archwire also know as 2nd order step bending long bend. 7mm. Box Joint. Stainless Steel.

[Request pre-order →](#)


N21150903

● IN STOCK

Individual torquing Pliers Set with key 12 cm

Specially designed to place torquing bend into the small sections of archwires (such as a bend for a single tooth) without distorting adjacent sections of the wire. Two piece set consist of one torquing plier and one Key.

[View on orthodontalusa.com →](#)


N21150905 (6mm)

● IN STOCK

Individual torquing Pliers Set (Male & Female) 6 mm Torque 12 cm

These Pliers allow the application of torque without distortion in other segments of the wire. Non-Slip grip tips facilitate the placement of labial or lingual torque. Tungsten Carbide Tips inserts, Rivet Joint. Stainless Steel.

[View on orthodontalusa.com →](#)


N21150900 (3.5mm)

● PRE-ORDER

Individual torquing Pliers Set (Male & Female) 3.5 mm Torque 12 cm

These Pliers allow the application of torque without distortion in other segments of the wire. Non-Slip grip tips facilitate the placement of labial or lingual torque. Tungsten Carbide Tips inserts, Rivet Joint. Stainless Steel.

[Request pre-order →](#)


N21150635

● PRE-ORDER

Individual Torquing Pliers Set. To make 6mm long torque 12.5cm

These Pliers allow the application of torque without distortion in other segments of the wire. Non-Slip grip, tapered tips facilitate the placement of labial or lingual torque. Box Joint. Stainless Steel.

[Request pre-order →](#)



N21150636

● PRE-ORDER

Individual Torquing Pliers Set. To make 3.5mm long torque 12.5cm

These Pliers allow the application of torque without distortion in other segments of the wire. Non-Slip grip, tapered tips facilitate the placement of labial or lingual torque. Box Joint. Stainless Steel.

[Request pre-order →](#)


N21150624

● PRE-ORDER

Individual Torque Forming Pliers Set, To Make 4mm Long Torque, 12.5cm

These Pliers allow the application of torque without distortion in other segments of the wire. Non-Slip grip tips facilitate the placement of labial or lingual torque. Angle marking gives confidence to bend till required angle,

[Request pre-order →](#)


N21150935

● PRE-ORDER

Aero clasp bending 12cm

Specialy designed Pliers for making bend to Aero clasp without putting extra force. Box Joint. Stainless Steel.

[Request pre-order →](#)


N21150936

● PRE-ORDER

Aero clasp bending star shape 12cm

Specialy designed Pliers for making Aero clasp on wires with easy & comfort Box Joint. Stainless Steel.

[Request pre-order →](#)


N21150913

● IN STOCK

Plaster Cutter 20cm

Specialy designed cutter for Plaster Cutting Rivet Joint. Stainless Steel.

[View on orthodontalusa.com →](#)



N21150585

IN STOCK

Elastic Separating Pliers 15.5cm

Angulated Beaks for accessibility, ease & accuracy. Grooved beaks assure positive hold for positioning elastics. Stainless Steel.

[View on orthodontalusa.com →](#)


N21150590

PRE-ORDER

Elastic Separating Pliers without stop 16cm

Angulated Beaks for accessibility, ease & accuracy. Grooved beaks assure positive hold for positioning elastics. Stainless Steel.

[Request pre-order →](#)


N21150587

PRE-ORDER

Elastic Separating Pliers Slim Long 17.5cm

Angulated Beaks for accessibility, ease & accuracy. Grooved beaks assure positive hold for positioning elastics. Stainless Steel.

[Request pre-order →](#)


N21150589

PRE-ORDER

Ligature Tying Pliers 17cm

For Tying Stainless Steel Ligature Wires. made of Stainless Steel.

[Request pre-order →](#)


N21150586

IN STOCK

Coon Style Ligature Tying Pliers 16.5cm

For Tying Stainless Steel Ligature Wires. made of Stainless Steel.

[View on orthodontalusa.com →](#)



N21150710

IN STOCK

Step Forming Pliers 12cm, Straight (90°)

Max: .022" (.56mm) Wires

Forms finishing step in archwires. Excellent pliers for bending utility archwires, also known as first order step bending. Rivet joint. Stainless steel. Kindly specify step size when ordering.

0.25 mm 0.50 mm 0.75 mm 1.00 mm 2.00 mm 3.00 mm

in stock pre-order



N21150723

PRE-ORDER

Step Forming Pliers 12cm, Angled (45°)

Max: .022" (.56mm) Wires

Forms finishing step in archwires. Excellent pliers for bending utility archwires, also known as first order step bending. Rivet joint. Stainless steel. Kindly specify step size when ordering.

0.25 mm 0.50 mm 0.75 mm 1.00 mm 2.00 mm

in stock pre-order



N21150730

PRE-ORDER

Step Forming Pliers 12cm, Angle Size Series

Max: .022" (.56mm) Wires

Forms finishing step in archwires. Excellent pliers for bending utility archwires, also known as first order step bending. Rivet joint. Stainless steel. Kindly specify angle size when ordering.

0.25 mm 0.50 mm 0.75 mm 1.00 mm 1.50 mm 2.00 mm

in stock pre-order

13

SECTION

Lingual Series Pliers

7

MODELS

6

IN STOCK

1

PRE-ORDER



N2LINGUL323

IN STOCK

Flush Cut Distal End Cutter Safety Hold. 15.5cm

Min: .016" (0.41mm) • Max: .021" (0.54mm) x .025" (.64mm) Hard wires

The Lingual series flush cut feature of this instrument allows the archwire to be cut flush with the end of the buccal tube, eliminating the possibility of patient discomfort resulting from archwire protrusion. Box Joint.

[View on orthodontalusa.com](https://orthodontalusa.com) →


N2LINGUL323S

PRE-ORDER

Flush Cut Distal End Cutter Safety Hold. 15.5cm

Min: .012" (0.31mm) • Max: .016" (0.41mm) x .025" (.64mm) Hard wires

The Lingual series flush cut feature of this instrument allows the archwire to be cut flush with the end of the buccal tube, eliminating the possibility of patient discomfort resulting from archwire protrusion. Box Joint.

[Request pre-order](#) →


N2LINGUL331

IN STOCK

Micro-Miniature Ligature Cutter 15.5cm

Max: .012" (0.31mm) Soft Wires

This Lingual series Miniature Ligature cutter small, narrow & sharp pointed tips allow easy access to the most difficult areas. Its feature is to ensure precision cutting of ligature wire and elastics. Box Joint. Stainless Steel

[View on orthodontalusa.com](https://orthodontalusa.com) →


N2LINGUL322

IN STOCK

Miniature Ligature Cutter 60 degree angled 16cm

Max: .012" (0.31mm) Soft Wires

This Lingual series Miniature Ligature cutter small, narrow, sharp pointed and 60 degree angled tips allow easy access to the most difficult areas of lower arch.

[View on orthodontalusa.com](https://orthodontalusa.com) →

★ CERTIFIED QUALITY

✓ ISO 9001:2008

✓ ISO 13485:2016

✓ CE CA006210



N2LINGUL377

IN STOCK

Weingart 60 degree angled 17cm

This Lingual series weingart pliers has the slim long shape serrated beak tips angled at 60 degree which allow easy access in the posterior region. Box Joint. Stainless Steel

[View on orthodontalusa.com →](#)


N2LINGUL380

IN STOCK

Bracket Removing Pliers angled 16cm

The Lingual series bracket remover has the narrow shape tips, wedge between the edges of the bracket base & tooth surface to remove direct bond brackets safely & comfortably. Box Joint. Stainless Steel

Order: sales@orthodontalusa.com


N2LINGUL802

IN STOCK

Niti Cinch Back Pliers 16.5cm

This Lingual series niti cinch back pliers has the serrated tips which hold the wire securely while bending the wire 90 degree distal to the buccal tube, avoiding patient discomfort. Box Joint.

[View on orthodontalusa.com →](#)

14

SECTION
Aligner Pliers

19

MODELS

4

IN STOCK

15

PRE-ORDER



N2119T015

● PRE-ORDER

Vertical Slot Forming Pliers 12 cm

The Vertical emphasizes over-correction of rotations by producing an indentation at the mesial or distal line angles on the facial and/or lingual of a specific tooth.

[Request pre-order →](#)


N2119T016

● PRE-ORDER

Horizontal Slot Forming Pliers 12 cm

The Horizontal instrument produces an indentation to accent individual root torque. It also serves a second purpose by enhancing retention of either clear aligners or retainers by placing an indent into a tooth's undercut or...

[Request pre-order →](#)


N2119T017

● PRE-ORDER

Tear Drop Aligner Forming Pliers 12 cm

specially designed to prevent distortion of the aligner without the need to heat the plier. The ultra slim tips provide easy access incisally, even with lower anterior aligners.

[Request pre-order →](#)


N2119T018

● PRE-ORDER

Hole Punch Aligner Forming Pliers 12.5 cm

specially designed to prevent distortion of the aligner without the need to heat the plier. The ultra slim tips provide easy access incisally, even with lower anterior aligners.

[Request pre-order →](#)

★ CERTIFIED QUALITY

✓ ISO 9001:2008

✓ ISO 13485:2016

✓ CE CA006210



N2119T011

IN STOCK

Bottom Forming Aligner Pliers 1mm, 13.5cm

Creates space in the appliance without cutting a window in the appliance.

[View on orthodontalusa.com →](#)


N2119T012

IN STOCK

Slot Forming Aligner Pliers 13cm

Micro-Ramp Thermopliers move anterior teeth laterally.

[View on orthodontalusa.com →](#)


N2119T013

IN STOCK

Circular Forming Aligner Pliers 13cm

create a platform which the elastic hook is extended

[View on orthodontalusa.com →](#)


N2119T014

IN STOCK

Slot Aligner Forming Plier 13cm

The Slot- forming plier take the second step in enabling elastics to be attached

[View on orthodontalusa.com →](#)


N21150770

PRE-ORDER

Aligner Pliers Slim Hole Punch 1.5mm, short 13cm

Designed to place a round hole in any location on the aligner with one simple squeeze. When placing a hole in the aligner it is not necessary to heat the plier or the aligner, saving time chairside.

[Request pre-order →](#)



N21150772

● PRE-ORDER

Aligner Pliers Hole Punch 3.00mm 13cm

specially designed to prevent distortion of the aligner without the need to heat the plier. The ultra slim tips provide easy access incisally, even with lower anterior aligners.

[Request pre-order →](#)


N21150771

● PRE-ORDER

Aligner Pliers Slim Hole Punch 4.00 mm 13cm

specially designed to prevent distortion of the aligner without the need to heat the plier. The ultra slim tips provide easy access incisally, even with lower anterior aligners.

[Request pre-order →](#)


N21150773

● PRE-ORDER

Aligner Pliers Hole Punch 6.00mm 13cm

specially designed to prevent distortion of the aligner without the need to heat the plier. The ultra slim tips provide easy access incisally, even with lower anterior aligners.

[Request pre-order →](#)


N21150775

● PRE-ORDER

Aligner Dimple Forming 1mm 14cm

is used on the gingival margin at the interdental papilla to create a circular indentation of 1mm for increased retention of aligners and retainers. This may be an essential requirement for teeth with short clinical crowns,

[Request pre-order →](#)


N21150776

● PRE-ORDER

Aligner Dimple Forming 3mm 13cm

is used on the gingival margin at the interdental papilla to create a circular indentation of 1mm for increased retention of aligners and retainers. This may be an essential requirement for teeth with short clinical crowns,

[Request pre-order →](#)



N21150777

● PRE-ORDER

Aligner Vertical Dimple Slot Froming 13.5cm

handles perfectly the material flexibility at the gingival margins, where achieving torque sometimes may be challenging. It creates a horizontal indentation of 0.75 by 1.5mm and creates a force point for individual lingual or...

[Request pre-order →](#)


N21150778

● PRE-ORDER

Aligner Horizontal Dimple Slot Froming 13.5cm

handles perfectly the material flexibility at the gingival margins, where achieving torque sometimes may be challenging. It creates a horizontal indentation of 0.75 by 1.5mm and creates a force point for individual lingual or...

[Request pre-order →](#)


N21150779

● PRE-ORDER

Aligner Opening Instruments 13 cm

Creates an open space on the aligner with a simple squeeze, Tooth movement can occur allowing the aligner not to seat properly.

[Request pre-order →](#)


N21150774

● PRE-ORDER

Aligner Pliers slot Punch Width 4.5mm / Length 10mm, 13.5cm

specially designed to prevent distortion of the aligner without the need to heat the plier. The ultra slim tips provide easy access incisally, even with lower anterior aligners.

[Request pre-order →](#)


N21150780

● PRE-ORDER

Dimple Remover 13cm

The flat surface of the eraser plier reduces or flattens the dimple to remove the pressure.

[Request pre-order →](#)

15

SECTION

Mathieu Needle Holders & Mosquito Forceps

29

MODELS

22

IN STOCK

7

PRE-ORDER



N21150465FT/C

IN STOCK

Mathieu Needle Holder 14cm. Fine beak. T/C

Finer serration T/C Tips for gripping O rings and needles while precision working. Fitted with T/C Tips inserts. Made of stainless steel.

[View on orthodontalusa.com →](#)


N21150466

IN STOCK

Mathieu Needle Holder 14cm. Regular beak with T/C Tips

Regular tips with wider serrated T/C Tips inserts for gripping well while working. Made of stainless steel.

[View on orthodontalusa.com →](#)


N21150468

IN STOCK

Mathieu Needle Holder 14cm. wide beak with T/C tips

Wide beak and fine serration T/C Tips inserts combination makes this needle holder more popular for working. Made of stainless steel.

[View on orthodontalusa.com →](#)


N21150553

IN STOCK

Mathieu Needle Holder 14cm. Offset ratchets. With T/C tips

T/C tips inserts and regular firm beak ensure non slipping grip for working precisely. Offset ratchet makes it more comfortable in working. Made of stainless steel.

Order: sales@orthodontalusa.com

★ CERTIFIED QUALITY

✓ ISO 9001:2008

✓ ISO 13485:2016

✓ CE CA006210



N21150555

IN STOCK

Mathieu Needle Holder Long 7" With T/C tips

T/C tips inserts and wide firm beak ensure non slipping grip for working precisely. Made of stainless steel.

Order: sales@orthodentalusa.com



N21150465FR

IN STOCK

Mathieu Needle Holder 14cm. Fine beak. With Ring

Finer serration for gripping O rings and needles while precision working. Utility of ring makes it easier for working. Made of stainless steel.

Order: sales@orthodentalusa.com



N21150465

IN STOCK

Mathieu Needle Holder 14cm. Regular Tips

Regular tips with wider serration for gripping well while working. Made of stainless steel.

[View on orthodentalusa.com](#) →



N21150465F

IN STOCK

Mathieu Needle Holder 14cm. Fine beak

Finer serration for gripping O rings and needles while precision working. Made of stainless steel.

[View on orthodentalusa.com](#) →



N21150465EXF

IN STOCK

Mathieu Needle Holder 14cm. Extra Fine beak

Super Fine serrated beak for finest working and gripping O rings and needles while precision working. Made of stainless steel.

[View on orthodentalusa.com](#) →



N21150467

IN STOCK

Mathieu Needle Holder 14cm. wide beak and fine serration

Wide beak and fine serration combination makes this needle holder more popular for working. Made of stainless steel.

[View on orthodontalusa.com →](#)


N21150461

IN STOCK

Mathieu Needle Holder 14cm. Hook Tip

This hook tip style needle holder has a notched tip for easy and assurance in placing elastic. Non slipping grip. Made of stainless steel.

[View on orthodontalusa.com →](#)


N21150462

IN STOCK

Mathieu Needle Holder 14cm. Hole Tip

This hole tip style needle holder features a groove at tip for non slip placement of elastic ligatures and other intraoral elastics. Made of stainless steel.

[View on orthodontalusa.com →](#)


N21150463

IN STOCK

Mathieu Needle Holder 14cm. Tips contact

This tips contact style needle holder has tightest grip on elastic ties. It is not recommended for ligature wires. Made of stainless steel.

[View on orthodontalusa.com →](#)


N21150469

IN STOCK

Mosquito Mathieu Needle Holder 14cm

Straight serration and fine beak ensure non slipping grip for working precisely. Made of stainless steel.

[View on orthodontalusa.com →](#)



N21150545

● PRE-ORDER

Mathieu Needle Holder TC Plain Handle 14cm

Narrow beak and fine serration T/C Tips inserts combination makes this needle holder more popular for working. Plain handles gives more comfort Made of stainless steel.

[Request pre-order →](#)


N21150546

● PRE-ORDER

Samaha Mathieu Needle Holder 14 cm

Straight serration and fine beak ensure non slipping grip for working precisely. Made of stainless steel.

[Request pre-order →](#)


N21150547

● PRE-ORDER

Samaha Mathieu Needle Holder Double Spring Action 14 cm

Straight serration and fine beak ensure non slipping grip for working precisely. Made of stainless steel.

[Request pre-order →](#)


N21150551

● IN STOCK

Mosquito Mathieu Needle Holder 14cm. Double Spring action. Straight serration

and fine beak ensure non slipping grip for working precisely. Made of stainless steel.

[View on orthodontalusa.com →](#)


N21150552

● IN STOCK

Mathieu Needle Holder 14cm. Offset

Cross serration and regular firm beak ensure non slipping grip for working precisely. Offset ratchet makes it more comfortable in working. Made of stainless steel.

[Order: sales@orthodontalusa.com](mailto:sales@orthodontalusa.com)



N21150554

IN STOCK

Mathieu Needle Holder 14 cm Mathieu Needle Holder 14cm. Offset ratchets, Double action Offset ratchets....

Cross serration and fine firm beak ensure non slipping grip for working precisely. Offset ratchet makes it more comfortable in working. Light pattern makes this Mathieu working soft and

[Order: sales@orthodentalusa.com](mailto:sales@orthodentalusa.com)


N21150752

PRE-ORDER

Mathieu Needle Holder 14cm, Pointed Beak

Cross serration and regular firm beak ensure non slipping grip for working precisely. Offset ratchet makes it more comfortable in working. Made of stainless steel.

[Request pre-order →](#)


NA-17140

IN STOCK

Mosquito Forceps 12.5cm Straight

Made of stainless steel.

[View on orthodentalusa.com →](http://orthodentalusa.com)


N21150753

PRE-ORDER

Mathieu Needle Holder 14cm, Plain Handles Pointed Fine Beak

Cross serration and regular firm beak ensure non slipping grip for working precisely. Offset ratchet makes it more comfortable in working. Made of stainless steel.

[Request pre-order →](#)


NA-17142H

IN STOCK

Mosquito Forceps 12.5cm Hook Tips

Made of stainless steel

[View on orthodentalusa.com →](http://orthodentalusa.com)



NA - 17142

IN STOCK

Mosquito Forceps 12.5cm Curved

Made of stainless steel.

[View on orthodontalusa.com →](#)


NA - 17140T

IN STOCK

Mosquito Forceps 12.5cm Tip Contact

Made of stainless steel.

[View on orthodontalusa.com →](#)


NA - 17140E

IN STOCK

Mosquito Forceps 12.5cm Hole Tips

Made of stainless steel.

[View on orthodontalusa.com →](#)


NA - 17144

PRE-ORDER

Mosquito Forceps Curved 12.5cm Light And Smooth Working

Made of stainless steel.

[Request pre-order →](#)


NA - 17143

PRE-ORDER

Mosquito Forceps Straight Fine 12.5cm Light And Smooth Working

Made of stainless steel.

[Request pre-order →](#)

16

SECTION

Tweezers & Cotton Pliers

16

MODELS

14

IN STOCK

2

PRE-ORDER



N21150450

IN STOCK

BONDING TWEEZERS 11.5cm. Made of stainless

steel, cross action. For placement of brackets with thin.025" plate on the rear side.

[View on orthodontalusa.com →](#)


N21150451

IN STOCK

BONDING TWEEZERS 14cm Heavy

Made of stainless steel, cross action. For placement of brackets.

[View on orthodontalusa.com →](#)


N21150452

IN STOCK

Bracket Placer/Slot Aligner Tweezers 13cm

Made of stainless steel, cross action. For placement & lineament of brackets.

[View on orthodontalusa.com →](#)


N21150453

IN STOCK

IMPROVED TWEEZERS 13cm. Made of stainless

steel, cross action. For placement of brackets with thin.025" plate on the rear side

[View on orthodontalusa.com →](#)

★ CERTIFIED QUALITY

✓ ISO 9001:2008

✓ ISO 13485:2016

✓ CE CA006210



N21150454

IN STOCK

Keat Buccal Tweezer 11.5cm

Made of stainless steel, For placement of buccal tubes.

[View on orthodontalusa.com →](#)


N21150455

IN STOCK

Buccal Tube Tweezer Cross action 13cm

Made of Stainless steel, cross action. For placement buccal tubes.

[View on orthodontalusa.com →](#)


N21150456

IN STOCK

Buccal Tube Tweezer Cross action with slot aligner 14.5cm

Made of stainless steel, For placement of buccal tubes and aligning them easily.

[View on orthodontalusa.com →](#)


N21150457

IN STOCK

BONDING Tweezers with slot aligner 13.5cm

Made of stainless steel, cross action. For placement of brackets and aligning them.

[View on orthodontalusa.com →](#)


N21150510

IN STOCK

Lingual Bracket Tweezers 14cm

Made of stainless steel, cross action. For placement of brackets.

[View on orthodontalusa.com →](#)



N22255/2

IN STOCK

College Cotton Pliers 15cm

Made of stainless steel, serrated jaws.

Order: sales@orthodentalusa.com



N22255/2WL

IN STOCK

College Cotton Pliers with locking 15cm

Made of stainless steel, serrated jaws.

Order: sales@orthodentalusa.com



N22265/2

IN STOCK

Meriam Cotton Pliers 16cm

Made of stainless steel, serrated jaws.

Order: sales@orthodentalusa.com



N21150458

IN STOCK

Soldering Tweezers 16.5cm

Made of stainless steel ,with Fiber Thumb Pad for protection, curved tips.

[View on orthodentalusa.com](#) →



N21150460

IN STOCK

ARTICULATING PAPER Forceps 15cm

Made of stainless steel, cross action firmly hold paper.

[View on orthodentalusa.com](#) →



N21150315

● PRE-ORDER

Bracket Holding and aligning (serrated) Tweezers extra long 16.5cm

Made of Stainless Steel

[Request pre-order →](#)


N21150459

● PRE-ORDER

Lingual Bracket and Tube Holder Cross Action 16cm

Made of Stainless Steel

[Request pre-order →](#)

17

SECTION

Scalers & Orthodontic Instruments

30

MODELS

20

IN STOCK

10

PRE-ORDER



N21150470

IN STOCK

Double-Ended Ligature Director 17.5cm

Handle made of Non-Magnet Hollow Stainless steel pipe. Corrosion Resistant for longer use. Points made of hardened Stainless Steel

[View on orthodontalusa.com →](#)


N21150471

IN STOCK

Ligature Director Scaler 17.5cm

Handle made of Non-magnet Hollow Stainless steel pipe. Corrosion Resistant for longer use. Points made of hardened Stainless steel

[View on orthodontalusa.com →](#)


N21150472

IN STOCK

Band Pusher Scaler 17cm

Handle made of Non-Magnet Hollow Stainless steel pipe. Corrosion Resistant for longer use. Points made of hardened Stainless steel

[View on orthodontalusa.com →](#)


N21150473

IN STOCK

Schure Band Seater 16.5cm

Round Handle. One part instrument made of hardened stainless steel

[View on orthodontalusa.com →](#)

★ CERTIFIED QUALITY

✓ ISO 9001:2008

✓ ISO 13485:2016

✓ CE CA006210



N21150476

IN STOCK

Schure Band Seater 16.5cm

Handle made of non-magnet Hollow stainless steel pipe. Corrosion resistant for longer use. Points Made of hardened stainless steel.

[View on orthodontalusa.com →](#)


N21150484

IN STOCK

Tucker and Pusher 17.5cm

Handle made of non-magnet octagonal Hollow stainless steel pipe. Corrosion resistant for longer use. Points Made of hardened stainless steel.

[View on orthodontalusa.com →](#)


N21150474

IN STOCK

TWIRL-ON Module 17cm

Handle made of Non-magnet Hollow stainless steel pipe. Corrosion Resistant for longer use. Points made of hardened stainless steel.

[View on orthodontalusa.com →](#)


N21150475

IN STOCK

Distal Bender Single End 15cm

for wire up to 22" Simple using made from durable stainless steel, bends tails of springs. Work on all Stainless steel

[View on orthodontalusa.com →](#)


N21150475D/E

IN STOCK

Distal Bender Double End 16.5cm

one side to accomodate wire upto .019" and other side to accommodate wires up to .030" simple using made from durable stainless steel, bends tails of spring. Works on all stainless steel wires.

[View on orthodontalusa.com →](#)



N21150478

IN STOCK

Ligature Director curved, fine narrow 17.5cm

Octagonal Handle. Narrow Groove. One part instrument made of hardened stainless steel

[View on orthodontalusa.com →](#)


N21150479

IN STOCK

Ligature Director curved, 17.5cm

Regular Made of hardened Stainless steel

[View on orthodontalusa.com →](#)


N21150482

IN STOCK

Ligature Director Double Ended, 18cm

Octagonal Handle. Wide groove. Firm Grip. One part instrument made of hardened stainless steel

[View on orthodontalusa.com →](#)


N21150483

PRE-ORDER

2D Opener 15cm

Round Handle. One part instrument made of hardened stainless steel

[Request pre-order →](#)


N21150535

IN STOCK

Stainless steel Spatula No. 8

Double Ended 18.5cm.

Order: sales@orthodontalusa.com



N21150536

IN STOCK

Stainless steel Spatula No.3

Double Ended 16.5cm

[View on orthodontalusa.com →](#)


N21150481

IN STOCK

Shortly Twister Holder 15.5cm. Designed for sample insertion and locking of shortly ligature

twists and shortly koby twists. Easy spinning body for tight fit. Made of stainless steel.

[View on orthodontalusa.com →](#)


N21150485

IN STOCK

Tension Mapping Tool Measures the force applied in specific orthodontic applications

Sterilizable. Measures upto 16oz for max force. Made of stainless steel

[View on orthodontalusa.com →](#)


N21150486

IN STOCK

Nickel Titanium High Heat Ice Instruments 11cm

Made of stainless steel

[View on orthodontalusa.com →](#)



N21150678

IN STOCK

Positive Positioner 16.5cm

Designed to hold brackets firmly and align them positively. Made of Aluminum.

[View on orthodontalusa.com →](#)



N21150680

IN STOCK

Cavity Filler Tube 13cm

Made of Stainless Steel

[View on orthodontalusa.com →](#)


N21150487

IN STOCK

SpinTek Instruments. Bracket Cap slider. (Ormco Type) 12.5cm

Made of Stainless steel. Capable of being sterilize.

[View on orthodontalusa.com →](#)


N21150504

PRE-ORDER

Ligature Director Double Point & Opener Hollow Handle Wide Round 16cm

Made of stainless steel

[Request pre-order →](#)


N21150507

PRE-ORDER

Band Pusher & Scaler Octagonal Wide Hollow Handle 17cm

Made of stainless steel

[Request pre-order →](#)


N21150508

PRE-ORDER

Band Pusher & Scaler Hollow Handle 17.5cm

Made of stainless steel

[Request pre-order →](#)

N21150509

● PRE-ORDER

Probe & Opener, Solid Handle 16cm

Made of stainless steel.


[Request pre-order →](#)

N21150525

● PRE-ORDER

Orthodontic probe slim short point Double Ended Solid Handle 17.5cm

Made of stainless steel


[Request pre-order →](#)

N21150843

● PRE-ORDER

Band Pusher & Scaler Hollow Handle 17.5cm

Made of stainless steel


[Request pre-order →](#)

N21150844

● PRE-ORDER

Ligature Director with Scaler Hollow Handle 18cm

Made of stainless steel


[Request pre-order →](#)

N21150845

● PRE-ORDER

Band Pusher with Ligature Director Hollow Handle 17cm

Made of stainless steel


[Request pre-order →](#)



N21150846

● PRE-ORDER

Band Pusher with Scaler Hollow Handle 17cm

Made of stainless steel

[Request pre-order →](#)

18

SECTION

Gauges & Mapping Tools

33

MODELS

11

IN STOCK

22

PRE-ORDER



N21150673

● PRE-ORDER

Bracket Mapping Tool Aluminum 15 cm

To work on: .018" Slots

[Request pre-order →](#)


N21150675

● IN STOCK

Bracket Mapping Tool Aluminum 15 cm

To work on: .022" Slots

[View on orthodontalusa.com →](#)


N21150675SS

● IN STOCK

Bracket Mapping Tool Stainless Steel 15 cm

To work on: .018" & .022" Slots

Order: sales@orthodontalusa.com


N21150674

● IN STOCK

Bracket Mapping Tool Fully Adjustable 15 cm

Made of Stainless steel. Capable of being sterilizable. To work on: .018" & .022" Slots

[View on orthodontalusa.com →](#)

★ CERTIFIED QUALITY

✓ ISO 9001:2008

✓ ISO 13485:2016

✓ CE CA006210

N21150703

● PRE-ORDER

Mapping Tool Stainless Steel 13.5cm

To Work on .022" Slots.


[Request pre-order →](#)

N21150701

● PRE-ORDER

Mapping Tool Stainless Steel 13.5cm

To Work on .022" Slots.



[Request pre-order →](#)

N21150702

● PRE-ORDER

Mapping Tool Stainless Steel 13.5cm

To Work on .022" Slots.


[Request pre-order →](#)

N21150699

● PRE-ORDER

Mapping Tool Stainless Steel 13cm

To Work on .022" Slots.


[Request pre-order →](#)

N21150697

● PRE-ORDER

Mapping Tool Stainless Steel 13cm

To Work on .018" Slots.


[Request pre-order →](#)

N21150693

● PRE-ORDER

Mapping Tool Stainless Steel 10.5cm

To Work on .022" Slots.


[Request pre-order →](#)

N21150691

● PRE-ORDER

Mapping Tool Stainless Steel 10.5cm

To Work on .022" Slots.


[Request pre-order →](#)

N21150698

● PRE-ORDER

Mapping Tool Stainless Steel 10.5cm

To Work on .022" Slots.


[Request pre-order →](#)

N21150704

● PRE-ORDER

Mapping Tool Stainless Steel 10.5cm

To Work on .022" Slots.


[Request pre-order →](#)

N21150696

● PRE-ORDER

Mapping Tool Stainless Steel 10.5cm

To Work on .022" Slots.


[Request pre-order →](#)




N21150695

● PRE-ORDER

Mapping Tool Stainless Steel 10.5cm

To Work on .022" Slots.


[Request pre-order →](#)


N21150690

● PRE-ORDER

Mapping Tool Stainless Steel 10.5cm


To Work on .022" Slots.

[Request pre-order →](#)


N21150685

● PRE-ORDER

Mapping Tool 2.0mm & 2.5mm work on .022" Slots. 11.5cm Mapping Tool 5.5mm & 6.0mm work on .018"...

[Request pre-order →](#)


N21150661

● PRE-ORDER

Mapping Tool 14cm

Set of 4 Pieces made of Aluminum

[Request pre-order →](#)


N21150670N

● IN STOCK

Bracket Mapping tool Star type

Made of Stainless steel. Capable of being sterilize.

[View on orthodontalusa.com →](#)



N21150669

● PRE-ORDER

Mapping Tool

Made of Stainless steel. Capable of being sterilize.

[Request pre-order →](#)


N21150668

● PRE-ORDER

Mapping Tool

Made of Stainless steel. Capable of being sterilize.

[Request pre-order →](#)


N21150662

● PRE-ORDER

Mapping Tool

Made of Aluminum Sizes: 0.10mm, 0.15mm, 0.20mm, 0.25mm, 0.30mm, 0.40mm, 0.50mm,

[Request pre-order →](#)


N21150671

● IN STOCK

Bracket Mapping tool Star type fully adjustable with lead

points made of stainless steel fully adjustable to accommodate

[View on orthodontalusa.com →](#)


N21150670

● IN STOCK

Bracket Mapping tool Star type

Made of Stainless steel. Capable of being sterilize.

[View on orthodontalusa.com →](#)

N21150505

● PRE-ORDER

Mapping Scale 15cm

Made of Stainless steel.


[Request pre-order →](#)

N21150677

● PRE-ORDER

Boley Mapping Tool 14cm

Made of Stainless steel. Capable of being sterilize.


[Request pre-order →](#)

N21150679

● IN STOCK

Breene Donk Mapping Ruler 13 cm

Made of Stainless steel. Capable of being sterilize.



Order: sales@orthodentalusa.com

N21150682

● IN STOCK

Dividing Mapping Tool 11cm

Made of Brass and Stainless Steel Working Points


[View on orthodentalusa.com →](#)

N21150762

● PRE-ORDER

Band Seater/Adaptor 14.5cm


[Request pre-order →](#)



N21150477

IN STOCK

Mershon Band Pusher 14.5cm

All sided angled serrated tips reduced risk of slippage, lightweight hollow handle gives easy working Made of Stainless steel. Capable of being sterilize.

[View on orthodontalusa.com →](#)


N21150480

IN STOCK

Band Setter fully stainless steel 15.5cm

Fully serrated tips reduced risk of slippage, lightweight hollow handle gives easy working Made of Stainless steel. Capable of being sterilize.

[View on orthodontalusa.com →](#)


N21150902

IN STOCK

Interproximal Stripping Instrument 9.5cm

Fully Stainless steel. Durable and Sterilizable. Surely holds lightening strips for controlled reduction.

[View on orthodontalusa.com →](#)


N21150947

PRE-ORDER

Stripper Handle 10.5cm

Made of Stainless steel.

[Request pre-order →](#)

19

SECTION

Spatulas, Knives & Scissors

17

MODELS

12

IN STOCK

5

PRE-ORDER



N21150500

IN STOCK

PLASTER Spatula 21.5cm

Made of stainless steel wooden handle, ridge blade.

[View on orthodontalusa.com →](#)


N21150512

IN STOCK

PLASTER Spatula 21cm

Made of stainless steel wooden handle, flexible blade.

[View on orthodontalusa.com →](#)


N21150520

IN STOCK

FAHNENSTOCK Wax Knife 18cm

Made of stainless steel with wooden handle

[View on orthodontalusa.com →](#)


N21150530

IN STOCK

GRITMAN Plaster Trimmer 18cm

Made of stainless steel, plaster knife with scraper

[View on orthodontalusa.com →](#)

★ CERTIFIED QUALITY

✓ ISO 9001:2008

✓ ISO 13485:2016

✓ CE CA006210



N21150532

IN STOCK

Plaster knife 17.5cm

Made of stainless steel, plaster knife with wooden handle.

[View on orthodontalusa.com →](#)


N21150534

IN STOCK

PLASTER Knife 16.5cm

Made of stainless steel, plaster knife with scraper with wooden handle.

[View on orthodontalusa.com →](#)


N21150640

IN STOCK

Fine Scissors 11cm

Straight 11.5cm Made of stainless steel

[View on orthodontalusa.com →](#)


N21150646

IN STOCK

Crown Scissors Str 11.5cm

Made of Stainless steel. One side serrated blade provides non slipping cutting ability. Fine beak. Capable of being sterilize.

[View on orthodontalusa.com →](#)


N21150700

IN STOCK

Dressing Scissors 14cm

Made of stainless steel, blunt/blunt, sharp/blunt s sharp/sharp, straight models are available

[View on orthodontalusa.com →](#)



N21150647

IN STOCK

Aligner Scissors Curved 11.5cm

Made of Stainless steel. One side serrated blade provides non slipping cutting ability. Fine beak. Capable of being sterilize.

[View on orthodontalusa.com →](#)


N21150650

IN STOCK

Utility Scissors 19cm

Made of stainless steel with molded plastic handles.

[View on orthodontalusa.com →](#)


N21150656

PRE-ORDER

Aligner Scissors Straight 9CM

Serrated blades for extra durable holding while cutting Capable of being sterilize, Made of Stainless steel.

[Request pre-order →](#)


N21150657

PRE-ORDER

Aligner Scissors Curved 9cm

Serrated blades for extra durable holding while cutting Capable of being sterilize, Made of Stainless steel.

[Request pre-order →](#)


N21150644

IN STOCK

Crown Scissors Straight 12CM

Made of Stainless steel. One side Serrated Blade Provides non slipping cutting ability fine beak capable of being sterilize.

[View on orthodontalusa.com →](#)



N21150645

● PRE-ORDER

Crown Scissors Curved 12CM

Made of Stainless steel. One side Serrated Blade Provides non slipping cutting ability fine beak capable of being sterilize.

[Request pre-order →](#)


N21150651

● PRE-ORDER

Aligner Scissors Straight 12.5cm

Serrated blades for extra durable holding while cutting Capable of being sterilize, Made of Stainless steel.

[Request pre-order →](#)


N21150652

● PRE-ORDER

Aligner Scissors Curved 12.5cm

Serrated blades for extra durable holding while cutting Capable of being sterilize, Made of Stainless steel.

[Request pre-order →](#)

20

SECTION

Photography & Intraoral Accessories

23

MODELS

3

IN STOCK

20

PRE-ORDER



N21150908

IN STOCK

Band Seater Stick, Autoclavable, Blue 14.5cm

Autoclavable.

[View on orthodontalusa.com →](#)


N21150518

PRE-ORDER

Cheek Retractor, Plastic, Slim & Curved, Set of 2 Pieces 15.5cm

Autoclavable.

[Request pre-order →](#)


N21150511

IN STOCK

Lip & Cheek Retractors Set 11.5cm

Autoclavable 137°C. Made of stainless steel.

[Order: sales@orthodontalusa.com](mailto:sales@orthodontalusa.com)


N21150521

PRE-ORDER

Black Contraster Set of 6 Pieces

Autoclavable. Stainless steel.

[Request pre-order →](#)

★ CERTIFIED QUALITY

✓ ISO 9001:2008

✓ ISO 13485:2016

✓ CE CA006210



N21150521A

Black Contraster, Shape A

Autoclavable. Stainless steel.

● PRE-ORDER

[Request pre-order →](#)



N21150521B

Black Contraster, Shape B

Autoclavable. Stainless steel.

● PRE-ORDER

[Request pre-order →](#)



N21150521C

Black Contraster, Shape C

Autoclavable. Stainless steel.

● PRE-ORDER

[Request pre-order →](#)



N21150521D

Black Contraster, Shape D

Autoclavable. Stainless steel.

● PRE-ORDER

[Request pre-order →](#)



N21150521E

Black Contraster, Shape E

Autoclavable. Stainless steel.

● PRE-ORDER

[Request pre-order →](#)



N21150521F

● PRE-ORDER

Black Contraster, Shape F

Autoclavable, Stainless steel.

[Request pre-order →](#)



N21150980

● PRE-ORDER

Adult Occlusal Child Occlusal Photographic Mirrors(shape A)

Stainless Steel. Autoclavable

[Request pre-order →](#)



N21150981

● PRE-ORDER

Buccal Lateral Photographic Mirrors (shape B)

Stainless Steel. Autoclavable

[Request pre-order →](#)



N21150982

● PRE-ORDER

Adult Occlusal Child Occlusal Photographic Mirrors(shape C)

Stainless Steel. Autoclavable

[Request pre-order →](#)



N21150983

● PRE-ORDER

Adult Buccal Adult Occlusal Photographic Mirrors(shape D)

Stainless Steel. Autoclavable

[Request pre-order →](#)

	N21150984 ● PRE-ORDER Lingual Upper Lingual Lower Photographic Mirrors(shape E) Stainless Steel. Autoclavable Request pre-order →
	N21150985 ● PRE-ORDER Child Buccal Adult Buccal Photographic Mirrors(shape F) Stainless Steel. Autoclavable Request pre-order →
	N21150986 ● PRE-ORDER Child Occlusal Child Buccal Photographic Mirrors(shape G) Stainless Steel. Autoclavable Request pre-order →
	N21150988 ● PRE-ORDER Occlusal-XL Adult (Shape I) Stainless Steel. Autoclavable Front Surface Request pre-order →
	N21150987 ● PRE-ORDER Buccal (Shape H) Stainless Steel. Autoclavable Request pre-order →



N21150683

Front Surface Mirror Size 9 with Handle

Made of Stainless Steel

IN STOCK

[View on orthodontalusa.com →](#)



N21150940

Bite Stick, Lining Handle 14.5cm

Autoclavable 137°C. Extra durable and firm strength. Available in five colors.

Blue
Red
Orange
Purple
Green

■ in stock
 ■ pre-order

[Request pre-order →](#)



N21150945

Bite Stick, Serrated Handle 14.5cm

Autoclavable 137°C. Extra durable and firm strength. Available in four colors.

Blue
Red
Purple
Orange

■ in stock
 ■ pre-order

[Request pre-order →](#)



N21150987

Bite Stick, Plain Handle 14.5cm

Autoclavable. Available in six colors.

Blue
Red
Purple
Orange
Green
White

■ in stock
 ■ pre-order

[Request pre-order →](#)

21

SECTION Sterilization, Trays & Kits

15

MODELS

12

IN STOCK

3

PRE-ORDER



N21150920

IN STOCK

Impression Trays PERFORATED Sizes: L1, L2, L3, L4, L5 & U1, U2, U3, U4, U5 Individual numbers as per

requirement are also available. Kindly mention while ordering Made of non Magnet Stainless Steel

[View on orthodontalusa.com →](#)


CT-B12-006

PRE-ORDER

Sterilization Cassettes Tray 28cm

For 12 Instruments, Box Type Made of Stainless Steel

[Request pre-order →](#)


N21150921

IN STOCK

Impression Trays SOLID Sizes: L1, L2, L3, L4, L5 & U1, U2, U3, U4, U5 Individual numbers as per

requirement are also available. Kindly mention while ordering Made of non Magnet Stainless Steel

[View on orthodontalusa.com →](#)


CT-B8-005

PRE-ORDER

Sterilization Cassettes Tray 21cm

For 8 Instruments, Box Type Made of Stainless Steel

[Request pre-order →](#)

★ CERTIFIED QUALITY

✓ ISO 9001:2008

✓ ISO 13485:2016

✓ CE CA006210



CT-B15-003

IN STOCK

Sterilization Cassettes Tray 20cm

For 15 Instruments, Box Type Made of Stainless Steel

[View on orthodontalusa.com →](#)


CT-B20-004

IN STOCK

Sterilization Cassettes Tray 29cm

For 20 Instruments, Box Type

[View on orthodontalusa.com →](#)


CT-B7-002

IN STOCK

Sterilization Cassettes Tray 21cm

For 7 Instruments, Box Type Made of Stainless Steel

[View on orthodontalusa.com →](#)


N21150705

IN STOCK

Scaler Tray 20cm

Made of Stainless Steel

[View on orthodontalusa.com →](#)


N21150706

IN STOCK

Kidney Tray. Medium Size 21.5cm

Made of Stainless Steel

[View on orthodontalusa.com →](#)



CT-B5-001

IN STOCK

Sterilization Cassettes Tray 18.5cm

For 5 Instruments, Box Type Made of Stainless Steel

[View on orthodontalusa.com →](#)


N21150602

IN STOCK

Pliers Hanging Rack SS with Holding Cup 25cm

Made of non magnetic Stainless steel with additional facility of cup to hold hand instruments

[View on orthodontalusa.com →](#)


N21150601

IN STOCK

Pliers Hanging Rack 25cm

With Acrylic Transparent Plates joining three non magnetic Stainless steel Rods

[View on orthodontalusa.com →](#)


N21150600

IN STOCK

Pliers Hanging Rack SS 22cm

Made of non magnetic Stainless steel

[View on orthodontalusa.com →](#)


KIT-1106

IN STOCK

Basic Dental Examination Kit, Regular Quality

Set of 3 pieces made of stainless steel: mirror #5, explorer and cotton plier.

[View on orthodontalusa.com →](#)



KIT-1107

Basic Dental Examination Kit, High Quality

Set of 3 pieces made of stainless steel.

● PRE-ORDER

[Request pre-order →](#)

22

SECTION

Additional In-Stock Instruments

26

MODELS

25

IN STOCK

1

PRE-ORDER



K-1106-MR

● IN STOCK

Replacement Mirror #5 (for Kit 1106)

Order: sales@orthodentalusa.com



N2-00-006

● IN STOCK

Mirror Handle, Hollow, with Mirror #5

Order: sales@orthodentalusa.com



N2025005

● IN STOCK

#6 Stainless Steel Explorer, Single-Ended

[View on orthodentalusa.com →](#)



N21150047

● IN STOCK

Ligature Director Scaler

[View on orthodentalusa.com →](#)

★ CERTIFIED QUALITY

✓ ISO 9001:2008

✓ ISO 13485:2016

✓ CE CA006210



N211503023

IN STOCK

Facebow Plier / Aderer Pliers, Serrated

Laboratory plier.

[View on orthodontalusa.com →](#)


N211503485

IN STOCK

Micro Ligature & Pin Cutter, 15° Angled 11.5cm

Maximum cutting capacity .012" soft wires.

Order: sales@orthodontalusa.com


N211503485B

IN STOCK

Ligature Cutter, Black

[View on orthodontalusa.com →](#)


N21150374-S

IN STOCK

Set x 5 Studs (Only) for Band Removing Pliers (Pliers NOT Included)

Not for autoclave.

[View on orthodontalusa.com →](#)


N21150389-TC

IN STOCK

Bird Beak Pliers with Cutter, Slim & Sloped Beak, T/C, Grooved

[View on orthodontalusa.com →](#)



N21150389 B

IN STOCK

Bird Beak Pliers with Cutter, Slim & Sloped Beak, No T/C, No Grooves

[View on orthodontalusa.com →](#)


N21150400/TC

IN STOCK

Universal Plier 15.5cm, T/C

[View on orthodontalusa.com →](#)


N21150416

IN STOCK

Cinch Back Pliers

[Order: sales@orthodontalusa.com](mailto:sales@orthodontalusa.com)


N21150419-L

IN STOCK

Hollow Chop Long Beak Arch Forming Plier, Long Tip

[View on orthodontalusa.com →](#)


N21150427S

IN STOCK

Replacement Screws for Adhesive Removing Pliers

[Order: sales@orthodontalusa.com](mailto:sales@orthodontalusa.com)



N21150427T

● IN STOCK

Carbide Blade for Adhesive Removing Pliers

[View on orthodontalusa.com →](#)



N21150427TSS

● IN STOCK

Stainless Steel Blade for Adhesive Removing Pliers

1 blade per bag.

Order: sales@orthodontalusa.com



N21150454L

● IN STOCK

Keat Tube Holding Tweezers

[View on orthodontalusa.com →](#)



N21150465EXF-GH

● IN STOCK

Mathieu Needle Holder 14cm, Super Fine Beak, Gold Handle

[View on orthodontalusa.com →](#)



N21150476N

● IN STOCK

Ligature Director & Pusher, Hollow Handle, Wide

[View on orthodontalusa.com →](#)



N21150551C

● PRE-ORDER

Mathieu Needle Holder, Double Spring, Curved

[View on orthodontalusa.com →](#)



N21150586T

● IN STOCK

Ligature Tying Pliers, T-Type

Order: sales@orthodontalusa.com



N21150688

● IN STOCK

Bracket Height Gauge, Adjustable Type 2mm–3.5mm

Settings: 2mm–3.5mm / 2.5mm–3mm.

[View on orthodontalusa.com →](#)



N21150708

● IN STOCK

Instrument & Accessory Organizer Tray

[View on orthodontalusa.com →](#)



N21150709

● IN STOCK

Silicone Handle Set, 5 pcs

Order: sales@orthodontalusa.com



N21150904T

● IN STOCK

Band Slitter Pliers with T/C Tips

Order: sales@orthodentalusa.com



N2JN0662

● IN STOCK

IPR Increase Gauge, Stainless Steel Set

[View on orthodentalusa.com](https://orthodentalusa.com) →

QUICK REFERENCE INDEX

● In stock ● Pre-order

REF	MODEL	PG	REF	MODEL	PG
N21150317	Distal End Cutter, Universal Cut 12cm	● 5	N21150360	Heavy Wire Cutter.15.5cm	● 20
N21150319	Distal End Cutter, Universal Cut. Long Neck	● 5	N21150380	Bracket Removing Plier 12.5cm angled	● 21
N21150320	Distal End Cutter, Universal Cut. Small Neck	● 5	N21150381	Bracket Removing Plier 13cm Straight	● 21
N21150321	Distal End Cutter, Universal Cut, Small Neck	● 5	N21150383	Bracket Removing Plier 12.5cm angled T/C Tip	● 21
DEN-002	Distal End Cutter, Universal Cut Box Joint 1	● 6	N21150577	Bracket Removing Pliers Angled Box Joint 13c	● 21
N21150326	Distal End Cutter, Universal Cut. Long Neck	● 6	N21150371	Band Removing Pliers long neck 13cm	● 22
N21150327	Distal End Cutter, Universal Cut	● 6	N21150373	Band removing pliers short neck 13cm	● 22
N21150328	Distal End Cutter, Universal Cut. Slim Patte	● 6	N21150374	Band Removing Pliers 13.5cm	● 22
N21150571	Distal End Cutter 12 cm	● 6	N21150812	Band Removing Pliers Long Neck 13.5cm	● 22
N21150565	Distal End Cutter Universal 12.5cm	● 7	N21150824	Band Removing Pliers Long beak with T/C Tip	● 22
N21150606	Distal End Cutter Universal 12.5cm	● 7	N21150578	Band Removing Pliers Box Joint 14cm	● 23
N21150956	Distal End Cutter Universal 11.5cm	● 7	TEF-201	Teflon Studs / Pads	● 23
N21150334	Distal End Cutter, Mini Head 13cm	● 7	PHI-105	High Heat Band Removing Pads	● 23
N21150335	Distal End Cutter, Box Joint 12.5cm	● 7	N21150378	Weingart Pliers 13cm	● 24
N21150317L	Distal End Cutter, Universal Cut. Long Handl	● 8	N21150377	Weingart Pliers 13cm	● 24
N21151000	Distal End Cutter	● 8	N21150379	Weingart Pliers 13.5cm	● 24
N21150606L	Distal End Cutter, Universal 15cm	● 8	N21150395	Weingart Pliers Box Joint 14 cm	● 24
DEN-001	Distal End Cutter, Universal Cut. Long Handl	● 8	N21150396	Slim Weingart Pliers Box Joint 14.5 cm	● 25
N21150318	Distal End Cutter, Universal Cut. Long Slim	● 8	N21150375	Slim Weingart Pliers 13.5cm	● 25
N21150325	Distal End Cutter, Universal Cut. Long Handl	● 9	N21150376	Slim Weingart Pliers 13.5cm	● 25
N205001	Distal End Cutter, Universal Cut. Heavy patt	● 9	N21150638	Slim Weingart Pliers Narrow	● 25
N21150322	Flush Cutter without Hold 45 degree angled 1	● 9	N21150863	Weingart Pliers (MINI) 13cm	● 25
N21150323	Flush Cutter with Hold 12cm	● 9	N21150397	Slim Weingart 90 degree 13.5cm	● 26
N21150323H	Flush Cutter with Hold 13.5cm	● 10	N21150957	Slim Weingart Pliers 13cm	● 26
N21150323F	Flush Cutter with Hold, Heavy pattern, Long	● 10	N21150390	Bird Beak Pliers 12cm	● 27
N21150606F	Flush Cutter with Hold 12.5 cm	● 10	N21150958	Bird Beak Pliers Fine 11.5cm	● 27
N21150324	Slim Flush Cutter with Hold 12.5cm	● 10	N21150576	Bird Beak Pliers Box Joint 12.5cm	● 27
N21150575	Flush Cutter without Hold 12.5cm	● 10	N21150813	Slim Sloped Bird Beak 12.5cm	● 27
N21150345	Flush Cutter with Hold 13 cm	● 11	N21150637	Slim Sloped Bird Beak 12.5cm	● 28
N21150330	Ligature & Pin Cutter straight 12cm	● 12	N21150388	Slim Bird Beak Pliers Tapered Beak 12.5cm	● 28
N21150330L	Ligature & Pin Cutter straight Long Handles1	● 12	N21150390H	Bird Beak Pliers Heavy and serrated jaw 12.5	● 28
N21150342	Ligature & Pin Cutter 15 degree angled 12cm	● 12	N21150822	Lab Bird Beak Pliers with Serrated Tips 14cm	● 28
N21150862	Ligature & Pin Cutter 15 degree Regular Beak	● 12	N21150823	Bird Beak Long Beak for Heavy Wires 12.5cm	● 28
N21150572	Ligature & Pin cutter 12 cm	● 13	N21150631	Light Wire Bending Pliers with grooves tips	● 29
N21150605	Ligature and Pin Cutter 12 cm	● 13	N21150392SB	Light Wire Bending Plier Small Beak 12.4cm	● 29
N21150815	Ligature Cutter Long Handles 7 degree 14cm	● 13	N21150392	Light Wire Bending Pliers Long 12.5cm	● 29
N21150816	Ligature Cutter 7 degree 12.5cm	● 13	N21150390W/C	Bird Beak Pliers with cutters 12cm	● 29
N21150817	Ligature Cutter 15 degree 12.5cm	● 13	N21150389	Bird Beak pliers with cutters 12.5cm Tapered	● 29
N21151002	Ligature & Pin Cutter wide beak 15 degree an	● 14	N21150391	Light wire bending Pliers with cutters 13 cm	● 30
N21150332	Mini Ligature & Pin Cutter 15 degree	● 14	N21150413	TWEED Arch Forming Pliers 12cm	● 31
N21150332L	Mini Ligature & Pin Cutter 15 degree	● 14	N21150417	TWEED Arch Forming Pliers 12.5cm	● 31
N21150333	Ligature & Pin Cutter wide beak 15 degree an	● 14	N21150303	Tweed Arch Bending Pliers	● 31
N21150352	Micro Ligature & Pin Cutter T/C 15 degree an	● 14	N21150412	Slim Arch Forming Pliers 12.5cm	● 31
N21150351	Micro Ligature & Pin Cutter 15 degree light	● 15	N21150411	Rectangular Arch Forming Pliers 12.5cm	● 32
N21150818	Micro Ligature Cutter 7 degree 12cm	● 15	N21150568	Arch forming pliers box joint 13cm	● 32
N21150354	Micro Ligature & Pin Cutter 15 degree angle	● 15	N21150410	Lingual Arch Forming Pliers 12cm	● 32
N21150603	Mini Ligature and pin Cutter 12cm	● 15	N21150405	Lingual Arch Forming Pliers 12cm	● 32
N21150329	Ligature Cutter 45 degree angled wide handle	● 15	N21150407	Arch Forming Pliers 12cm	● 32
N21150336	Ligature & Pin cutter, Box Joint 12.5 cm	● 16	N21150382	Nance Pliers 12cm	● 33
N211505002	Micro Ligature & Pin Cutter 7 degree angle 1	● 16	N21150864	Nance Pliers 12.5cm	● 33
N211505003	Ligature & Pin Cutter straight 13cm	● 16	N21150580	Nance Pliers Box Joint 12.5 cm	● 33
N21150331	Mini Ligature & Pin Cutter Str, Slim	● 16	N21150639	Nance Pliers Double Sided Loop Forming 12cm	● 33
N21150355	Ligature & Pin Cutter str, Light Slim Model	● 16	N21150761	Nance Double Sided Pliers 12cm	● 33
N21150566	Ligature and pin Cutter 12cm	● 17	N21150304	Optical Loop Forming Pliers Wide Handles 12.	● 34
N21150353	Hard Wire Cutter Straight 11.5cm	● 18	N21150911	Optical Pliers 12.5cm	● 34
N21150353A	Hard Wire Cutter 15 degree 11.5cm	● 18	N21150915	Loop Forming Pliers, Small Beak 12cm	● 35
N21150353L	Hard Wire Cutter Str Long Beak 12cm	● 18	N21150916	Loop Forming Pliers, Long Beak 12.5cm	● 35
N21150353H	Hard Wire Cutter 15 degree Long	● 18	N21150955	Ruhland Loop Forming Pliers 12.5cm	● 35
N21150353B	Hard Wire Cutter Black, 15 degree	● 19	N21150385	Three Jaw Pliers 12.5cm	● 35
N21150560	Hard Wire Cutters 13cm	● 19	N21150579	Three Prong Pliers Box Joint 13 cm	● 36
N21150819	Hard Wire Cutter str 12.5cm	● 19	N21150386	Aderer Pliers 13cm	● 36
N21150820	Hard Wire Cutter 15 degree 12.5cm	● 19	N21150305	Three Jaw Pliers Rounded 12cm	● 36
N21150356	Heavy Duty Cutter, Multi-action	● 19	N21150370	NiTi Wire Bender (3 Jaw) Slim	● 36
N21150357	Heavy Duty Cutter, Multi-action. Carbon Stee	● 20	N21150387	Mini Aderer Pliers 11.5cm	● 36
N21150358	Heavy Duty Cutter, Multi-action Mechanism 9"	● 20	N21153023	Lab Heavy Three Prong / Heavy Aderer 14.5cm	● 37
N21150359	Heavy Duty Cutter, Multi-action Mechanism 7"	● 20	N21150422	Tweed Pliers 12.5cm	● 37

QUICK REFERENCE INDEX

● In stock ● Pre-order

REF	MODEL	PG	REF	MODEL	PG
N21150420	Tweed Pliers 12.5cm	● 37	N21150636	Individual Torquing Pliers Set. To make 3.5m	● 53
N21150423	Tweed forming Pliers 13.5 cm	● 37	N21150624	Individual Torque Forming Pliers Set, To Mak	● 53
N21150424	Omega Loop Forming Pliers, Tweed Type 12.5cm	● 37	N21150935	Aero clasp bending 12cm	● 53
N21150425	Young Loop Bending Pliers 12.5cm	● 38	N21150936	Aero clasp bending star shape 12cm	● 53
N21150339	Young Pliers Box Joint 13.5 cm	● 38	N21150913	Plaster Cutter 20cm	● 53
N21150435	Auxiliary Spring Forming Pliers 13cm	● 38	N21150585	Elastic Separating Pliers 15.5cm	● 54
N21150421	Loop Forming Pliers with cutters. Tweed Type	● 38	N21150590	Elastic Separating Pliers without stop 16cm	● 54
N21150426	Kim Pliers	● 38	N21150587	Elastic Separating Pliers Slim Long 17.5cm	● 54
N21150393	Jaraback Pliers 13cm	● 39	N21150589	Ligature Tying Pliers 17cm	● 54
N21150338	Jaraback Pliers Box joint (Slim) 14 cm	● 40	N21150586	Coon Style Ligature Tying Pliers 16.5cm	● 54
N21150393W/C	Jaraback Pliers with cutters 13.5cm	● 40	N21150710	Step Forming Pliers 12cm, Straight (90°)	● 55
N21150394	Jaraback Pliers with T/C cutters 13.5cm	● 40	N21150723	Step Forming Pliers 12cm, Angled (45°)	● 55
N21150408	Arch Forming Pliers 12cm	● 40	N21150730	Step Forming Pliers 12cm, Angle Size Series	● 55
N21150414	Hollow Chop Arch Forming Pliers 12cm	● 41	N2LINGUL323	Flush Cut Distal End Cutter Safety Hold. 15.	● 56
N21150409	Hollow Chop Arch Forming pliers with 3 Groov	● 41	N2LINGUL323S	Flush Cut Distal End Cutter Safety Hold. 15.	● 56
N21150419	Arch Forming Pliers 13.5cm	● 41	N2LINGUL331	Micro-Miniature Ligature Cutter 15.5cm	● 56
N21150436	Wire Bending Pliers wide curve 12.5cm	● 41	N2LINGUL322	Miniature Ligature Cutter 60 degree angled 1	● 56
N21150570	Hollow Chop Pliers with Grooves	● 41	N2LINGUL377	Weingart 60 degree angled 17cm	● 57
N21150559	Hollow Chop Pliers Plain Box Joint 12.5cm	● 42	N2LINGUL380	Bracket Removing Pliers angled 16cm	● 57
N21150432	How Utility Pliers Straight 13.5cm	● 42	N2LINGUL802	Niti Cinch Back Pliers 16.5cm	● 57
N21150430	How Utility Pliers angled 13.5cm	● 42	N2119T015	Vertical Slot Forming Pliers 12 cm	● 58
N21150431	How Utility Pliers angled 13cm	● 42	N2119T016	Horizontal Slot Forming Pliers 12 cm	● 58
N21150427	Adhesive Removing Pliers 13.5cm	● 42	N2119T017	Tear Drop Aligner Forming Pliers 12 cm	● 58
N21150959	Adhesive Removing Pliers 13cm	● 43	N2119T018	Hole Punch Aligner Forming Pliers 12.5 cm	● 58
N21150415	Band Crimping Pliers 12.5cm	● 43	N2119T011	Botton Forming Aligner Pliers 1mm, 13.5cm	● 59
N21150418	Crown Crimping Pliers 12cm	● 43	N2119T012	Slot Forming Aligner Pliers 13cm	● 59
N21150814	Tube Crimping Plier (BB finish 1/3 Polish)	● 43	N2119T013	Circular Forming Aligner Pliers 13cm	● 59
N21150402	Wire Bending Plier with 1 groove 12.5cm	● 43	N2119T014	Slot Aligner Forming Plier 13cm	● 59
N21150401	Wire Bending Plier with 1 groove 12.5cm with	● 44	N21150770	Aligner Pliers Slim Hole Punch 1.5mm, short	● 59
N21150399	Universal Pliers Flat Pattern 15.5cm	● 45	N21150772	Aligner Pliers Hole Punch 3.00mm 13cm	● 60
N21150400	Universal Pliers 16cm	● 45	N21150771	Aligner Pliers Slim Hole Punch 4.00 mm 13cm	● 60
N21150404	Wald Sachs Universal Pliers 15cm	● 45	N21150773	Aligner Pliers Hole Punch 6.00mm 13cm	● 60
N21150064	Adams Pliers No. 64. 14.5cm	● 45	N21150775	Aligner Dimple Forming 1mm 14cm	● 60
N21150065	Adams Pliers No. 65. 15cm	● 46	N21150776	Aligner Dimple Forming 3mm 13cm	● 60
N21150625	Adams Pliers 12.5cm	● 46	N21150777	Aligner Vertical Dimple Slot Froming 13.5cm	● 61
N21150626	Universal Pliers Box Joint, 14cm	● 46	N21150778	Aligner Horizontal Dimple Slot Froming 13.5c	● 61
N21150801	Ligature Forming Pliers 13cm	● 46	N21150779	Aligner Opening Instruments 13 cm	● 61
N21150804	Ligature Forming Pliers Improved Design 13cm	● 46	N21150774	Aligner Pliers slot Punch Width 4.5mm / Leng	● 61
N21150632	Four Step Loop Forming Pliers 12.5cm	● 47	N21150780	Dimple Remover 13cm	● 61
N21150906	Band Slitter Pliers 12cm	● 47	N21150465FT/C	Mathieu Needle Holder 14cm. Fine beak. T/C	● 62
N21150950	Renold Pliers 12.5cm	● 47	N21150466	Mathieu Needle Holder 14cm. Regular beak wit	● 62
N21150917	Pesso Pliers 13cm	● 47	N21150468	Mathieu Needle Holder 14cm. wide beak with T	● 62
N21150406	Johnson Contouring Pliers 13cm	● 47	N21150553	Mathieu Needle Holder 14cm. Offset ratchets.	● 62
N21150912	Hook Crimping Pliers	● 48	N21150555	Mathieu Needle Holder Long 7" With T/C tips	● 63
N21150573	Hook Crimping Pliers 11.5 cm	● 48	N21150465FR	Mathieu Needle Holder 14cm. Fine beak. With	● 63
N21150901	Hook Crimping Pliers 12.5cm	● 48	N21150465	Mathieu Needle Holder 14cm. Regular Tips	● 63
N21150569	Hook Crimping Pliers Box joint 13cm	● 48	N21150465F	Mathieu Needle Holder 14cm. Fine beak	● 63
N21150904	Hook Crimping Pliers Cvd 12cm	● 48	N21150465EXF	Mathieu Needle Holder 14cm. Extra Fine beak	● 63
N21150907	Band Contouring Pliers Spring Action 12.5cm	● 49	N21150467	Mathieu Needle Holder 14cm. wide beak and fi	● 64
N21150910	Stop "V" Bend Pliers 12.5cm	● 49	N21150461	Mathieu Needle Holder 14cm. Hook Tip	● 64
N21150890	Stop V bend Pliers 12.5cm	● 49	N21150462	Mathieu Needle Holder 14cm. Hole Tip	● 64
N21150574	NITI Cinch Back Pliers Fine Beak 12.5cm	● 49	N21150463	Mathieu Needle Holder 14cm. Tips contact	● 64
N21150802	Niti Cinch Back Pliers. Slim Beak 14cm	● 50	N21150469	Mosquito Mathieu Needle Holder 14cm	● 64
N21150810	Niti Cinch Back Pliers 13.5cm. Short Handle	● 50	N21150545	Mathieu Needle Holder TC Plain Handle 14cm	● 65
N21150809	Niti Cinch Back Pliers Short Head Long Handl	● 50	N21150546	Samaha Mathieu Needle Holder 14 cm	● 65
N21150437	Ceramic Bracket Remover 45 degree angled 12c	● 50	N21150547	Samaha Mathieu Needle Holder Double Spring A	● 65
N21150875	Ceramic Bracket remover with groove STR 13cm	● 51	N21150551	Mosquito Mathieu Needle Holder 14cm. Double	● 65
N21150876	Ceramic Bracket Removing plier Angled 13cm	● 51	N21150552	Mathieu Needle Holder 14cm. Offset	● 65
N21150930	Cap Removing Pliers 12.5cm	● 51	N21150554	Mathieu Needle Holder 14 cm Mathieu Needle H	● 66
N21150930B	REPLACEMENT BLADE	● 51	N21150752	Mathieu Needle Holder 14cm, Pointed Beak	● 66
N21150633	Double Step Bending Pliers " Z" Bending shor	● 51	NA-17140	Mosquito Forceps 12.5cm Straight	● 66
N21150634	Double Step Bending Pliers " Z" Bending long	● 52	N21150753	Mathieu Needle Holder 14cm, Plain Handles Po	● 66
N21150903	Individual torquing Pliers Set with key 12 c	● 52	NA-17142H	Mosquito Forceps 12.5cm Hook Tips	● 66
N21150905 (6mm)	Individual torquing Pliers Set (Male & Femal	● 52	NA-17142	Mosquito Forceps 12.5cm Curved	● 67
N21150900 (3.5mm)	Individual torquing Pliers Set (Male & Femal	● 52	NA-17140T	Mosquito Forceps 12.5cm Tip Contact	● 67
N21150635	Individual Torquing Pliers Set. To make 6mm	● 52	NA-17140E	Mosquito Forceps 12.5cm Hole Tips	● 67

QUICK REFERENCE INDEX

● In stock ● Pre-order

REF	MODEL	PG	REF	MODEL	PG
NA-17144	Mosquito Forceps Curved 12.5cm Light And Smo	● 67	N21150670N	Bracket Mapping tool Star type	● 82
NA-17143	Mosquito Forceps Straight Fine 12.5cm Light	● 67	N21150669	Mapping Tool	● 83
N21150450	BONDING TWEEZERS 11.5cm. Made of stainless	● 68	N21150668	Mapping Tool	● 83
N21150451	BONDING TWEEZERS 14cm Heavy	● 68	N21150662	Mapping Tool	● 83
N21150452	Bracket Placer/Slot Aligner Tweezers 13cm	● 68	N21150671	Bracket Mapping tool Star type fully adjusta	● 83
N21150453	IMPROVED TWEEZERS 13cm. Made of stainless	● 68	N21150670	Bracket Mapping tool Star type	● 83
N21150454	Keat Buccal Tweezer 11.5cm	● 69	N21150505	Mapping Scale 15cm	● 84
N21150455	Buccal Tube Tweezer Cross action 13cm	● 69	N21150677	Boley Mapping Tool 14cm	● 84
N21150456	Buccal Tube Tweezer Cross action with slot a	● 69	N21150679	Breene Donk Mapping Ruler 13 cm	● 84
N21150457	BONDING Tweezers with slot aligner 13.5cm	● 69	N21150682	Dividing Mapping Tool 11cm	● 84
N21150510	Lingual Bracket Tweezers 14cm	● 69	N21150762	Band Seater/Adaptor 14.5cm	● 84
N22255/2	College Cotton Pliers 15cm	● 70	N21150477	Mershon Band Pusher 14.5cm	● 85
N22255/2WL	College Cotton Pliers with locking 15cm	● 70	N21150480	Band Setter fully stainless steel 15.5cm	● 85
N22265/2	Meriam Cotton Pliers 16cm	● 70	N21150902	Interproximal Stripping Instrument 9.5cm	● 85
N21150458	Soldering Tweezers 16.5cm	● 70	N21150947	Stripper Handle 10.5cm	● 85
N21150460	ARTICULATING PAPER Forceps 15cm	● 70	N21150500	PLASTER Spatula 21.5cm	● 86
N21150315	Bracket Holding and aligning (serrated) Twe	● 71	N21150512	PLASTER Spatula 21cm	● 86
N21150459	Lingual Bracket and Tube Holder Cross Action	● 71	N21150520	FAHNENSTOCK Wax Knife 18cm	● 86
N21150470	Double-Ended Ligature Director 17.5cm	● 72	N21150530	GRITMAN Plaster Trimmer 18cm	● 86
N21150471	Ligature Director Scaler 17.5cm	● 72	N21150532	Plaster knife 17.5cm	● 87
N21150472	Band Pusher Scaler 17cm	● 72	N21150534	PLASTER Knife 16.5cm	● 87
N21150473	Schure Band Seater 16.5cm	● 72	N21150640	Fine Scissors 11cm	● 87
N21150476	Schure Band Seater 16.5cm	● 73	N21150646	Crown Scissors Str 11.5cm	● 87
N21150484	Tucker and Pusher 17.5cm	● 73	N21150700	Dressing Scissors 14cm	● 87
N21150474	TWIRL-ON Module 17cm	● 73	N21150647	Aligner Scissors Curved 11.5cm	● 88
N21150475	Distal Bender Single End 15cm	● 73	N21150650	Utility Scissors 19cm	● 88
N21150475D/E	Distal Bender Double End 16.5cm	● 73	N21150656	Aligner Scissors Straight 9CM	● 88
N21150478	Ligature Director curved, fine narrow 17.5cm	● 74	N21150657	Aligner Scissors Curved 9cm	● 88
N21150479	Ligature Director curved, 17.5cm	● 74	N21150644	Crown Scissors Straight 12CM	● 88
N21150482	Ligature Director Double Ended, 18cm	● 74	N21150645	Crown Scissors Curved 12CM	● 89
N21150483	2D Opener 15cm	● 74	N21150651	Aligner Scissors Straight 12.5cm	● 89
N21150535	Stainless steel Spatula No. 8	● 74	N21150652	Aligner Scissors Curved 12.5cm	● 89
N21150536	Stainless steel Spatula No.3	● 75	N21150908	Band Seater Stick, Autoclavable, Blue 14.5cm	● 90
N21150481	Shortly Twister Holder 15.5cm. Designed for	● 75	N21150518	Cheek Retractor, Plastic, Slim & Curved, Set	● 90
N21150485	Tension Mapping Tool Measures the force appl	● 75	N21150511	Lip & Cheek Retractors Set 11.5cm	● 90
N21150486	Nickel Titanium High Heat Ice Instruments 11	● 75	N21150521	Black Contraster Set of 6 Pieces	● 90
N21150678	Positive Positioner 16.5cm	● 75	N21150521A	Black Contraster, Shape A	● 91
N21150680	Cavity Filler Tube 13cm	● 76	N21150521B	Black Contraster, Shape B	● 91
N21150487	SpinTek Instruments. Bracket Cap slider. (Or	● 76	N21150521C	Black Contraster, Shape C	● 91
N21150504	Ligature Director Double Point & Opener Holl	● 76	N21150521D	Black Contraster, Shape D	● 91
N21150507	Band Pusher & Scaler Octagonal Wide Hollow H	● 76	N21150521E	Black Contraster, Shape E	● 91
N21150508	Band Pusher & Scaler Hollow Handle 17.5cm	● 76	N21150521F	Black Contraster, Shape F	● 92
N21150509	Probe & Opener, Solid Handle 16cm	● 77	N21150980	Adult Occlusal Child Occlusal Photographic M	● 92
N21150525	Orthodontic probe slim short point Double En	● 77	N21150981	Buccal Lateral Photographic Mirrors (shape B	● 92
N21150843	Band Pusher & Scaler Hollow Handle 17.5cm	● 77	N21150982	Adult Occlusal Child Occlusal Photographic M	● 92
N21150844	Ligature Director with Scaler Hollow Handle	● 77	N21150983	Adult Buccal Adult Occlusal Photographic Mir	● 92
N21150845	Band Pusher with Ligature Director Hollow Ha	● 77	N21150984	Lingual Upper Lingual Lower Photographic Mir	● 93
N21150846	Band Pusher with Scaler Hollow Handle 17cm	● 78	N21150985	Child Buccal Adult Buccal Photographic Mirro	● 93
N21150673	Bracket Mapping Tool Aluminum 15 cm	● 79	N21150986	Child Occlusal Child Buccal Photographic Mir	● 93
N21150675	Bracket Mapping Tool Aluminum 15 cm	● 79	N21150988	Occlusal-XL Adult (Shape I)	● 93
N21150675SS	Bracket Mapping Tool Stainless Steel 15 cm	● 79	N21150987	Buccal (Shape H)	● 94
N21150674	Bracket Mapping Tool Fully Adjustable 15 cm	● 79	N21150683	Front Surface Mirror Size 9 with Handle	● 94
N21150703	Mapping Tool Stainless Steel 13.5cm	● 80	N21150940	Bite Stick, Lining Handle 14.5cm	● 94
N21150701	Mapping Tool Stainless Steel 13.5cm	● 80	N21150945	Bite Stick, Serrated Handle 14.5cm	● 94
N21150702	Mapping Tool Stainless Steel 13.5cm	● 80	N21150987	Bite Stick, Plain Handle 14.5cm	● 94
N21150699	Mapping Tool Stainless Steel 13cm	● 80	N21150920	Impression Trays PERFORATED Sizes: L1, L2, L	● 95
N21150697	Mapping Tool Stainless Steel 13cm	● 80	CT-B12-006	Sterilization Cassettes Tray 28cm	● 95
N21150693	Mapping Tool Stainless Steel 10.5cm	● 81	N21150921	Impression Trays SOLID Sizes: L1, L2, L3, L4	● 95
N21150691	Mapping Tool Stainless Steel 10.5cm	● 81	CT-B8-005	Sterilization Cassettes Tray 21cm	● 95
N21150698	Mapping Tool Stainless Steel 10.5cm	● 81	CT-B15-003	Sterilization Cassettes Tray 20cm	● 96
N21150704	Mapping Tool Stainless Steel 10.5cm	● 81	CT-B20-004	Sterilization Cassettes Tray 29cm	● 96
N21150696	Mapping Tool Stainless Steel 10.5cm	● 81	CT-B7-002	Sterilization Cassettes Tray 21cm	● 96
N21150695	Mapping Tool Stainless Steel 10.5cm	● 82	N21150705	Scaler Tray 20cm	● 96
N21150690	Mapping Tool Stainless Steel 10.5cm	● 82	N21150706	Kidney Tray. Medium Size 21.5cm	● 96
N21150685	Mapping Tool 2.0mm & 2.5mm work on.022" Slot	● 82	CT-B5-001	Sterilization Cassettes Tray 18.5cm	● 97
N21150661	Mapping Tool 14cm	● 82	N21150602	Pliers Hanging Rack SS with Holding Cup 25cm	● 97

QUICK REFERENCE INDEX

REF	MODEL	PG
N21150601	Pliers Hanging Rack 25cm	97
N21150600	Pliers Hanging Rack SS 22cm	97
KIT-1106	Basic Dental Examination Kit, Regular Qualit	97
KIT-1107	Basic Dental Examination Kit, High Quality	98
K-1106-MR	Replacement Mirror #5 (for Kit 1106)	99
N2-00-006	Mirror Handle, Hollow, with Mirror #5	99
N2025005	#6 Stainless Steel Explorer, Single-Ended	99
N21150047	Ligature Director Scaler	99
N211503023	Facebow Plier / Aderer Pliers, Serrated	100
N21150348S	Micro Ligature & Pin Cutter, 15° Angled 11.5	100
N21150348SB	Ligature Cutter, Black	100
N21150374-S	Set x 5 Studs (Only) for Band Removing Plier	100
N21150389-TC	Bird Beak Pliers with Cutter, Slim & Sloped	100
N21150389 B	Bird Beak Pliers with Cutter, Slim & Sloped	101
N21150400/TC	Universal Plier 15.5cm, T/C	101
N21150416	Cinch Back Pliers	101
N21150419-L	Hollow Chop Long Beak Arch Forming Plier, Lo	101
N21150427S	Replacement Screws for Adhesive Removing Pli	101
N21150427T	Carbide Blade for Adhesive Removing Pliers	102
N21150427TSS	Stainless Steel Blade for Adhesive Removing	102
N21150454L	Keat Tube Holding Tweezers	102
N21150465EXF-GH	Mathieu Needle Holder 14cm, Super Fine Beak,	102
N21150476N	Ligature Director & Pusher, Hollow Handle, W	102
N21150551C	Mathieu Needle Holder, Double Spring, Curved	103
N21150586T	Ligature Tying Pliers, T-Type	103
N21150688	Bracket Height Gauge, Adjustable Type 2mm-3.	103
N21150708	Instrument & Accessory Organizer Tray	103
N21150709	Silicone Handle Set, 5 pcs	103
N21150904T	Band Slitter Pliers with T/C Tips	104
N2JN0662	IPR Increase Gauge, Stainless Steel Set	104

● In stock ● Pre-order

REF	MODEL	PG
-----	-------	----

PLIER CARE & MAINTENANCE

Make your instruments last for years.

Your pliers are an investment. Most performance problems come from sterilization and cleaning mistakes — these simple rules keep them cutting and bending like new.

★ CERTIFIED QUALITY

✓ ISO 9001:2008

✓ ISO 13485:2016

✓ CE CA006210

DO

- Pre-clean with surgical milk or distilled water.
- Use distilled water in the ultrasonic cleaner and change the solution daily.
- Keep pliers separated from each other in the ultrasonic cleaner.
- Clean and sterilize with instruments in the **open position**.
- Lubricate hinged instruments before sterilization with lubricants made for precision joints.
- Match the plier to the wire: check each model's cutting capacity.

DON'T

- Don't use soap and water, or glutaraldehyde solutions.
- Don't disinfect with phenols.
- Don't scrub with hard brushes — they scratch the finish.
- Don't over-sterilize; it can dull the polish.
- Don't use cold sterilization — it can chemically attack instruments.
- Don't cut wire heavier than the stated capacity; it damages the inserted tips.

Sterilization

135°C / 275°F

5 minutes · steam autoclave

121°C / 250°F

20 minutes · steam autoclave

Max 177°C / 350°F

never heat above

Warranty

Our pliers are covered against defects in material and manufacturing workmanship under normal use for **one year** from the date of purchase. Defective pliers will be repaired or replaced at our discretion. Misuse, negligence, or lack of proper maintenance and care voids this warranty. Regular sharpening is recommended to maintain cutting ability — contact us for service.

For use by dental professionals only. Clean, inspect and sterilize before each use.

TIPS

Choosing the right plier

Cutting wire

Distal end cutters with safety hold keep the cut wire end from becoming a hazard. Flush cutters trim archwires flush with the buccal tube. Ligature & pin cutters handle soft wire up to .012"; hard-wire and heavy-duty cutters are built for rigid and heavy archwires.

Bending & forming

Bird beak, light wire and Jaraback pliers make precise bends; Tweed, loop-forming and Nance pliers shape loops and arches; step pliers and Z-bend pliers create offsets and root-tip bends.

Removing

Bracket removing pliers lift bonded brackets safely; band removing pliers work with replaceable high-heat pads; ceramic bracket removers and adhesive removing pliers protect enamel.

Holding & placing

Mathieu needle holders place ligatures and elastics; tweezers and bonding tweezers position brackets and tubes; band pushers and ligature directors seat bands and tuck wire.

Aligner care

Aligner forming pliers create precise slots, tear-drop and hole-punch features that prevent distortion of the aligner.

Need help choosing?

Call +1 (954) 655-8487 or email sales@orthodentalusa.com — we will match the right instrument to your technique.



Let's equip your practice.

Can we help you? Our team answers availability questions, quotes and pre-orders.

WEBSITE

orthodentalusa.com

PHONE / WHATSAPP

+1 (954) 655-8487

EMAIL

sales@orthodentalusa.com

ADDRESS

1500 Weston Rd Suite 204, Weston, FL
33326, USA

★ CERTIFIED QUALITY

✓ ISO 9001:2008

✓ ISO 13485:2016

✓ CE CA006210